

## MD Biosciences

Division of Morwell Diagnostics GmbH  
Gewerbestr. 9 8132 Egg b. Zürich Switzerland  
Tel: +41-1 986 26 28 Fax: +41-1 986 26 30

For USA/ Canada : (651) 641 1770 Fax: (651) 641 1773

For Israel Toll Free: 1-800-200 MDB Fax: (8) 940 9973

For United Kingdom (UK) : 0800 0197 200

e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

# **Cell Proliferation kit with XTT Reagent** *for 1000 assays*

**Cat. No.:** 409005 (for 1000 Assays)

**Storage:** -20°C

## **1. Introduction**

Cell proliferation assays are widely used in cell biology for the study of growth factors, cytokines and media components, for the screening of cytotoxic agents and for lymphocyte activation. The need for a reliable, sensitive and quantitative assay that would enable analysis of a large number of samples led to the development of methods, such as:

- use of radioactive thymidine to label DNA in live cells
- use of Brdu to label DNA in live cells (as a substitute for radioactive thymidine)

The above methods have a number of disadvantages, including: use of radioactive materials and relatively complex techniques. The use of tetrazolium salts, such as MTT, commenced in the 1950s and is based on the fact that live cells reduce tetrazolium salts into colored formazan compounds. The biochemical procedure is based on the activity of mitochondria enzymes which are inactivated shortly after cell death. This method was found to be very efficient in assessing the viability of cells. A colorimetric method based on the tetrazolium salt, XTT, was first described by P.A. Scudiero in 1988. Whilst the use of MTT produced a non-soluble formazan compound which necessitated dissolving the dye in order to measure it, the use of XTT produces a soluble dye. The use of XTT greatly simplifies the procedure of measuring proliferation, and is, therefore, an excellent solution to the quantitating of cells and their viability without using radioactive isotopes. This kit was developed to assay cell proliferation in reaction to different growth factors, cytokines and nutrient components. In addition, it is suitable for assaying cytotoxicity of materials such as TNF or other growth inhibitors. XTT can be used as a non-radioactive substitute for cytotoxic tests based on the release of <sup>51</sup>Cr from cells with no less sensitivity.



## MD Biosciences

Division of Morwell Diagnostics GmbH

Gewerbestr. 9 8132 Egg b. Zürich Switzerland

Tel: +41-1 986 26 28 Fax: +41-1 986 26 30

For USA/ Canada : (651) 641 1770 Fax: (651) 641 1773

For Israel Toll Free: 1-800-200 MDB Fax: (8) 940 9973

For United Kingdom (UK) : 0800 0197 200

e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

The advantages of using this kit can be summarized with the following attributes:

- *easy-to-use*: there is no requirement for additional reagents and/or the cell washing procedures
- *speed*: multiwell plates and an ELISA reader can be used for reading
- *sensitivity*: can be assayed even in low cell concentrations
- *accuracy*: dye absorbance is proportional to the number of cells in each well
- *safety*: there is no need for radioactive isotopes

## 2. **Kit Components**

### 2.1 *XTT Reagent (10x5ml)*

A sterile solution containing the XTT reagent. The solution should be stored frozen and should not be exposed to light. To avoid repeated re-freezing, dividing the solution into a number of vials after defrosting the original vial is recommended.

Note: *if sediment is present in the solution, heat the solution to 37°C and swirl gently until a clear solution is obtained.*

### 2.2 *Activation Reagent (2x0.5ml)*

A sterile solution containing PMS (N-methyl dibenzopyrazine methyl sulfate). The solution should be stored frozen and should not be exposed to light. To avoid repeated re-freezing, dividing the solution into a number of vials after defrosting the original vial is recommended.

Note: *if sediment is present in the solution, heat the solution to 37°C and swirl gently until a clear solution is obtained.*

## **MD Biosciences**

Division of Morwell Diagnostics GmbH  
Gewerbestr. 9 8132 Egg b. Zürich Switzerland  
Tel: +41-1 986 26 28 Fax: +41-1 986 26 30

For USA/ Canada : (651) 641 1770 Fax: (651) 641 1773

For Israel Toll Free: 1-800-200 MDB Fax: (8) 940 9973

For United Kingdom (UK) : 0800 0197 200

e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

### **Assay Principles**

The assay is based on the ability of metabolic active cells to reduce the tetrazolium salt XTT to orange colored compounds of formazan. The dye formed is water soluble and the dye intensity can be read at a given wavelength with a spectrophotometer. The intensity of the dye is proportional to the number of metabolic active cells. The use of multiwell plates and an ELISA reader enables testing a large number of samples and obtaining easy and rapid results. The test procedure includes cultivation of cells in a 96-well plate, adding the XTT reagent and incubation for 2-24 hours. During incubation an orange color is formed, the intensity of which can be measured with a spectrophotometer, in this instance with an ELISA reader. The greater the number of active cells in the well, the greater the activity of mitochondria enzymes, and the higher the concentration of the dye formed, which can then be measured and quantitated.

### **4. Procedure**

- 4.1 The cells should be cultivated in a flat 96-well plate. To each well add 100 $\mu$ l of growth media. The cells should be incubated in a CO<sub>2</sub> incubator at 37°C. In most cases cells can be used to assay proliferation after 24-96 hours.
- 4.2 Defrost the XTT reagent solution and the activation solution immediately prior to use in a 37°C bath. Swirl gently until clear solutions are obtained.
- 4.3 To prepare a reaction solution sufficient for one plate (96 wells), add 0.1ml activation solution to 5ml XTT reagent.
- 4.4 Add 50 $\mu$ l of the reaction solution to each well and incubate the plate in an incubator for 2-24 hours (usually, 2-5 hours are sufficient).
- 4.5 Shake the plate gently to evenly distribute the dye in the wells.
- 4.6 Measure the absorbance of the samples with a spectrophotometer (ELISA reader) at a wavelength of 450-500 nanometer. In order to measure reference absorbance (to measure non-specific readings), use a wavelength of 630-690 nanometer.

## MD Biosciences

Division of Morwell Diagnostics GmbH  
Gewerbestr. 9 8132 Egg b. Zürich Switzerland  
Tel: +41-1 986 26 28 Fax: +41-1 986 26 30

For USA/ Canada : (651) 641 1770 Fax: (651) 641 1773

For Israel Toll Free: 1-800-200 MDB Fax: (8) 940 9973

For United Kingdom (UK) : 0800 0197 200

e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

### 5. Notes

- 5.1 Defrost and prepare the reaction mixture only immediately prior to use.
- 5.2 Since the test is extremely sensitive, it is possible to use a low concentration of cells in the wells (approximately 5000 cells per well). Since there are cell types which show low metabolic activity, such as lymphocytes, kartinocytes and melanocytes, it is recommended to increase the concentration of cells to  $2.5 \times 10^5$  cells per well, in order to obtain development of formazan color within a reasonable period of time.
- 5.3 Incubation time with the reaction mixture varies according to the type and concentration of the cells. Therefore, it is advisable to perform an initial test by reading the absorbance at various time lapses, i.e. after 4, 6, 8, 12 hours using the same plate.
- 5.4 Prior to reading the absorbance with a spectrophotometer, the plate should be gently shaken in order to evenly distribute the dye in the wells.
- 5.5 If the volume of the media in each well is larger than 100 $\mu$ l, add a larger amount of reaction mixture by the same increment (i.e. 100 $\mu$ l reaction mixture to 200 $\mu$ l growth media).

### 6. Summary of Proliferation Assay Using XTT Reagent: Flow Chart

- Defrost XTT reagent and activation reagent (37°C)
- Prepare reaction mixture (0.1ml activation reagent and 5ml XTT reagent for one plate)
- Add 50 $\mu$ l reaction mixture for each well (containing 100 $\mu$ l media)
- Incubate at 37°C for 2-24 hours (in most cases incubation for 2-5 hours is sufficient)
- Measure absorbance at a wavelength of 450-500 nanometer (reference absorbance: at a wavelength of 630-690 nanometer)

### 7. References

- \* Hansen, M.B., et al, (1989), *J. Immunol. Meth.* 119, 203-210
- \* Jost, L.M., et al, (1992), *J. Immunol. Meth.* 147, 153-165
- \* Roehm, N.W., et al, (1991), *J. Immunol. Meth.* 142, 257-265
- \* Scudiero, P.A., et al, (1988), *Cancer Res.* 48, 4827-4833
- \* Tada, H., et al, (1986), *J. Immunol. Meth.* 93, 157-165
- \* Weislow, O.S., et al, (1989), *J. Natl Cancer Inst.* 81, 577-586

**MD Biosciences**  
Division of Morwell Diagnostics GmbH  
Gewerbestr. 9 8132 Egg b. Zürich Switzerland  
Tel: +41-1 986 26 28 ♦ Fax: +41-1 986 26 30  
For USA/ Canada : (651) 641 1770 ♦ Fax: (651) 641 1773  
For Israel Toll Free: 1-800-200 MDB ♦ Fax: (8) 940 9973  
For United Kingdom (UK) : 0800 0197 200  
e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

## **Material Safety Data Sheet (MSDS)**

### **CELL PROFILARATION KIT**

### **(XTT BASED)**

**Cat. No.:** 409005

**Section 1. ----- Chemical identification -----**

Name: Cell Proliferation Kit (XTT based)

**Section 2. ----- Composition/information on ingredients -----**

Cas #: 111072-31-2  
Mf: c22h16n7nao13s2

**Section 3. ----- Hazards identification -----**

Label precautionary statements.  
Caution: Avoid contact and inhalation.

**Section 4. ----- First-aid measures -----**

If swallowed, wash out mouth with water provided person is conscious. Call a physician. If inhaled, remove to fresh air. If not breathing give artificial respiration. If breathing is difficult, give oxygen. In case of contact, immediately wash skin with soap and copious amounts of water. In case of contact, immediately flush eyes with copious amounts of water for at least 15 minutes.

**Section 5. ----- Fire fighting measures -----**

Extinguishing media with water spray, carbon dioxide, dry chemical powder or appropriate foam.  
Special firefighting procedures: wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.  
Unusual fire and explosions hazards: emits toxic fumes under fire conditions.

## **MD Biosciences**

Division of Morwell Diagnostics GmbH

Gewerbestr. 9 8132 Egg b. Zürich Switzerland

Tel: +41-1 986 26 28 ♦ Fax: +41-1 986 26 30

For USA/ Canada : (651) 641 1770 ♦ Fax: (651) 641 1773

For Israel Toll Free: 1-800-200 MDB ♦ Fax: (8) 940 9973

For United Kingdom (UK) : 0800 0197 200

e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

### **Section 6. - - - - - Accidental release measures- - - - -**

Wear respirator, chemical safety goggles, rubber boots and heavy rubber gloves. Sweep up, place in a bag and hold for waste disposal. Avoid raising dust. Ventilate area and wash spill site after material pickup is complete. Evacuate area.

### **Section 7. - - - - - Handling and storage- - - - -**

Refer to section 8.

### **Section 8. - - - - - Exposure controls/personal protection- - - - -**

Niosh/msha-approved respirator.  
Compatible chemical-resistant gloves.  
Chemical safety goggles.  
Safety shower and eye bath.  
Mechanical exhaust required.  
Use only in a chemical fume hood. Avoid inhalation.  
Avoid contact with eyes, skin and clothing.  
Avoid prolonged or repeated exposure.  
Wash thoroughly after handling.  
Wash contaminated clothing before re-use.  
Keep tightly closed.  
Refrigerate.  
Store at 2-8°C.

### **Section 9. - - - - - Physical and chemical properties - - - - -**

Appearance and odor: liquid  
Physical properties  
Solubility: water -z1076

### **Section 10. - - - - - Stability and reactivity - - - - -**

*Stability:* stable under ordinary conditions of use and storage.

*Incompatibilities:*

Strong oxidizing agents  
Strong acids  
Strong bases

*Hazardous combustion or decomposition products:*

Carbon monoxide, carbon dioxide  
Nitrogen oxides

**MD Biosciences**

Division of Morwell Diagnostics GmbH  
Gewerbestr. 9 8132 Egg b. Zürich Switzerland  
Tel: +41-1 986 26 28 ♦ Fax: +41-1 986 26 30

For USA/ Canada : (651) 641 1770 ♦ Fax: (651) 641 1773

For Israel Toll Free: 1-800-200 MDB ♦ Fax: (8) 940 9973

For United Kingdom (UK) : 0800 0197 200

e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

Sulfur oxides

Sodium/sodium oxides

*Hazardous polymerization:* not expected to occur.

**Section 11. ----- Toxicological information -----**

*Acute effects:*

May cause skin irritation.

May be harmful if absorbed through the skin.

May cause eye irritation.

Material may be irritating to mucous membranes and upper respiratory tract.

May be harmful if swallowed.

May be harmful by inhalation, ingestion, or skin absorption.

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

**Section 12. ----- Ecological information -----**

Data not yet available.

**Section 13. ----- Disposal considerations -----**

Cleanup and disposal in accordance to local environmental regulations.

**Section 14. ----- Transport information -----**

Contact Biological Industries for transportation information.

**Section 15. ----- Regulatory information -----**

*European information:*

Caution: substance not yet fully tested.

**Section 16. ----- Other information -----**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Biological industries shall not be held liable for any damage resulting from handling or from contact with