

**Moss, Inc.**  
P.O. Box 189  
Pasadena, MD 21123  
1-800-932-6677  
410-768-3442  
FAX 410-768-3971

### CERTIFICATE OF ANALYSIS

PRODUCT: PE-IgG conjugate, goat anti-human IgA, alpha-specific  
CODE: GTIA-001  
LOT NUMBER: 02038061  
CONCENTRATION: 1 mg/mL  
EXPIRATION: AUG 2007  
ORIGIN: Goat anti-human IgA alpha, affinity purified, U.S.A. origin  
R-phycoerythrin, U.S.A. origin  
BUFFER: 10 mM sodium phosphate, 140 mM sodium chloride, pH 7.3,  
stabilizer 15 mg/mL BSA  
PRESERVATIVE: 0.05% (w/v) sodium azide.  
STORAGE: 2-8°C  
APPEARANCE: Pass  
UV/VIS: Pass  
FLUORESCENCE: Pass  
SHIPPING: Ice Packs

Alvydas Ozinskas  
Manager QC/QA

# **FLUORESCENT CONJUGATE SOLUTION PE-IgG, GOAT ANTI-HUMAN IgA, ALPHA-SPECIFIC MOSS INC. PRODUCT NO. GTIA-001**

## **INTRODUCTION**

Moss Inc. Phycoerythrin-IgG (PE-IgG) conjugates excel in diagnostic, molecular and cellular fluorescence assays based on antibody detection. Moss PE-IgG conjugates are manufactured reproducibly in homogeneous, liquid stable form and are suitable for use on various immunoassay, flow cytometry, and multiplexing platforms such as the Luminex and microarrays. Moss PE-IgG conjugation technology produces conjugates that result in exceptional signal-to-noise ratios, high titers, and the conjugates can also be customized to maximize performance for specific platform applications.

The goat IgG, anti-human IgA, alpha-specific is affinity purified and binds the heavy chain or alpha-chain of human IgA. R-Phycoerythrin is a pink-colored protein purified from seaweed. The intense R-PE absorption maximum at 566 nm ( $E_{566\text{ nm}} = 1.96 \times 10^6 \text{ M}^{-1} \text{ cm}^{-1}$ ) and the strong relative maxima near 545 nm and 498 nm provide multiple opportunities to select excitation wavelengths. The emission maximum is at 578 nm with a high quantum yield. Product No. GTIA-001, available from **Moss Inc.**, is a highly sensitive single component reagent that is ready to use for the quantitative detection of human IgA bound to a solid phase such as a microsphere or microarray, a biological cell, or in free solution. **GTIA-001 is stable for 28 months when refrigerated** and is not appreciably sensitive to normal laboratory light over the course of a typical usage and detection cycle. It should be refrigerated dark when not in use, and exposure to sunlight should be avoided.

## **USAGE SYNOPSIS**

Upon completion of a procedure requiring the detection of human IgA, the appropriate dilution of GTIA-001 is added. Dilutions of PE-anti-hlgA can be made using Moss Inc. PE diluent, Product No. PECD-100, -500, and -1000. If the required dilution is unknown, dilutions of GTIA-001 can be made beginning with a 1:50 dilution. Appropriate fluorescence excitation and emission wavelengths and associated parameters can be selected according to instrument capabilities. Excitation wavelengths selected near the absorption maxima of 566 nm, 545 nm, and 498 nm result in higher fluorescence emission intensities.

## **REAGENT PROVIDED**

**GTIA-001 SOLUTION:** Contains 1 mg/mL PE-IgG, goat anti-human IgA, alpha-specific in buffer containing 10 mM sodium phosphate, 140 mM sodium chloride, pH 7.3. The solution is stabilized with 15 mg/mL BSA and preserved with 0.05% sodium azide.

Store dark at refrigerator temperatures, 2-8°C.

**Protect from exposure to sunlight.**

## **AVAILABLE PE-Anti-hlgA DILUENT**

**PECD-100 SOLUTION:** Contains a proprietary phosphate buffered saline solution, pH 7.3 with 0.05% sodium azide preservative.

Store at refrigerator temperatures, 2-8°C.



## QC USING FRET METHOD

The Moss Inc. FRET test was developed for internal use to test PE-IgG conjugates, and of course any suitable methodology may be used to test PE-IgG conjugates. PE-anti-hlgA serves as the energy transfer donor in the test. The energy transfer acceptor conjugate is Cy5-human IgA (Cy5-hlgA) and is prepared by Moss Inc. The FRET test is a homogeneous assay, which means that there are no separation or wash steps. The reagents are mixed together in microtiter plate wells, and then the fluorescence is measured after 15 minutes.

### FRET TEST PROCEDURE USING 96 WELL BLACK MICROTITER PLATES

The assay is conducted in a black microtiter plate. The plate is divided into an upper half (Rows A-D) and a lower half (Rows E-H). Each column (1-12) has n=2 donor alone replicates (each well in rows A and B has PE-anti-hlgA and buffer, 100 µL final volume). Each column also has n=2 donor plus acceptor replicates (each well in Rows C and D has PE-anti-hlgA and Cy5-hlgA, 100 µL final volume). One or more columns can be used for the buffer alone control (no donor or acceptor) and for the acceptor alone control (no donor) with 100 µL final volume for each well.

1. Prepare 0.01 mg/mL dilutions each of PE-anti-hlgA and Cy5-hlgA using a buffer of 10 mM PBS, pH 7.3, 1.5 mg/mL BSA. It is convenient to use a second microtiter plate to prepare the PE-anti-hlgA dilutions (300 µL each) for use with a multi-channel pipettor.
2. Add 50 µL of buffer to Rows A, B E, and F.
3. For the upper half of the plate, add 50 µL of one set of PE-anti-hlgA donor conjugates to each well in rows A, B, C and D. A different PE-anti-hlgA donor conjugate can be used for each column, and in this case n=2. We usually use the same PE-anti-hlgA conjugate for 4 columns, so n=8 in this case.
4. To the lower half of the plate, add 50 µL of the next set of PE-anti-hlgA donor conjugates to each well in rows E, F, G and H.
5. Add 50 µL of Cy5-hlgA acceptor to each well in Rows C, D, G and H.
6. Mix plate on orbital shaker at 700 RPM for 5-7 seconds.
7. Read plate after 15 minutes at ambient temperature using a fluorescence microtiter plate reader.

### **Data Processing:**

1. Average the 15 minute readings for each PE-anti-hlgA conjugate, separately for the donor alone (Fo), and for the donor plus acceptor (F). The CVs are generally 1-3%.
2. Correct for background by subtracting the fluorescence readings for the buffer alone control from Fo, and by subtracting the fluorescence readings for the buffer plus acceptor control from F.
3. Calculate F/Fo for each PE-anti-hlgA conjugate, using F and Fo values that were corrected for background. Normalize the data if required. Prepare a table, graph, or chart.

### **Significance of results:**

1. No binding results in no energy transfer and F/Fo = 1.00, theoretically.
2. A binding interaction between donor and acceptor results in fluorescence resonance energy transfer and F/Fo < 1.00.
3. The smaller the F/Fo value, the stronger the binding interaction.
4. Variations of the assay are possible by adding binding ligands to generate an analyte titration curve.
5. In general, for a given PE-anti-hlgA conjugate, the range between replicates for F/Fo = 0.03.

The method described is comparable to that used at **Moss Inc.** Variations in test parameters require standardization by the user.

MOSS DOCUMENT NO. GTIA-001, MARCH 2007. FOR TECHNICAL ASSISTANCE OR FOR INFORMATION ABOUT MOSS, INC. PRODUCTS,  
PLEASE VISIT OUR WEBSITE AT [www.MossSubstrates.com](http://www.MossSubstrates.com) OR CALL 800-932-6677. OUTSIDE THE USA CALL 410-768-3442 OR FAX 410-768-3971.