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[Reagent for food testing]

FASTKIT SLIM SOYBEANS <<Instruction Manual>>

*Please read this manual before using this product.

Cat# NPH-NFS006

[Introduction]

The Food Sanitation Law obliges that 7 specified allergenic ingredients (egg, milk, wheat, buckwheat, peanut, shrimp and crab), which have a high risk of inducing food allergy, should be listed on food labels. The Food Sanitation Law also recommends that 18 ingredients, including soybeans, should be listed on the label in a way similar to the specified allergenic ingredients.

“Providing accurate information by listing of ingredients on food labels” and “preventing contamination with the ingredients which are not listed on the labels” are required to prevent food allergy. Moreover, as a measure to prevent contamination, routine management at the manufacturing sites, including verification of the process of washing the machines and devices used for production, is important.

This product is a kit for detecting soybean protein by immunochromatography. Tests can be conducted rapidly and simply by means of the kit.

[Product Features]

- 1) The simple one-step operation enables easy assessment, and anybody can perform the test easily.
- 2) Because there is no need for special test equipment and the test gives rapid results (in 15 minutes), it is best for routine management at the production sites.

[Kit contents]

1) Components

- A: Test strip..... 2-test × 10 packs
B: Dilution Buffer..... 50mL×1 bottle
C: Extraction buffer (1/10 concentration).... 100mL×1 bottle
D: Instruction manual..... 1 sheet
E: Plastic pouched bag..... 1 bag

2) Ingredients


- (1) Reagent-containing section
Gold colloid-labeled anti-soybean protein antibody (rabbits)
(2) Detecting section
Anti-soybean protein antibody (rabbits)
Anti-rabbit immunoglobulin antibody (goat)

[Application]

- 1) Detection of soybean protein in foods or solutions


[Illustration of Test strip and the Principle of assay]

1) Illustration of Test strip

- a. Sample solution drop section (*Be careful not to touch this section with your finger.*)
- b. Reagent-containing section
- c. Detecting section (*Be careful not to scratch this section and touch this section with your finger.*)
- d. Absorbent pad
- e. Measurement items listing position
- f. Test line appearance position (*Approx. 30mm from the sample solution drop section.*)
- g. Control line appearance position (*Approx. 38mm from the sample solution drop section.*)

2) Principle of assay

When a sample solution is dropped onto the sample solution drop section of the test strip, the gold colloid-labeled anti-soybean protein antibody (2) in the reagent-containing section dissolves and forms complexes with soybean protein (1). These complexes move to the detecting section by capillary attraction and are trapped by the anti-soybean protein antibody (3) that is fixed in the test line appearance position. This results in the appearance of a reddish purple line of gold colloid. This reddish purple line can be detected by visual inspection and used to judge the presence or absence of soybean

protein in the sample solution.

The excess gold-labeled antibodies, regardless of the presence or absence of soybean protein in the sample solution, travel further through the detecting section and are trapped by the anti-goat immunoglobulin rabbit antibody (4) fixed at the control line appearance position, where they form a second reddish purple line. The presence of this line indicates that the sample solution has reached the detecting section.

[Preparation of the sample solution 1 (detection in food)]

1) Required Equipment and Instruments

Grinder (food cutter), homogenizer (Millser), centrifuge (3,000×g or more, 4°C is recommended), sedimentation tubes, etc., filter paper, funnel, measuring cylinder, beaker, micropipette, etc.

2) Preparation of Reagent

- Dilution buffer: Leave it at room temperature (20°C-25°C) before use.
- Extraction buffer (1/10 concentration): Dilute the solution 1:10 with purified water before use.

Note 1) If any sediments are observed in the extraction buffer (1/10 concentration), leave the solution at room temperature until the sediments dissolve, and then, dilute it with purified water.

Note 2) Since the dilution buffer and the extraction buffer (1/10 concentration) of the FASTKIT Slim Series have the same composition, they can be used for any type of tests.

3) Extracting operation (example of operation in ordinary foods)

- (1) Grind the food to be tested (sample) of each package unit to be a more homogenous state using a grinder or food cutter (Note 1).
- (2) To 2 g of homogenized sample, add 38 mL of the extraction buffer, which is dispensed in advance, and repeat the extracting operation 30 to 60 seconds × 3 times using a homogenizer, etc. (Note 2).
- (3) Centrifuge the above sample at 3,000×g or more, 4°C for 20 minutes, and filter the supernatant (Note 3).
- (4) Dilute the filtrate with the dilution buffer (1:10) and use it as the sample solution (Note 4).

Note 1) To prevent contamination via equipments, thoroughly wash the equipments before use. In particular, grinders and homogenizers should be washed completely after processing of each sample. (Wash them with neutral detergent and soak them in alkaline detergent over night, or perform ultrasonic cleaning with alkaline detergent.)

Note 2) Determine the pH during extraction. If necessary, adjust the pH to the neutral (pH 6.0 to 8.0) if necessary.

Note 3) If the sample solution contains insoluble matters, such as lipids, in a large amount, accurate results may not be obtained. Insoluble matters should be eliminated as much as possible.

Note 4) Store the sample solution at 4°C. The shelf life of the sample solutions varies depending on the foods. Perform the test as soon as possible.

[Preparation of the sample solution 2 (Smear test, etc.)]

1) Required Equipment and Instruments

Test tube, etc., swab, micropipettes, etc.

2) Preparation of reagents and devices (Note 1)

- Smearing solution and swab: Prepare 0.9% sodium chloride solution (physiological saline) as the smearing solution, and dispense a certain portion into a test tube, etc. Immerse a swab in the dispensed smearing solution (Note 2, 3).

Note 1) A commercially available smearing kit for microbiological tests can be used. However, please beware of the composition of the solution contained.

Note 2) Since there is a risk of contamination of the manufacturing machine and devices with the smearing solution, do not use the dilution buffer or the extraction buffer (1/10 concentration) provided in the kit.

Note 3) The sensitivity of the smear test may vary depending on the amount of the dispensed smearing solution. The minimum volume of the dispensed smearing solution should be 1 mL. Determine the optimal volume by yourself.

3) Smearing operation (example of ordinary operation)

- (1) Specify the smearing sites of the target machines and tools (Note 1).
- (2) Rub the designated smearing sites in (1) with a swab which has been moisturized in advance (Note 2).
- (3) Thoroughly wash the swab with the smearing solution to suspend the soybean protein adhered to the swab (Note 3).
- (4) Use the suspended smearing solution as the sample solution (Note 4).

Note 1) It is recommended that the smearing sites should be the sites where stains are likely to remain, or the sites which are difficult to wash.

Note 2) It is recommended that the entire surface of the smearing site should be rubbed.

Note 3) After you have washed the swab, you should thoroughly squeeze out the smearing solution absorbed by the swab, using the wall of the test tube, etc.


Note 4) If the sample solution contains insoluble matters in a large amount, accurate results may not be obtained. Insoluble matters should be eliminated as much as possible by centrifugation or filtration.

Note 5) When physiological saline is used for dispensing the sample solution, a very subtle line may be observed as time goes by. Make assessment 15 minutes after the start of the test.

[Operating Procedures for Testing]

- (1) Leave the test strip in the aluminum pouch at room temperature and remove from the pouch immediately before use (Notes 1, 2).
- (2) With an oil-based marker pen, write the name of the test sample or the number of the subject under test on the absorbent pad of the test strip removed from the pouch.
- (3) Place the test strip carefully on the flat stand and drop a 100 μ L-portion of the sample solution onto the sample solution drop section (see the figure on the left). Otherwise, dispense a 150 μ L portion of the sample solution into a test tube and attach the test strip to the test tube so that the sample solution drop section of the test strip is immersed in the sample solution (see the figure on the right)(Notes 3, 4).

Note 1) Do not remove the test strip from the aluminum pouch until it has returned to room temperature, otherwise incorrect test results may be obtained as a result of moisture absorption. Test strips which are not used should be placed in the plastic pouch bag again together with a desiccant, and should be preserved in a refrigerator.


Note 2) Be careful not to scratch the sample solution drop section or detecting section and do not touch them with your fingers. When handling the test strip, make sure that you hold the absorbent pad.

Note 3) Make sure that you change the pipette or chip to drop or dispense the sample solution for every sample solution.

Note 4) Make sure that the 100- μ L portion of sample solution does not overflow the test strip when dropping it. If necessary, drop the solution in two or more portions

[Judgment of the Test Results]

- (1) The test results is judged as positive when a reddish purple line is observed at the test line appearance position and at the control line appearance position 15 minutes after the start of the test.
- (2) Judge the test results as negative when no reddish purple line is observed at the test line appearance position, but a line is observed at the control line appearance position.

Note 1) Be sure to use the results obtained 15 minutes after the start of the test as the test results. Deviation from the 15-minute results may occur as time goes by.

Note 2) Judge the test results as positive if a reddish purple line is observed, regardless of the shading of the reddish purple line.

Note 3) Retest in cases where no reddish purple line is observed at the control appearance position. It is likely that there is something abnormal in the development of the sample solution in such cases.

Note 4) This kit is a reagent designed to identify soybean protein qualitatively. It cannot determine the content of soybean protein. Use FASTKIT ELISA Ver.II Soybeans if you need to determine the content of soybean protein.

[Performance]

1) Sensitivity Test

The results of tests conducted in accordance with instructions for “the preparation of the sample solution” and “the operating procedures” described in this instruction manual will be positive when the concentration of soybean protein in the sample solution is 25 ng/mL or more.

2) Repeatability Tests

When positive and negative sample solutions of soybean protein were simultaneously tested three times each, all positive sample solutions exhibited positive results and all negative sample solutions show negative results.

3) False Positive / False Negative Results

- (1) Cross-reactivity with the 7 specified ingredients (egg, milk, wheat, buckwheat, peanut, shrimp and crab) has not been observed.
- (2) It may exhibit false-positive responses to some types of beans (ex. pinto beans, tiger beans, large kintoki beans, red kidney beans, green kidney beans, purple speckled kidney beans, garbanzo beans)
- (3) False positive results may be obtained due to non-specific reactions in the presence of viscous food or very high of concentrations of

protein. In such cases, dilute the sample to an appropriate concentration before conducting the test.

- (4) Pressurized and heated food packaged in containers, canned food, and food in retort pouch packaging may exhibit false negative results.
- (5) Soy sauce, soy oil, natto, soybean peptide, soybean isoflavone, and soybean lecithin may exhibit false-negative results because the soybean protein content is small, or because soybean protein is degraded during the process of production.
- (6) When the sample solution contains a high concentration of soybean protein, an obscure line or no line may be observed in the smear test. In such cases, dilute the sample to an appropriate concentration before conducting the test.
- (7) See the “List of Foods Exhibiting False Positive and False Negative Results” on the website of Nippon Meat Packers Inc. R&D center for the foods exhibiting false positive or false negative results.

[Precautions in Using the Kit]

1) Precautions in Handling the Kit

- (1) Read the instruction manual carefully before use. Use the kit in accordance with the test method described in this manual.
- (2) Do not use a kit whose use-by date has passed. The expiry date is indicated on the label on the external package of the kit and on the aluminum pouch of the test strip.
- (3) The kit is a reagent designed to detect soybean protein in foods or solutions. It is not a reagent for diagnosing the presence or absence of food allergy to soybeans. Correlation between the results obtained from the test with this kit and the occurrence of allergic symptoms has not been confirmed.
- (4) Globally judge the presence or absence of soybean protein not only on the basis of the results obtained with this kit, but also using other methods, such as confirmation of raw materials and ingredients and production records.
- (5) Ask their manufacturers or distributors about how to use the equipments and instruments to be used for the test with this kit.
- (6) This instruction manual is intended as a guideline for those in charge of testing. Verify your own operating procedures for the kit and the appropriateness of its use for each particular food.
- (7) Product specifications may be changed without notice.

2) Precautions Regarding Risk Prevention

- (1) Beware so that the reagents of this kit or sample solutions do not adhere to skin, mucosa or clothes.
- (2) If you accidentally get any reagent or sample solution in your eyes or mouth, take emergency measures, such as immediately washing away the solution with tap water, and then seek medical attention, if necessary.

3) Precautions Regarding Disposal of Waste Materials

- (1) Discard the kits, test samples, and surplus sample solutions in strict compliance with your local waste-disposal regulations and with full consideration of environmental sanitation.

[Storage Method and Use-By Date]

- (1) Storage method: Refrigerate at 2–8 °C and shade from the light. Avoid freezing.
- (2) Use-by date: 12 months from the date of manufacture. Stated on the external package and on each component.

[References]

- 1) Ministry of Health, Labour and Welfare: Labels for foods containing allergenic substances (*Shokukihatsu* No. 2 *Shokukanhatsu* No.46 dated March 21, 2001, final revision, *Shokuankihatsu* No. 012201 *Shokukanhatsu* No. 0122002 dated January 22, 2009)
- 2) Ministry of Health, Labour and Welfare: Enforcement of the ministerial ordinance partially amending the Ordinance for Enforcement of the Food Sanitation Act (*Shokuanhatsu* No. 0603001 dated June 3, 2008)
- 3) Ministry of Health, Labour and Welfare: Interim report (outline) by the Allergy Labeling Commission, Food Labeling Study Team (October 29, 2001)
- 4) Ministry of Health, Labour and Welfare: Methods of the test of foods containing allergenic substances (final revision, *Shokuanhatsu* No. 0122001 dated January 22, 2009)

Manufactured by Nippon Meat Packers, Inc. R&D Center

COSMO BIO CO., LTD.

[JAPAN]

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

[Outside Japan]

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com