


Bradford ULTRA

Detergent Compatible

Coomassie-based protein quantitation method



Recommended Protocol and Notes on Usage

Before Use

Mix the Bradford ULTRA Reagent solution immediately before use by gently inverting the bottle several times (Do not shake the bottle to mix the solution). Remove the amount of reagent needed and equilibrate it to room temperature before use. It is good practice to mix the Bradford ULTRA Reagent before dispensing and to mix each tube or plate immediately before measuring absorbances.

Protein Standards

The intensity of the colour formation is protein dependent and is correlated to the basicity of the protein, in particular the number of positive charges on the protein (lysine, arginine and histidine), and hydrophobicity of the protein. A suitable standard will have a similar mol% positively charged residues as the target protein. For most proteins the content of basic amino acids ranges from 10 to 17 mol%¹. BSA is the most commonly used reference standard in Bradford methods. Some useful standards in this range are given in the table below.

	mol% positive residues	Supplier
Hen egg white lysozyme	13.9	Sigma (# L3790)
Bovine serum albumin (BSA)	16.5	Sigma (# A7638)
Bovine Immunoglobulin (IgG)	11.3	Sigma (# I5506)
Bovine beta lactoglobulin	11.8	Sigma (# L0130)

To calculate the mol% basic residues of your target protein visit:

<http://us.expasy.org/tools/protparam.html>

¹ Tal, M., et al., *J. Biol. Chem.* **260**, 9976, (1985)

Introduction

The Bradford ULTRA kit is a quick and ready-to-use Coomassie-binding, colorimetric method for total protein quantitation in an environment containing up to 1% detergent (1% high protein range, 0.1% low protein range).

Novexin's Bradford ULTRA is an improvement over classical Bradford formulations that cannot tolerate detergent contamination of the protein samples. In addition, the Bradford ULTRA reagent shows excellent linearity for a defined range of protein concentrations and shows significantly less protein-to-protein variation than is observed with other Bradford-type Coomassie formulations.

When Coomassie dye binds protein in an acidic medium, an immediate shift in absorption maximum occurs from 465 nm to 595 nm with a concomitant color change from brown to blue. Protein concentrations are estimated by reference to absorbances obtained for a series of standard protein dilutions, which are assayed alongside the samples with unknown protein concentrations.

Contents

500 ml reagent, containing Coomassie dyes, ethanol, phosphoric acid and solubilizing agents in water. (**Caution:** Phosphoric acid is a corrosive liquid.)

Storage

Upon receipt store at + 4°C. Discard any reagents that show discoloration or evidence of microbial contamination.

Protocol

- 1) Make a dilution series of the chosen model protein in the range:
0.1 mg/ml – 1.5 mg/ml (high protein range) OR
1 µg/ml – 25 µg/ml (low protein range)
- 2) Mix the samples, standards and a blank (buffer, no protein) with Bradford ULTRA reagent.

For 0.1 mg/ml – 1.5 mg/ml protein (high range)

Sample / Reagent ratio: 1 / 15

For microtiter plate : 20 µl sample + 300 µl Bradford ULTRA

For cuvette : 100 µl sample + 1.5 ml Bradford ULTRA

For 1 µg/ml – 25 µg/ml protein (low range)

Sample / Reagent ratio: 1 / 1

For microtiter plate : 150 µl sample + 150 µl Bradford ULTRA


For cuvette : 750 µl sample + 750 µl Bradford ULTRA

NOTE:

Preferably all samples, standards and blanks are prepared in triplicates.



- 3) Read absorbance at 595 nm.
- 4) Subtract the average 595 nm measurement for the blank from the 595 nm measurements of all other individual standards and unknown samples. Plot the average blank-corrected 595 nm measurement for each standard vs. concentration. Use the slope of this standard curve to estimate the protein concentration of the unknown samples.

Detergent Compatibility



Comparison of Novexin's "Bradford ULTRA" Assay with Pierce's "Coomassie Protein Assay". The graph shows the average blank corrected A595 measurements for detergents.

Novexin's Bradford ULTRA assay is not affected by the detergents. Classic formulations are subject to the high background noise at 595nm in the presence of detergents.



Standard curves obtained with Bradford ULTRA are unaffected by the presence of detergents. Standard curves obtained with classical Bradford formulation are significantly affected by the presence of detergents resulting in loss of sensitivity and inaccurate results.

Contact

For technical support email:

scientist@novexin.com

Novexin Ltd
Babraham Bioincubator
Babraham
Cambridge CB2 4AT

Tel: +44 (0)1223 496500
Fax: +44 (0)1223 496741