TRAP Staining Kit

Catalog No. AK04F

July 1, 2024

Introduction

The TRAP Staining Kit (Cat.No.AK04F) is used for the staining of Tartrate-Resistant Acid Phosphatase in osteoclasts. Bone mass is controlled by the balance between the activity of osteoblasts bone formation and the activity of osteoclasts. Alkaline phosphatase (ALP: ALP staining kit, Cat.No.AK20) and tartrate-resistant acid phosphate are used as markers for osteoblasts and osteoclasts, respectively.

Components

Component	Quantity	Storage
Tartrate-containing Buffer	50mL	4°C (Do Not Freeze)
Chromogenic Substrate	10vials	4°C

One kit contains reagents for 10 × 96-well plates

Additional Materials Required

- 10% Formalin Neutral Buffer Solution
- Phosphate Buffered Saline (PBS)
- Distilled or deionized water

Protocol (96-well plate format)

1. Remove culture medium. Wash each well once with 100uL of PBS.
2. Add 50uL of the 10% Formalin Neutral Buffer Solution to each well and fix for 5 minutes at room temperature.
3. Wash each well with 250uL of deionized water. (× 3 times)
4. Dissolve 1 vial of Chromogenic Substrate with 5mL of Tartrate-containing Buffer.
5. Add 50uL of Chromogenic Substrates to each well.
6. Incubate at 37°C for 20-60 minutes. Adjust incubation time until stained TRAP is clearly showing the result in figure 1.
7. Wash with deionized water to stop the reaction.

Note: Excess incubation will be cause of over stained.


Figure 1: TRAP staining of Osteoclast

References

- (1) Sunao Takesita, Keisuke Kaji, Akira Kudo. Identification and Characterization of the New Osteoclast Progenitor with Macrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. *JOURNAL OF BONE AND MINERAL RESEARCH* Volume 15, Number 8 (2000). 1477-1488.
- (2) BEN A. A. SCHEVEN, JONE S. MILNE, SIMON P. ROBINS. A SEQUENTIAL CULTURE APPROACH TO STUDY OSTEOCLAST DIFFERENTIATION FROM NONADHERENT PORCINE BONE MARROW CELLS. *In Vitro Cell. Dev. Biol.* July-August (1998). *Animal* 34: 568-577.
- (3) Martha J. Somerman, Janice E. Berry, Zhila Khalkhali-Ellis, Philip Osdoby, Robert U. Simpson. Enhanced Expression of α_v Integrin Subunit and Osteopontin during Differentiation of HL-60 Cells along Monocytic Pathway. *EXPERIMENTAL CELL RESEARCH* 216, 335-341(1995).
- (4) Ichiro Itonaga, Afsie Sabokbar, Susan D. Neale, Nicholas A. Athanasou. 1,25-Dihydroxyvitamin D₃ and Prostaglandin E₂ Act Directly on Circulating Human Osteoclast Precursors. *Biochemical and Biophysical Research Communications* 264, 590-595(1999).
- (5) Hiroshi Takayanagi, Kouetsu Ogasawara, Shigeaki Hida, Tomoki Chiba, Shigeo Murata, Kojiro Sato, Akinori Takaoka, Taeko Yokochi, Hiromi Oda, Keiji Tanaka, Kozo Nakamura and Tadatsugu Taniguchi. T-cell-mediated regulation of osteoclastogenesis by signaling cross-talk between RANKL and IFN- γ . *Nature*, Vol. 408, No.6812, 600-605, 30 November 2000.
- (6) Morinobu, A., Biao, W., Tanaka, S., Horiuchi, M., Jun, L., Tsuji, G., Sakai, Y., Kurosaka, M., Kumagai, S. (-)-Epigallocatechin-3-Gallate Suppresses Osteoclast Differentiation and Ameliorates Experimental Arthritis in Mice. *Arthritis Rheum.* 58, 2012-2018 (2008)
- (7) Hase, H., Kanno, Y., Kojima, H., Sakurai, D., Kobata, T. Coculture of Osteoclast Precursors With Rheumatoid Synovial Fibroblasts Induces Osteoclastogenesis via Transforming Growth Factor beta-Mediated Down-Regulation of Osteoprotegerin. *Arthritis Rheum.* 58, 3356-3365 (2008)
- (8) Cao, H., Kuboyama, N. A Biodegradable Porous Composite Scaffold of PGA/beta-TCP for Bone Tissue Engineering. *Bone.* 46, 386-395 (2010)
- (9) Okada, Y., Hamada, N., Kim, Y., Takahashi, Y., Sasaguri, K., Ozono, S., Sato, S. Blockade of Sympathetic beta-receptors Inhibits Porphyromonas Gingivalis-induced Alveolar Bone Loss in an Experimental Rat Periodontitis Model. *Arch. Oral Biol.* 55, 502-508 (2010)
- (10) Okamoto, A., Ohnishi, T., Bandow, K., Kakimoto, K., Chiba, N., Maeda, A., Fukunaga, T., Miyawaki, S., Matsuguchi, T. Reduction of Orthodontic Tooth Movement by Experimentally Induced Periodontal Inflammation in Mice. *Eur. J. Oral Sci.* 117, 238-247 (2009)

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- (1) Sunao Takesita, Keisuke Kaji, Akira Kudo. Identification and Characterization of the New Osteoclast Progenitor with Macrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. JOURNAL OF BONE AND MINERAL RESEARCH Volume 15, Number 8 (2000). 1477-1488.
- (2) BEN A. A. SCHEVEN, JONE S. MILNE, SIMON P. ROBINS. A SEQUENTIAL CULTURE APPROACH TO STUDY OSTEOCLAST DIFFERENTIATION FROM NONADHERENT PORCINE BONE MARROW CELLS. In Vitro Cell. Dev. Biol. July-August (1998). Animal 34: 568-577.
- (3) Martha J. Somerman, Janice E. Berry, Zhila Khalkhali-Ellis, Philip Osdoby, Robert U. Simpson. Enhanced Expression of α_v Integrin Subunit and Osteopontin during Differentiation of HL-60 Cells along Monocytic Pathway. EXPERIMENTAL CELL RESEARCH 216, 335-341(1995).
- (4) Ichiro Itonaga, Afsie Sabokbar, Susan D. Neale, Nicholas A. Athanasou. 1,25-Dihydroxyvitamin D₃ and Prostaglandin E₂ Act Directly on Circulating Human Osteoclast Precursors. Biochemical and Biophysical Research Communications 264, 590-595(1999).
- (5) Hiroshi Takayanagi, Kouetsu Ogasawara, Shigeaki Hida, Tomoki Chiba, Shigeo Murata, Kojiro Sato, Akinori Takaoka, Taeko Yokochi, Hiromi Oda, Keiji Tanaka, Kozo Nakamura and Tadatsugu Taniguchi. T-cell-mediated regulation of osteoclastogenesis by signaling cross-talk between RANKL and IFN- γ . Nature, Vol. 408, No.6812, 600-605, 30 November 2000.
- (6) Morinobu, A., Biao, W., Tanaka, S., Horiuchi, M., Jun, L., Tsuji, G., Sakai, Y., Kurosaka, M., Kumagai, S. (-)-Epigallocatechin-3-Gallate Suppresses Osteoclast Differentiation and Ameliorates Experimental Arthritis in Mice. Arthritis Rheum. 58, 2012-2018 (2008)
- (7) Hase, H., Kanno, Y., Kojima, H., Sakurai, D., Kobata, T. Coculture of Osteoclast Precursors With Rheumatoid Synovial Fibroblasts Induces Osteoclastogenesis via Transforming Growth Factor beta-Mediated Down-Regulation of Osteoprotegerin. Arthritis Rheum. 58, 3356-3365 (2008)
- (8) Cao, H., Kuboyama, N. A Biodegradable Porous Composite Scaffold of PGA/beta-TCP for Bone Tissue Engineering. Bone. 46, 386-395 (2010)
- (9) Okada, Y., Hamada, N., Kim, Y., Takahashi, Y., Sasaguri, K., Ozono, S., Sato, S. Blockade of Sympathetic beta-receptors Inhibits Porphyromonas Gingivalis-induced Alveolar Bone Loss in an Experimental Rat Periodontitis Model. Arch. Oral Biol. 55, 502-508 (2010)
- (10) Okamoto, A., Ohnishi, T., Bandow, K., Kakimoto, K., Chiba, N., Maeda, A., Fukunaga, T., Miyawaki, S., Matsuguchi, T. Reduction of Orthodontic Tooth Movement by Experimentally Induced Periodontal Inflammation in Mice. Eur. J. Oral Sci. 117, 238-247 (2009)