# DNA Quantity Kit

Catalog No. AK06

March 31, 2020

## Introduction

The DNA Quantity Kit (Cat.No.PMC-AK06-COS) is designed to quantify DNA directly from cell culture. DNA can be quantified without purification process from cell culture as Color Development Reagent binds to cellular double stranded DNA and emits blue fluorescence at 458 nm.

## Components

### Component Quantity Storage

Color Development Reagent 10 mL 4°C

Dilution Buffer 150 mL × 2 4°C

DNA Standard (100 g/ ml) 2 mL 4°C

## Additional Materials Required

- Purified water
- Fluorometer
- Phosphate Buffered Saline (PBS)


## Preparation of Standard Curve

### 1. Standard sample preparation

Dilute the stock DNA Standard (100 mg/mL) to 50, 25, 12.5 and 0 mg/mL with purified water.

Diluted sample standards can be stored frozen at -20°C.

2. Transfer each 50 ul of standard sample to new tubes.
3. Add 1ml Dilution Buffer to each standard sample and mix well.
4. Add 50 ul of Color Development Reagent to each tube.
5. Mix thoroughly.
6. Measure fluorescence. (Excitation wavelength at 356nm, emission wavelength at 458nm)
7. Calibration curve (DNA)

## Sample Analysis (24well plate format)

1. Remove culture medium and rinse the culture plate with Phosphate Buffered Saline (PBS).
2. Remove PBS and add 500 ul of dilution buffer to each well.
3. Sonicate cells until completely homogenated.
4. Transfer 50 ul of homogenate sample to new tube.
5. Add 1ml of Dilution Buffer to each tube and mix well.
6. Add 50 ul of Color Developer to each tube.
7. Mix thoroughly.
8. Measure fluorescence. (Excitation wavelength at 356nm, emission wavelength at 458nm)

## References

- (1) Kalgaonkar, S., Nishioka, H., Gross, H. B., Fujii, H. Keen, C. L., Hackman, R. M. Bioactivity of a Flavanol-rich Lychee Fruit Extract in Adipocytes and Its Effects on Oxidant Defense and Indices of Metabolic Syndrome in Animal Models. *Phytother. Res.* 24, 1223-1228 (2010)
- (2) Okura, H., Komoda, H., Saga, A., Kakuta-Yamamoto, A., Hamada, Y., Fumimoto, Y., Lee, CM., Ichinose, A., Sawa, Y., Matsuyama, A. Properties of Hepatocyte-like Cell Clusters from Human Adipose Tissue-Derived Mesenchymal Stem Cells. *Tissue Eng. Part C Methods.* 16, 761-70 (2010)
- (3) Saito, Y., Goto, M., Maya, K., Ogawa, N., Fujimori, K., Kurokawa, Y., Satomi, S. The Influence of Brain Death on Tissue Factor Expression in the Pancreatic Tissues and Isolated Islets in Rats. *Transplant. Proc.* 41, 307-310 (2009)

### Notice to Purchaser

This product is to be used for Research Purposes Only. It is not to be used for Drug or Diagnostic Purposes, nor is it intended for Human Use. Primary cell products may not be resold, modified for resale, or used to manufacture commercial products without the express written consent of Primary Cell Co., Ltd.

Except as otherwise expressly set forth in this use manual, Primary Cell does not make any representation or warranties or conditions of any kind, either express or implied, with respect to the products, or information disclosed hereunder, including, but not limited to, the implied warranties of merchantability, fit for a particular purpose, or noninfringement of the intellectual property rights of third parties.

**COSMO BIO Co., LTD.**

**[JAPAN]**

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: 81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

**COSMO BIO USA**

**[Outside Japan]**

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobioussa.com](mailto:info@cosmobioussa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobioussa.com](http://www.cosmobioussa.com)

**For research use only. Not for clinical diagnosis.**

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COSMO BIO CO., LTD.