DNA Quantity Kit

Catalog No. AK06

March 2, 2026

Introduction

The DNA Quantity Kit (Cat.No.PMC-AK06-COS) is designed to quantify DNA directly from cell culture. DNA can be quantified without purification process from cell culture as Color Development Reagent binds to cellular double stranded DNA and emits blue fluorescence at 458 nm.

Components

Component Quantity Storage

Color Development Reagent 10 mL 4°C

Dilution Buffer 150 mL × 2 4°C

DNA Standard (100 µg/ml) 2 mL 4°C


Additional Materials Required

- Purified water
- Fluorometer
- Phosphate Buffered Saline (PBS)

Preparation of Standard Curve

1. Standard sample preparation
Dilute the stock DNA Standard (100 µg/mL) to 50, 25, 12.5 and 0 µg/mL with purified water.
Diluted sample standards can be stored frozen at -20°C.
2. Transfer each 50 µl of standard sample to new tubes.
3. Add 1ml Dilution Buffer to each standard sample and mix well.
4. Add 50 µl of Color Development Reagent to each tube.
5. Mix thoroughly.
6. Measure fluorescence. (Excitation wavelength at 356nm, emission wavelength at 458nm)

Calibration curve (DNA)

Sample Analysis (24well plate format)

1. Remove culture medium and rinse the culture plate with Phosphate Buffered Saline (PBS).
2. Remove PBS and add 500 µl of dilution buffer to each well.
3. Sonicate cells until completely homogenated.
4. Transfer 50 µl of homogenate sample to new tube.
5. Add 1ml of Dilution Buffer to each tube and mix well.
6. Add 50 µl of Color Developer to each tube.
7. Mix thoroughly.
8. Measure fluorescence. (Excitation wavelength at 356nm, emission wavelength at 458nm)

References

- (1) Kalgaonkar, S., Nishioka, H., Gross, H. B., Fujii, H. Keen, C. L., Hackman, R. M. Bioactivity of a Flavanol-rich Lychee Fruit Extract in Adipocytes and Its Effects on Oxidant Defense and Indices of Metabolic Syndrome in Animal Models. *Phytother. Res.* 24, 1223-1228 (2010)
- (2) Okura, H., Komoda, H., Saga, A., Kakuta-Yamamoto, A., Hamada, Y., Fumimoto, Y., Lee, CM., Ichinose, A., Sawa, Y., Matsuyama, A. Properties of Hepatocyte-like Cell Clusters from Human Adipose Tissue-Derived Mesenchymal Stem Cells. *Tissue Eng. Part C Methods.* 16, 761-70 (2010)
- (3) Saito, Y., Goto, M., Maya, K., Ogawa, N., Fujimori, K., Kurokawa, Y., Satomi, S. The Influence of Brain Death on Tissue Factor Expression in the Pancreatic Tissues and Isolated Islets in Rats. *Transplant. Proc.* 41, 307-310 (2009)

For research use only. Not for clinical diagnosis.

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