# POLARIC™ PLT-500c6

Cat No. AK12

March 31, 2020

POLARIC™ PLT-500c6 induces dramatic changes in the emission spectra depending on hydrophobic variation of the microenvironment of cell organelles.

## I-1. Contents

Name	Unit size	Shipping	Storage	Caution
POLARIC™ PLT-500c6	10 µg/tube 5 tubes	Room temperature	Protect from light >4°C	<i>Always wear gloves and safety glasses when working these materials.</i>

※ For staining 5 x 96-wells plates

※ Only for ethanol resolution

## I-2. Advantages

- Extremely low cytotoxicity
- Minimum cell damage by excitation wavelengths around 460 - 520nm.
- Fade-resistant fluorescence
- Dramatically changes the emission spectra depending on hydrophobicity / hydrophilicity of the microenvironment of cell organelles.

**Fig1. Emission spectra change of solvatochromic fluorophore**


## II. Experimental Protocol

1. After brief centrifugation, resolve the red dye pellet in 3 $\mu$ L of ethanol and transfer dye solution into 10mL of medium (staining solution). Staining solution can be stored at 4°C for 2 weeks. Protect from light.
2. Culture the cells onto a glass bottom dish (nonluminescent glass).
3. Remove the culture medium from the culture dish and add the same volume of **pre-warmed** staining solution.
4. Incubate at 37°C/ 5% CO<sub>2</sub> for 10min – 2hrs. Staining conditions should be optimized for your each cell type.
5. Wash the cells 3 times with culture medium after staining.
6. Visualize the cells using a fluorescence microscope at Ex 460 nm – 520 nm and Em 520 nm – 700 nm.

## Storage


Store POLARIC™ PLT-500c6 reagents at room temperature. Protect from moisture and light. Use only fresh solutions. Store staining solution at 4°C and use solution within two weeks.

## IV. Applications

**Fig 2. Rat Visceral Adipocyte**

A. differential interference contrast microscopy □ □  
B. POLARIC™ PLT-500c6 stained cell  
Mitochondria produces orange fluorescence, cell membrane produces green fluorescence, ER produces yellow fluorescence.

**Fig 3. Rat Cardiomyocyte double staining by Hoechst33342 and POLARIC™ PLT-500c6**

POLARIC™ PLT-500c6 can stain rat cardiomyocyte and non-cardiomyocyte with different colors. Cardiomyocytes stained with POLARIC™ PLT-500c6 can live and beat over for 2 months while retaining stain.

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COSMO BIO CO., LTD.