Rat astrocyte

Catalog No. AST01C

April 2, 2020

Summary

Astrocytes is one of the glial cells in the central nervous system that play key roles in developing the nervous system and the blood-brain barrier and in maintaining extracellular fluid homeostasis.

This product is cryopreserved primary astrocytes isolated from the fetal brain of SD rat. It can be used in the research for antioxidant, anti-inflammatory and gliosis and in the analysis for response to various cytokines.

This product was developed with the guidance of Dr. Yamanoha, Faculty of Science and Engineering, Soka University.

Precaution

- Read the manual carefully before beginning the culture.
- Cell culture should be performed under aseptic conditions. The biosafety level is [Level 1].
- Use the optimized culture medium (selling separately).

Warranty

Cosmo Bio warrants that product (cells) will be viable until the expiration date and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.

Component

Product Name	Quantity	Amount	Storage conditions	Expiration date
Rat astrocytes (cryopreserved cells)	1×10 ⁶ cells/vial	1	Liquid nitrogen	6 months

※ Store the cells in liquid nitrogen if not immediately used after receiving them.

※ Cells were confirmed positive for GFAP.

Cell origin

Cerebrum of SD rat embryo (E18-E20)

Optimized Culture Medium (selling separately)

Product name	Code	Volume	Storage conditions	Expiration date
Astrocyte Culture medium	ASTM	250 mL	- 20 °C	Described on bottle
			4°C after thawing.	3 months after thawing

Culture medium component: DMEM/F12, serum, etc.

Protocol of Cell thawing and seeding

※ This product cannot be passaged.

※ This product is shipped in dry ice packages. Begin culture of cells immediately upon receipt. Store the cells in liquid nitrogen for long term storage.

Materials Not Supplied

- Astrocyte Culture medium (cat#ASTM)

1. Thaw the frozen cryovial in 37°C water bath for 2 minutes.
2. Transfer the thawed cell suspension to a 50 mL conical tube containing 10 mL of Astrocyte Culture medium (cat#ASTM) and mix them. Then transfer 1 mL of culture medium in the same conical tube to the cryovial and collect the remaining cells.
3. Seed the cell at the density of 0.5 - 1×10^4 cells/cm² and culture in 37° C, 5% CO₂ incubator.

※ It is recommended to seed the cells without centrifugal washing. It may worsen cell adhesion and proliferation, because the cell damage by centrifuge after thawing is greater than that by residual cell freezing solution.

4. Replacement the medium with the pre-warmed culture medium at 37° C on the day after beginning culture. Then medium replacement should be done twice or three times a week.

※ Astrocytes reach 80 to 100% confluent 3 to 5 days after beginning culture.

※ When cultured in glass slides or dishes, seed the cells at the density of $2-5 \times 10^4$ cells/cm² because astrocytes have poor adhesion to the glass surface. The cultured cells reach 80 to 100% confluent for about 1 week.

Frequently asked questions

Q. Can the astrocytes be passaged?

A. No, passage is not recommended. Passage result in the predominance of highly proliferating non-astrocytic cells because this product is a primary cell taken from an organ.

Q. Is the purity of astrocyte measured?

A. The purity of astrocyte has not been determined. However, immunostaining for GFAP is performed and it is confirmed that more than 90% of cells are GFAP positive by visual inspection.

Technical information


Figure 1. Cell morphology

A. Anti-GFAP

B. Nuclear stained

C. Phase contrast

GFAP (Glial Fibrillary Acidic Protein) is an intermediate filament with a molecular weight of about 50kDa that is specifically found in astrocytes. Its expression level is increased in neurological diseases such as brain injury, dementia, prion disease, and multiple sclerosis.

Reference

- 1) Miller, R. H., Ffrench-Constant, C., and Raff, M. C. (1989). The macroglial cells of the rat optic nerve. *Annu. Rev. Neurosci.* 12, 517-534.

COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN

Phone: +81-3-5632-9610

FAX: +81-3-5632-9619

URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA

email: info@cosmobiousa.com

Phone/FAX: (+1) 760-431-4600

URL: www.cosmobiousa.com

For research use only. Not for clinical diagnosis.

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COSMO BIO CO., LTD.