For research use only. Not for clinical diagnosis.

Astrocyte, Mouse

Catalog No. AST02C

August 9, 2024

I. Product Overview

Astrocyte is one of the glial cells in the central nervous system, which plays an important role in the construction of the nervous system, homeostasis of extracellular fluid, and formation of the blood-brain barrier.

This product is a frozen astrocyte isolated from fetal brain of C57BL/6N mouse. It can be used for experiments such as anti-oxidation, anti-inflammation, response to various cytokines, and gliosis.

II. Precautions before use

Please review this manual before use.

This product should be operated under aseptic conditions. The biosafety level is [Level 1].

Please use the dedicated medium (sold separately) for culturing this product.

III. Product Warranty

The product is warranted against growth failure after the start of culture only when cells are properly stored in liquid nitrogen and cultured according to the manual using the dedicated medium and reagents.

The warranty period is within 6 months from the date of receipt of the product.

Please note that the warranty does not apply to any modifications to the medium or method of use, or to the use of re-frozen cells.

IV. Product Composition

Composition	capacity	Number of bottles	Storage Temperature	Expiration date
Mouse astrocytes (frozen cells)	1×10 ⁶ cells/vial	1 vial	Liquid nitrogen storage	6 months

Please store frozen cells in liquid nitrogen if not used immediately after receipt.

Cells are confirmed GFAP-positive.

V. Cell origin

Cerebrum of fetal C57BL/6N mice (E16-E17)

VI. Dedicated medium (sold separately)

Product name	Catalog No.	Volume	Storage Temperature	Expiration Date
Culture medium for astrocytes	ASTM	250 mL	Stored at -20°C (Store at 4°C after thawing)	Bottle description (stored at -20°C) 3 months after thawing (stored at 4°C)

Main components of culture medium: DMEM/F12, serum, others

VII. Procedure

*This product cannot be passaged.

Cell thawing and seeding

*Frozen cells for this product are shipped in dry ice packaging. Cells should be cultured immediately upon receipt. For long-term storage, store in liquid nitrogen.

[Preparation]

Culture medium for astrocytes: ASTM

1. Thaw the vial of frozen cells by warming the vial in 37°C water for 2 minutes.
2. Transfer the thawed cell suspension to a 50 mL centrifuge tube containing 10 mL of culture medium for astrocytes, mix the thawed cell suspension with the culture medium in the centrifuge tube, and then remove 1 mL of the culture medium from the centrifuge tube and collect the cell suspension by prewashing the vials.
3. Seed the above cell suspension at a density of 0.5 to 1.0×10^4 cells/cm² and culture in a 37°C incubator in the presence of 5% CO₂.

*We recommend seeding without centrifugal washing because the damage to cells caused by centrifugal washing after thawing is greater than the damage caused by residual frozen cell suspension, which may result in poor cell adhesion and proliferation.

4. The medium should be changed once the day after seeding, and 2 or 3 times a week thereafter, using medium warmed to 37°C.

*The cells will be 80% to 100% confluent 3 to 5 days after seeding.

*Because astrocytes do not adhere well to glass surfaces, when cultured on glass slides or glass dishes, seed cells at a density of $2.0\text{-}5.0 \times 10^4$ cells/cm² and culture for about 1 week until 80-100% confluence is achieved.

VIII. Frequently Asked Questions

Q. Can the cells be passaged?

A. This product is a primary cell isolated from organs, so we do not recommend passaging, as the proliferative non-astrocyte cells will predominate in the culture.

Q. Is the purity of astrocytes measured?

A. We do not measure the purity of astrocytes, but we have confirmed that almost 90% or more of the cells are GFAP-positive by visual confirmation in immunostaining.

IX. Technical Information




Fig. 1. Cell morphology
A. Anti-GFAP antibody staining
B. Nuclear staining C. Bright field

GFAP (Glial Fibrillary Acidic Protein) is an intermediate-sized filament with a molecular weight of approximately 50 kDa found specifically in astrocytes, and its expression is known to increase in neurological diseases such as brain damage, dementia, prion diseases and multiple sclerosis.

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