# Bone Marrow-derived mesenchymal stem cell

(Rat Bone Marrow-derived)

Catalog No. MSB01C

August 26, 2020

## Principle

Mesenchymal Stem Cells (MSC), being one kind of somatic stem cells, can be isolated from adipose tissue, bone marrow, umbilical cord, dental pulp or the like. In 2006, Dominici et al. have defined, as Mesenchymal Stem Cells, the cells satisfying the following three criteria: (1) MSC must be adherent to plastic culture vessels; (2) MSC must express CD105, CD73 and CD90 as positive markers and lack expression of CD45, CD34, CD11b, CD79a or CD19 and HLA-Class II(DR) as negative markers; and (3) MSC must have differentiation potency to osteoblasts, adipocytes and chondroblasts<sup>1)</sup>. Mesenchymal Stem Cells are known to have excellent characteristics of, in addition to excellent autologous growth, having pluripotency not only to bone, cartilage, and fat but also to liver or the like<sup>2)3)5)7)</sup>, lacking the expression of MHC Class II<sup>3)</sup>, having immunomodulatory capabilities<sup>4)</sup>, collecting on a diseased region to secrete various cytokines and growth factors<sup>5)</sup>, being effective in promoting tissue regeneration and repair<sup>6)7)</sup>, and the like.

This product is Mesenchymal Stem Cells (BMSC / BSC: Bone Marrow Delivered Mesenchymal Stem Cells / Bone Marrow Delivered Stromal Cells) and is prepared by collecting, from the cell groups (product name: Bone Marrow Cell Culture Kit BMC01) isolated from bone marrow of adult SD rats, the MSC fractions being positive in CD44, CD73, CD90 and CD105 and negative in CD14, CD31 and CD45, and then cryopreserving the resultant material.

## Warranty

**Cosmo Bio warrants that product (cells) will be viable until the expiration date, and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.**

## Kit Components

Components	Size	Quantity	Storage	Expiration date
Rat Bone Marrow-derived Mesenchymal Stem Cell (frozen cells)	$0.5 \times 10^6$ cells/vial	1	Liquid nitrogen vapor phase	6 months

※This product is delivered with dry ice package. Immediately store the received cells in liquid nitrogen.

## Optimized Culture Medium (selling separately)

Catalog No.	Product Name	Size	Quantity	Storage Conditions	Expiration date
MSB-GM	Bone Marrow-derived Mesenchymal Stem Cell (BMSC) Growth Medium	500 mL	1	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
		250 µL	2	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 2weeks after thawing (stored at 4°C)

(Manufacturer Abbreviation: PMC)

## Other Items Not Packaged

### Trypsin

Product Name	Code No.	Size
Trypsin EDTA Solution C (0.05%), EDTA (0.02%)	03-053-1B	100 mL

(Biological Industries Ltd. Manufacturer Abbreviation: BLG)

### Cell Freezing Medium

Product Name	Code No.	Size
COS banker	COS-CFM01	120 mL

(Manufacturer Abbreviation: KOJ)

**Collagen Coating Solution** to be used during differentiation. Unnecessary upon amplifying MSC in the growth medium

Product Name	Code No.	Size
Collagen Coating Solution	SCO	100 mL

(Manufacturer Abbreviation: PMC)

## Adipocytes Staining Kit

Product Name	Code No.	Components	Volume
Lipid Assay Kit	AK09F	Oil Red O Solution Solvent for Oil Red O Extraction	150 mL × 2 200 mL × 2

(Manufacturer Abbreviation: PMC)

## Related Products

### List of Single Articles for Adipogenic Differentiation Medium (Serum-Containing Medium)

Product Name	Size	Quantity	Code No.	Storage	Shelf Life
Mesenchymal Stem Cell (MSC) Adipogenic Differentiation Medium Set	MSC-ADGM 125 mL MSC-ADDM 100 mL MSC-ADMM 125 mL	1 set	MSC- ADM	-20°C (4°C after thawing)	- Written on the bottle (stored at - 20°C) 3months after thawing (stored at 4°C)
Growth Medium (Mesenchymal Stem Cell (MSC) Adipogenic Differentiation Medium Set)	500 mL	1	MSC- ADGM		
Differentiation Medium (Mesenchymal Stem Cell (MSC) Adipogenic Differentiation Medium Set)	500 mL	1	MSC- ADDM		
Maintenance Medium (Mesenchymal Stem Cell (MSC) Adipogenic Differentiation Medium Set)	500 mL	1	MSC- ADMM		

(Manufacturer Abbreviation: PMC)

## Bone-marrow Derived Mesenchymal Stem Cell (BMSC)

### Chondrogenic Differentiation Medium

Components	Size	Quantity	Code No.	Storage	Shelf Life
Chondrogenic Differentiation Medium	50 mL	1	MSC-CHB	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
Supplement	250 µL	1			- Written on the bottle (stored at -20°C) 1months after thawing (stored at 4°C)

(Manufacturer Abbreviation: PMC)

**Other Mesenchymal Stem Cell** Products are partially listed in the Table. Refer to the web site for details.

Product Name	Code No.	Size	Quantity	Storage	Expiration date
Adipose tissue-derived mesenchymal stem cell	MSA01C	$1 \times 10^6$ cells/vial	1	Liquid nitrogen vapor phase	6 months

(Manufacturer Abbreviation: PMC)

### Cell Culture Method, cultured in 10 cm dish

#### Preparation of Growth Medium

1. Thaw the growth medium at room temperature and transfer 100 mL thereof to a new medium bottle. The remaining medium can be stored for 3 months at 4°C.
2. Thaw the supplement when used, and add it to the growth medium. Then, 100 mL of growth medium can be prepared by one tube of 250 µL supplement. The medium after addition of the supplement can be stored for two weeks at 4°C.

#### Thawing and Seeding of Cells

1. Thaw the frozen cells in warm water at 37°C for 80 to 90 seconds (to a degree at which about 1 mm of ice lump remains.).  
※ The frozen cells in this kit are delivered in dry ice packaging. When the culture is not made immediately, immediately keep the cells at -80°C or in a liquid nitrogen container. When the cells are stored for a long time of six months or more, store the cells in liquid nitrogen.
2. For contamination prevention, wipe off moisture attached outside the vial with an alcohol swab, and then open the vial, and transfer the thawed cell liquid to a 15 mL conical tube having 10 mL of growth medium returned to room temperature therein. Gently mix the liquid, and then centrifuge it at  $100 \times g$  for five minutes.
3. Remove the supernatant by suction, add 9 mL of growth medium thereto and disperse the cells thereinto ( $0.55 \times 10^5$  cells/mL).

4. Put 3 mL of cell suspension and 7 mL of growth medium in a 10 cm dish, and gently shake the resultant mixture to uniformly disperse the cells thereinto. Incubate the cells in the presence of 5% CO<sub>2</sub> in an incubator at 37°C (one vial of frozen cells corresponds to three 10 cm dishes.).
5. Replace the medium with the growth medium kept warm at 37°C on the following day from seeding (the first day of culture).
6. Conduct subculture and differentiation experiments on the cells or cryopreserve the cells at a 70 to 80% confluent state (the 2<sup>nd</sup> to 4<sup>th</sup> day of culture).

#### **Cell Collection and Subculture (in the case of 10 cm dish)**

1. Remove the old medium from the cell culture vessel at the 70 to 80% of confluent state (the 2<sup>nd</sup> to 4<sup>th</sup> day of culture), and rinse the cells using HBSS or PBS (-) once or twice.
2. Put 3 mL of 0.05% trypsin-EDTA in the dish and leave it at 37°C for one to two minutes while observing the state.
3. Gently tap the culture vessel horizontally, and then confirm peeling thereof with a microscope. If not, leave it at 37°C further for one minute.
4. Add 10 mL of growth medium thereto, deactivate trypsin, and collect the cells into the centrifugal tubes by pipetting.
5. Centrifuge the tubes at 100 × g for five minutes, and then remove the supernatant, and add a suitable amount of medium to suspend the cells thereinto to disperse the cells.
6. Count the number of cells, and then adjust the cell density necessary to a level for the objective experiments to conduct seeding and differentiation induction experiments.

#### **Precautions:**

- Subculture twice or more is not recommended because growth rate or differentiation efficiency is significantly reduced.
- Trypsin-EDTA processing for 10 minutes or more causes deterioration of the state of cells.

### **Preparation of Frozen Stock**

1. Perform the steps 1 to 4 in the method of cell collection and subculture.
2. Centrifuge the tubes at  $100 \times g$  for 5 minutes at room temperature, and then remove the supernatant, and add COS banker (KOJ, Code No. OCS-CFM01) thereto to be  $1 \times 10^6$  cells/mL, and suspend the cells thereinto.
3. Put the resultant suspension in a cell freezing tube at 1 mL/tube, freeze the suspension at  $-80^\circ\text{C}$  by using a cell freezing container (BM Equipment Co, Ltd., Code No.: BCS-136 or equivalent), and store the tubes.
4. Leave the tubes to stand at  $-80^\circ\text{C}$  for one day or more, and transfer the tubes into liquid nitrogen.

### **Differentiation into Adipocytes by Adipogenic Differentiation Medium**

Precaution: After adipogenic differentiation, the cells are significantly easily peeled. In such a case, collagen coating is recommended.

#### **Collagen Coating.**

1. Put the collagen coating solution (PMC, Code No.: SCO) in a culture vessel to be used at a degree of being covered, leave the resulting material to stand for one hour or more (overnight for the glass bottom), and then remove the solution.
2. Rinse the vessel with PBS twice, and then put PBS therein, and leave the resultant material to stand for one or more hours, and then use the resultant material. If the material is not used immediately, fill the vessel with PBS, which can be kept for several days at room temperature.

#### **Adipogenic Differentiation**

1. Incubate the cells in the growth medium in advance to a 70 to 80% confluent state according to the method of thawing and seeding of cells (cell density:  $1.0 \times 10^4$  cells/cm<sup>2</sup>).
  2. Replace the medium with the growth medium of the adipogenic differentiation set, kept warm at  $37^\circ\text{C}$ , and incubate the cells for one to two days until the cells reach a 100% confluent state.
  3. Replace the medium with the differentiation medium of the adipogenic differentiation set, kept warm at  $37^\circ\text{C}$ , and incubate the cells for three days.
  4. Replace the medium with the maintenance medium of the adipogenic differentiation set, kept warm at  $37^\circ\text{C}$ , and incubate the cells for two days.
  5. Repeat the step 3 again. In this stage, small fat globules begin to emerge.
  6. Replace the medium with the maintenance medium, and observe the growth of the fat globules.
- × Accumulated adipocyte can be stained with the Lipid Assay Kit (Fig. 2).

## Differentiation into Chondrocytes by Chondrogenic Differentiation Medium

### Preparation of Medium

Thaw the chondrogenic differentiation medium (50 mL) at room temperature before use, and add one vial of thawed supplement.

### Chondrogenic Differentiation

1. Seed the cells in a 10 cm dish or the like in advance according to the method of Cell Collection and Subculture, collect the cells which reached the 70 to 80% confluent state, and count the number of cells.
2. Centrifuge the tubes at  $100 \times g$  for five minutes at room temperature, remove the supernatant, suspend the remaining solid with the chondrogenic differentiation medium, and adjust the cell density to be  $5.0 \times 10^5$  cells/mL. Dispense the resultant material by 0.5 mL for each 15 mL tube. ( $2.5 \times 10^5$  cells per tube).
3. Centrifuge the tubes at  $100 \times g$  for five minutes at room temperature, and then lightly loosen the cap of each tube in a clean bench so as to exchange gases, and incubate the cells in the presence of 5%  $\text{CO}_2$  in the incubator at  $37^\circ\text{C}$ . On that occasion, do neither suck the supernatant nor re-suspend the pellets.
4. Leave the pellets to stand for 24 hours, and replace the medium with the chondrogenic differentiation medium to which the supplement is added, every 2 to 3 days. When the pellets are attached to the tube during medium exchange, gently tap the tube to drop the pellets.
5. Chondrocytic cell mass is formed in two to four weeks.
6. The formed chondrocytic cell aggregates are fixed by formalin, embedded with paraffin, and then can be stained with Alcian Blue or the like (Fig. 3).

### Cell Morphology


Fig. 1 Cells on the third day of culture (Culturing in BMSC growth medium)


Fig. 2 Adipogenically differentiated cells (Staining with Lipid Assay Kit)

Fig. 3 Chondrogenically differentiated cells (Staining with Alcian Blue)

## References

1. Dominici M et al. Minimal criteria for defining multipotent mesenchymal stromal cells. The International Society for Cellular Therapy position statement. *Cytotherapy*. 2006; 8(4):315-7.
2. Uccelli A et al. Mesenchymal stem cells in health and disease, *Nat. Rev. Immunol.*, 8, 726, 2008.
3. Le Blanc K. Immunobiology of human mesenchymal stem cells and future use in hematopoietic stem cell transplantation. *Biol Blood Marrow Transplant*. 2005 May; 11(5):321-34.
4. Alma J N et al. Immunomodulatory properties of mesenchymal stromal cells. *Blood*. 2007 Nov 15; 110(10): 3499-506
5. Banas A et al. IFATS collection: In vivo therapeutic potential of human adipose tissue mesenchymal stem cells after transplantation into mice with liver injury. *Stem Cells*. 2008 Oct; 26(10): 2705-12.
6. Wollert.K.C. et al. *Lancet*, 364: 141, 2004
7. Sato.Y.et al. Human mesenchymal stem cells xenografted directly to rat liver are differentiated into human hepatocytes without fusion. *Blood*, 106, 756, 2005

**COSMO BIO Co., LTD.**

**【JAPAN】**

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

**COSMO BIO USA**

**【Outside Japan】**

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobioussa.com](mailto:info@cosmobioussa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobioussa.com](http://www.cosmobioussa.com)

**For research use only. Not for clinical diagnosis.**

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