Manual for handling the Optimized Medium and Collagen Coating for AMSC and BMSC

October. 17, 2025

Please read this user manual thoroughly before use. This product is for research use only.

I. Product Overview

These products are composition optimized culture medium for AMSC and BMSC, which contain serum and antibiotics.

II. Storage

These products will be shipped frozen. Please store below -20°C immediately after receipt. They should be stored at 4°C after thawing, and used as soon as possible.

The expiration date for these products will be noted on the label, which applies only for when they are unopened. After thawing, they should be stored at 4°C , and used as soon as possible. It is highly recommended to aliquot out to avoid freeze-thaw cycles. It should be noted that repeated freeze-thaw of these products will damage or cause quality loss of these products.

III. Precautions and Procedure for use

Please thaw these products in the refrigerator, so that they will thaw slowly. Never warm above 38°C . If necessary, aliquot out necessary amount, then warm it at 37°C .

This product cannot be warranted if it is left for a while under warm condition, or repeatedly warmed up, because of the quality loss.

Optimized Culture Medium for BMSC

Catalog No.	Product Name	Size	Quantity	Storage Conditions	Expiration date
MSB-GM	Bone Marrow-derived Mesenchymal Stem Cell (BMSC) Growth Medium	500 mL	1	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
		250 μL	2	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 2weeks after thawing (stored at 4°C)

Preparation of Growth Medium

1. Thaw the growth medium at room temperature and transfer 100 mL thereof to a new medium bottle. The remaining medium can be stored for 3 months at 4°C .
2. Thaw the supplement when used, and add it to the growth medium. Then, 100 mL of growth medium can be prepared by one tube of 250 μL supplement. The medium after addition of the supplement can be stored for two weeks at 4°C .

COSMO BIO CO., LTD.

Inspiration for Life Science

Optimized Culture Medium for AMSC

Catalog No.	Product Name	Size	Quantity	Storage Conditions	Expiration date
MSA-GM	Adipose Tissue-derived Mesenchymal Stem Cell (AMSC) Growth Medium	200 mL	1	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
		50 µL	2	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 2weeks after thawing (stored at 4°C)

Preparation of Growth Medium

1. Thaw the growth medium at room temperature and transfer 100 mL thereof to a new medium bottle. The remaining medium can be stored for 3 months at 4°C.
2. Thaw the supplement when used, and add it to the growth medium. Then, 100 mL of growth medium can be prepared by one tube of 50 µL supplement. The medium after addition of the supplement can be stored for two weeks at 4°C.

Collagen Coating Solution to be used during differentiation. Unnecessary upon amplifying MSC in the growth medium

Product Name	Code No.	Size
Collagen Coating Solution	SCO	100 mL

Collagen Coating.

1. Put the collagen coating solution (PMC, Code No.: SCO) in a culture vessel to be used at a degree of being covered, leave the resulting material to stand for one hour or more (overnight for the glass bottom), and then remove the solution.
2. Rinse the vessel with PBS twice, and then put PBS therein, and leave the resultant material to stand for one or more hours, and then use the resultant material. If the material is not used immediately, fill the vessel with PBS, which can be kept for several days at room temperature.

COSMO BIO CO., LTD.

Inspiration for Life Science

List of Single Articles for Adipogenic Differentiation Medium for AMSC and BMSC(Serum-Containing Medium)

Product Name	Size	Quantity	Code No.	Storage	Shelf Life
Differentiation Medium (Mesenchymal Stem Cell (MSC) Adipogenic Differentiation Medium Set)	500 mL	1	MSC-ADDM	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
Maintenance Medium (Mesenchymal Stem Cell (MSC) Adipogenic Differentiation Medium Set)	500 mL	1	MSC-ADMM		

Bone-marrow Derived Mesenchymal Stem Cell (BMSC) Chondrogenic Differentiation Medium

Components	Size	Quantity	Code No.	Storage	Shelf Life
Chondrogenic Differentiation Medium	50 mL	1	MSC-CHB	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
Supplement	250 µL	1			- Written on the bottle (stored at -20°C) 1months after thawing (stored at 4°C)

Preparation of Medium

1. Thaw the chondrogenic differentiation medium (50 mL) at room temperature before use, and add one vial of thawed supplement.

COSMO BIO CO., LTD.

Inspiration for Life Science

Adipose tissue-derived Mesenchymal Stem Cell (AMSC)

Chondrogenic Differentiation Medium

Components	Size	Quantity	Code No.	Storage	Shelf Life
Chondrogenic Differentiation Medium	50 mL	1	MSC-CHA	-20°C (4°C after thawing)	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
Supplement	500 µL	2			- Written on the bottle (stored at -20°C) 1months after thawing (stored at 4°C)

Preparation of Medium

1. Thaw the chondrogenic differentiation medium (50 mL) at room temperature before use, and add two vial of thawed supplement.

For research use only. Not for clinical diagnosis.

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COSMO BIO Co., LTD.