For research use only. Not for clinical diagnosis.

Osteogenesis Culture Kit (Mouse)

[Catalog No. OGC11]

July 1, 2024

I. Product Overview

There are hematopoietic stem cells and bone marrow stromal cells in bone marrow. Bone marrow stromal cells contain undifferentiated mesenchymal stem cells that can differentiate into a variety of cell types such as osteoblasts chondrocytes adipocytes and so on.

This product contains cryopreserved cells isolated from mouse bone marrow and two types of culture medium. The cells in this product can be grown using Growth Medium, and then can be differentiated into mature osteoblasts, which form calcified nodules, using Osteogenesis Medium.

II. Precautions Before Use

Please make sure to review this manual before using the product.

Product should be used under [aseptic operation]. The biosafety level is [Level 1].

Please use the dedicated medium supplied or sold separately for culturing this product.

III. Product Warranty

We guarantee against growth failure after the start of culture only if cells have been properly stored in liquid nitrogen, and cultured according to the manual using the dedicated medium and reagents.

The warranty period is [within 6 months of receiving the product].

Please note that the warranty does not apply if there have been changes to the medium or the method of use, if cells have been re-frozen cells have been used.

IV. Product Components

Product Name (Components)	Size	Quantity	Storage	Expiration
Mouse Bone Marrow Stromal Cells	1×10 ⁶ cells/vial	1	Liquid nitrogen vapor storage	6 months
Growth Medium	125 mL	1	-20°C (Frozen) 4°C (after thawing)	Labeled on the bottle (when stored at -20°C) 3 months after thawing (when stored at 4°C)
Osteogenesis Medium	250 mL	1		

(This product is shipped with dry ice package. If the product is not to be used immediately upon receipt, please store it in liquid nitrogen.)

V. Cell Origin

Cells are derived from the bone marrow of ICR mice.

VI. Dedicated Media (Serum Containing Media) (sold separately)

Product Name	Catalog No.	Size	Quantity	Storage	Expiration
Growth Medium	OGCMG	250 mL	1	-20°C (Frozen) 4°C (after thawing)	Labeled on the bottle (when stored at -20°C) 3 months after thawing (when stored at 4°C)
Osteogenesis Medium	OGCMO	250 mL	1		

(Medium main ingredients: α-MEM, serum, antibiotics, others.)

VII. Related Products

Collagen Coating Solution

Product Name	Catalog No.	Size	Quantity
Collagen Coating Solution	SCO	100 mL	1

Osteocytes Staining Kit

Product Name	Catalog No.	Components	Size	Quantity
Alkaline Phosphatase Staining Kit	AK20	Substrate-containing Buffer	50 mL × 1	1 Set
		Chromogenic Substrate	10 vials (For 5 mL)	
Calcified Nodule Staining Kit	AK21	Buffer (100 ×)	60 mL × 1	1 Set
		Staining Substrate	10 vials (For 10 mL)	

VIII. Instructions For Use

※ This product [is not capable of being subcultured].

IX. Protocol

Collagen Coating

[Requirements]

- Culture vessels
- Collagen coating solution (Cat No.: SCO)
- PBS(-) or HBSS(-)
- Sterile pipettes, conical tubes, and other culture equipment

1. Add enough collagen coating solution to cover the surface of the culture vessel (approximately 0.3–0.5 mL/well for a 24-well plate); thereafter, immediately remove the solution. For glass-bottom vessels, let the solution stand overnight before removing it.
2. Wash twice with PBS(-) or HBSS(-), then add PBS(-) or HBSS(-) and let stand for at least one hour before use. If not using immediately, fill the vessel with PBS(-) or HBSS(-) and store at room temperature for several days.

Cell Thawing and Seeding

※ For culturing in 24-well plates.

[Requirements]

- Mouse Bone Marrow Stromal Cells
- Growth Medium
- Osteogenesis Medium
- Collagen-coated 24-well plate
- Sterile pipettes, conical tubes, and other culture equipment

1. Remove the cryovial containing the frozen cells from liquid nitrogen storage, and immediately place it in a 37 °C water bath. Quickly thaw the cells (within 90–120 seconds) by gently swirling the vial in the 37 °C water bath until only a small bit of ice remains in the vial.
2. Transfer the vial to a clean bench. Before opening, wipe the outside of the vial with 70% ethanol.
3. Transfer the thawed cell solution to a 15-mL conical tube containing 10 mL of pre-warmed (37 °C) growth medium. After gentle mixing, centrifuge the solution at 600 ×g for 5 minutes at 4 °C.
4. Remove and discard the supernatant, then add 5 mL of pre-warmed (37 °C) growth medium and resuspend the cells.
5. Dispense 0.5 mL of cell suspension to each well of 24-well plate, and culture the cells in a 37 °C incubator with 5% CO₂.

- ※ In the case of using other types of culture plate, adjust cell density to 2.5×10^4 cells/ cm^2 . Seeding at a lower cell density than specified will adversely affect cell proliferation.
6. Change the medium in the following day, and every other day thereafter until the cells become confluent.
 - ※ Cells become confluent within 3-4 days.
 7. Once cells become confluent, replace the medium with pre-warmed (37 °C) Osteogenesis Medium to promote the formation of calcified nodules. Replace the medium every 2 or 3 days.
 - ※ Approximately 4-7 days after medium is changed into Osteogenesis Medium, cells whiten.
 - ※ Adipocytes appear from a part of cells.

Cell Subculturing

- ※ This product is not recommended for subculturing as the seeded cells are difficult to detach.

X. Technical Information



Figure 1: Cell Morphology

A. Cells on the 1st day of culture / B. Confluent / C. Cells on the 3 weeks of culture in osteogenesis medium (cultured in a 35 mm dish)

✧ Application example

Cells are cultured according to protocol and subjected to activity staining of alkaline phosphatase or calcified nodules.

The activity of alkaline phosphatase is detected by Alkaline Phosphatase Staining kit (Cat No.: AK20) and the calcium deposit is detected by Calcified nodule Staining kit (Cat No.: AK21).

XI. References

- 1) SL Cheng *et al.* Differentiation of human bone marrow osteogenic stromal cells in vitro: induction of the osteoblast phenotype by dexamethasone. *Endocrinology* 134(1) p277-286 (1994)
- 2) Ohgushi *et al.* In vitro bone formation by rat marrow cell culture. *Journal of Biomedical Materials Research Part A* 32(3) p333-340 (1996)

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