# Osteoclasts , Human

Catalog No. OSC15C

July 1, 2024

## I. Product Overview

Bone metabolism is composed of balanced osteogenesis and bone resorption. Research studies have shown that bone marrow cells can be differentiated into osteoclasts using M-CSF (Macrophage Colony Stimulating Factor) and RANKL (Receptor Activator of NF kappa B Ligand).

This product consists of human-derived osteoclast precursors (frozen). When cultured in a medium containing M-CSF and RANKL, which are sold separately, osteoclasts are formed.

It is useful for experiments on osteoclast formation, studies on bone resorption function, and research on conditions such as osteoporosis, rheumatoid arthritis, and Paget's disease.

\*The human cells used as raw materials in the cells have been confirmed to meet the following criteria:

- Informed consent has been properly obtained.
- Donor information confidentiality is maintained.
- The donor is negative for major infectious diseases (HIV, HBV, HCV, etc.).

\*These cells are shipped in dry ice packaging. Upon receipt, please promptly store the cells in liquid nitrogen. If liquid nitrogen storage is not possible, please proceed with cell culture immediately.

## II. Precautions Before Use

Please make sure to review this manual before using the product.

Product should be used under [aseptic operation]. The biosafety level is [Level 2].

Please use the dedicated medium sold separately for culturing this product.

## III. Product Warranty

**We guarantee against growth failure after the start of culture only if cells have been properly stored in liquid nitrogen, and cultured according to the manual using the dedicated medium and reagents.**

**The warranty period is [within 6 months of receiving the product].**

**Please note that the warranty does not apply if there have been changes to the medium or the method of use, if re-frozen cells have been used.**

## IV. Components

Product Name	Size	Quantity	Storage Conditions	Expiration date
Osteoclast precursor, Human, cryopreserved	1.5x10 <sup>6</sup> cells/vial	1 vial	Liquid Nitrogen vapor phase	6 months

\*Shipping: dry ice

## V. Dedicated Media (Serum Containing Media) (sold separately)

Product Name	Catalog No.	Size	Quantity	Storage Conditions	Expiration date
Osteoclast Wash Medium	OSCMW	100 mL	1	-20°C Freezer	- Written on the bottle
Osteoclast Culture Medium (w/ RANKL and M-CSF)	OSCMHB	30 mL	1	-20°C Freezer	- Written on the bottle

**Culture Medium components:  $\alpha$ -MEM, FBS, antibiotics, etc.**

\*If the test period is long, please aliquot the medium after thawing and store it frozen. Re-freezing after thawing is only allowed once. Avoid further freezing and thawing as it may cause a deterioration in quality.

## VI. Related Products (sold separately)

Product Name	Catalog No.	Size	Storage	Expiration date
RANKL <sup>*1</sup>	AK30	10 $\mu$ g	-70°C (frozen)	6 months
M-CSF <sup>*2</sup>	AK39	20 $\mu$ g	-70°C (frozen)	6 months

\*1 RANKL is identical to the substance included in human osteoclast culture medium. It is expressed and purified as a GST fusion protein from human-derived soluble RANKL.

\*2 M-CSF is identical to the substance included in human osteoclast culture medium.

## VII. General Information

Organism	<i>Homo sapiens</i> , human
Tissue	Bone Marrow
Cultural Properties	Adherent
Biosafety Level	2
Virus Check	HIV-1(-), HBV (-), HCV (-)

## VIII. Materials required but not provided

- ✧ Osteoclast Wash Medium (Cat No. OSCMW)
- ✧ Human Osteoclast Culture Medium (Cat No. OSCMHB)
- ✧ Variable volume pipettes
- ✧ Culture vessels
- ✧ 15 ml centrifuge tube

## IX. Protocol

### A) Cultured with the 96-well culture plate

- 1) Carefully remove the cryovial from liquid nitrogen and thaw cells in a water bath at 37°C with gently shaking.
- 2) Transfer thawed cells into a 15 ml centrifuge tube containing 10 ml of Osteoclast Wash Medium (Cat No. OSCMW) and Transfer 1mL of culture medium in the same conical tube back to the cryovial and pour the contents back to 15mL conical tube.
- 3) Centrifuge for 5 minutes at 4°C at 200 x g for 5 minutes.
- 4) Remove the supernatant, and re-suspend the cell pellet in 2.5 to 5.0 ml of Human Osteoclast Culture Medium (Cat No. OSCMHB).
- 5) Dispense 100  $\mu$ L of the cell suspension to each well of 96-well culture plate, and incubate the flask at 37°C under 5% CO<sub>2</sub> and 100% humidity.
- 6) Three days later, Replace the medium with fresh culture medium.
- 7) Change the medium every other day.  
※Osteoclasts will begin to fuse and form osteoclasts after 4 or 7 days of incubation

## Technical information

1. For experiments on osteoclastogenesis regulation factors, add the regulatory factors to the culture medium and culture in the wells for 7 days. Afterward, perform tartrate-resistant acid phosphatase (TRAP) staining to measure the number of osteoclasts. The TRAP staining kit (Cat No. AK04F), available separately, is a convenient tool for staining tartrate-resistant acid phosphatase, a marker enzyme for osteoclasts (Fig 2).
2. Pit Formation Assay Method (Staining Method)  
After culturing osteoclast precursors on dentin slices for 2 to 4 weeks, collect the slices and ultrasonicate them in 5 ml of 1 M ammonia solution to destroy the cells. After treatment, remove the slices from the ammonia solution, directly stain them with Mayer's hematoxylin solution for 1 minute, rinse with water, dry, and then measure the total area of resorption pits formed by osteoclasts (Figure 3-1, 3-2).
3. Pit Formation Assay Method (SEM Method)  
Perform a series of operations on dentin slices similar to the Pit Formation Assay Method (Staining Method), and observe the surface of the dried slices using a scanning electron microscope (SEM) to measure the area of resorption pits. For clearer images of resorption pits, SEM observation is recommended (Figure 4).


Fig.1 Osteoclast induced of differentiation with M-CSF/RANKL


Fig. 2 TRAP staining of human osteoclast

Fig.3-1 Pit on the slice of ivory.  
(HE staining)


Fig.3-2 Pit on the slice of ivory.  
(HE staining)


Fig.4 Pit on the slice of ivory.  
(SEM picture)

## X. References

- 1) Takeshita et al. (2000) Identification and Characterization of the New Osteoclast Progenitor with Macrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. *J. Bone Miner. Res.* 15, 1477-1488.
- 2) Scheven et al. (1998) A sequential culture approach to study osteoclast differentiation from nonadherent porcine bone marrow cells. *In Vitro Cell Dev Biol Anim* 34, 568-577.
- 3) Martha et al. (1995) Enhanced Expression of  $\alpha$ V Integrin Subunit and Osteopontin during Differentiation of HL-60 Cells along the Monocytic Pathway. *Exp. Cell Res.* 216, 335-341.
- 4) Itonaga et al. (1999) 1,25-Dihydroxyvitamin D3 and Prostaglandin E2 Act Directly on Circulating Human Osteoclast Precursors. *Biochem. Biophys. Res. Commun.* 242, 703-709.

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1. Takeshita *et al.*(2000) Identification and Characterization of the New Osteoclast Progenitor with Macrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. J. Bone Miner. Res. 15, 1477-1488.
2. Scheven *et al.* (1998) A sequential culture approach to study osteoclast differentiation from nonadherent porcine bone marrow cells. In Vitro Cell Dev Biol Anim 34, 568–577.
3. Martha *et al.* (1995) Enhanced Expression of  $\alpha$ v Integrin Subunit and Osteopontin during Differentiation of HL-60 Cells along the Monocytic Pathway. Exp. Cell Res. 216, 335-341.
4. Itonaga *et al.*(1999) 1,25-Dihydroxyvitamin D<sub>3</sub> and Prostaglandin E<sub>2</sub> Act Directly on Circulating Human Osteoclast Precursors. Biochem. Biophys. Res. Commun. 242, 703-709.