Rep-HepG2 cells

【Code No.REP-HEPG2C】

For research use only. Not for clinical diagnosis.

March 4, 2025

I. Product Overview

Human hepatocytes are one of the essential experimental tools for pharmacokinetic studies. However, primary human hepatocytes, which express high hepatic function, do not proliferate in normal culture methods, and their hepatic function, including drug metabolism, is known to rapidly decline during the culture period. Furthermore, the limited availability of primary human hepatocytes tends to make them expensive, limiting their use in drug discovery research. On the other hand, hepatocellular carcinoma-derived cell lines, such as the HepG2 cell line, which is widely used due to its easy availability, have low metabolic activity, which is a problem. In order to solve these problems, new cell culture methods and the creation of hepatocyte substitutes for primary human hepatocytes that can proliferate are being investigated.

Although this product is a hepatocellular carcinoma-derived cell, it highly expresses drug-metabolizing enzyme activities (phase I and phase II reactions), which are the main functions of hepatocytes, and is expected to be used as a human hepatocyte model for basic research and drug screening studies.

※This product was developed under the technical guidance of Dr. Luc Gailhouste, Department of Molecular Cell Therapy, National Cancer Center Research Institute.

!Caution!

See attached “Please Read BEFORE USE”

All the products are required to perform during aseptic processing.

Use the specific media (sold separately) to culture this product.

This product has been confirmed negative for HIV-1, HIV-2, HBV, and HCV:

Use the optional medium for incubation.

～Safety of human-derived cells～

The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. However, not all infectious disease tests are performed, and the possibility of false negatives or infection cannot be completely excluded. Therefore, when handling, please take sufficient measures such as using a safety cabinet, wearing infection-protective gloves, goggles, and masks, etc., before conducting experiments.

This material and its container must be disposed of according to approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

Warranty

We guarantee against growth failure after the start of culture only when cells are properly stored in liquid nitrogen and cultured according to the manual using the dedicated medium and reagents.

Our company, Product Support (e-mail: bio-products@cosmobio.co.jp) for more information on how to get in touch with us.

The warranty period is [within 6 months from receipt of product].

Please note that the warranty does not apply if any changes are made to the culture medium or method of use, or if re-frozen cells are used.

II. Component

Components	Size	Quantity	Storage	Expiration date
Rep-HepG2 cells (frozen vial)	3.0×10 ⁶ cells/vial	1vial	Liquid Nitrogen	6month

If you do not use the frozen cells immediately after receipt, please store them in liquid nitrogen.

III. Cell Origin

Human Hepatocellular Carcinoma

IV. Specific Medium (sold separately)

Product name	Catalog no.	Componen	Volume	Quantity	Storage	Validity period
Rep-HepG2 Maintenance Medium	REP-HEPG2M	Medium	125mL	1	Stored at -20°C (After Thawing Stored at 4°C)	Indicated on bottle (Stored at -20°C) After Thawing 3 month (Stored at 4°C)
		Supplement	125μL	1	Stored at -20°C	

Main components of culture medium: DMEM, serum, others

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Cell thawing and seeding

※The following is the protocol for culturing in 24 well culture plate.

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- Collagen coated 24well culture plate.
- Rep-HepG2 Maintenance Medium (**Prepare just prior to use**)

※Prepare required culture media volum. (Add 10μL of supplement per 10ml of medium)

- Sterile pipettes, conical tubes, and other culture equipment

1. Remove the cryovial containing the frozen cells from liquid nitrogen storage and immediately place it into a 37°C water bath. Quickly thaw the cells by gently swirling the vial in the 37°C

water bath until there is just a small bit of ice left in the vial for 90~120 sec.

2. Transfer the vial it into a laminar flow hood. Before opening, wipe the outside of the vial with 70% ethanol.
3. Transfer the thawed cells dropwise into the 10mL pre-warmed Rep-HepG2 Maintenance Medium in 15mL centrifuge tube. Mix gently and centrifuge the cell suspension 170 x g for 3 min 4°C.
4. After the centrifugation, aspirate supernatant and resuspend the cells in Rep-HepG2 Maintenance Medium in 15mL
5. Dispense 0.5 mL of cell suspension to each well of 24-well plate and incubate the plate at 37°C in a 5% CO₂ humidified incubator (24well culture pate can usually be established from one vial containing .)
6. . Medium replacement should be performed next day, thereafter medium should be changed every 2 to 3 days.
7. Wait for 3-5 days before assaying for preferred assay.

Technical note


Figure 1. Phase contrast micrograph of Rep-HepG2 cells the day after culturing


Figure 2. Drug metabolizing enzyme activity of primary human hepatocytes as 1,
Ratio of enzyme activity of Rep-HepG2 cells (day 4 culture) and HepG2 cells

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