**For research use only. Not for clinical diagnosis.**


# Rat Subcutaneous White Adipocyte

**【 Catalog No. SAC01C 】**



July 1, 2024

## I. Product Overview


This product is a primary cell product that contains a frozen cell population, including precursor adipocytes, obtained from the subcutaneous adipose tissue of rats. It can be used for research across the entire process from precursor adipocytes to mature adipocytes. Additionally, adiponectin in the culture supernatant can be easily measured using a commercially available kit. Use this product for elucidating the mechanisms of metabolic syndrome and for exploring therapeutic drugs for such conditions.

Subcutaneous adipocytes  
(Cells on the 9th day of culture)

Visceral adipocytes  
(Cells on the 9th day of culture)

Epididymal adipocytes  
(Cells on the 9th day of culture)

## II. Precautions Before Use

Please make sure to review this manual before using the product.

Product should be used under [aseptic operation]. The biosafety level is [Level 1].

Please use the dedicated medium sold separately for culturing this product.

## III. Product Warranty

**We guarantee against growth failure after the start of culture only if cells have been properly stored in liquid nitrogen, and cultured according to the manual using the dedicated medium and reagents.**

**The warranty period is [within 6 months of receiving the product].**

**Please note that the warranty does not apply if there have been changes to the medium or the method of use, if re-frozen cells have been used.**

#### IV. Components

Product Name	Size	Quantity	Storage Conditions	Expiration date
Rat Subcutaneous White Adipocyte, Cryopreserved	3 x 10 <sup>6</sup> cells/vial	1 vial	Liquid nitrogen storage	6months

\*Shipping: dry ice

#### V. Cell Origin

Derived from abdominal subcutaneous white adipose tissue of adult SD rat

#### VI. Dedicated Media (Serum Containing Media) (sold separately)

Product name	Catalog No.	Size	Storage Conditions	Expiration date
Subcutaneous Adipocyte Culture Medium	SACMR	250mL	-20°C Freezer	- Written on the bottle

#### VII. Materials required but not provided

Variable volume pipettes

15 ml centrifuge tube

24-well culture plate

Subcutaneous Adipocyte Culture Medium (Cat No. SACMR)

#### VIII. Protocols

- 1) Thaw the Subcutaneous Adipocyte Culture Medium (Cat No. SACMR) at 4°C
- 2) Thaw vials of frozen cells by warming them in 37°C warm water for 2 minutes.
- 3) Transfer thawed cell solution to a 15 mL centrifuge tube containing 10 mL of Subcutaneous Adipocyte Culture Medium (Cat No. SACMR). Mix gently and centrifuge at 4°C at 200 x g for 5 minutes.
- 4) Remove the supernatant and re-suspend cells in 10 mL of Subcutaneous Adipocyte Culture Medium (Cat No. SACMR) and centrifuge at 4°C at 200 x g for 5 minutes.
- 5) Remove the supernatant and re-suspend the cell pellet in 12.5 mL of Subcutaneous Adipocyte Culture Medium (Cat No. SACMR)
- 6) Dispense 0.5 mL of cell suspension to each well of 24-well plate and incubate the plate at 37°C in a 5% CO<sub>2</sub> humidified incubator.
- 7) CO<sub>2</sub> humidified incubator.
- 8) Gently add 0.5 ml of Subcutaneous Adipocyte Culture Medium (Cat No. SACMR) to each well of 24-well plate on the following day.
- 9) Carefully change the medium 1 ml per well on the second day after seeding.
- 10) Change the medium every other day. Be careful to not disturb the cell layer.
  - I. Approximately 3-4 days of culture, preadipocyte culture becomes confluent.
  - II. Approximately 5 days of culture, cells turn to mature adipocytes.
  - III. Approximately 8 days of culture, cells become hypertrophic and start detaching from the bottom of the well.

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