For research use only. Not for clinical diagnosis.

Rat Schwann Cell Line (IFRS1)

[Catalog No. SWN-IFRS1C]

July 1, 2024

I. Product Overview

Schwann cells are the principal glia cells of the peripheral nervous system (PNS) and are involved in several aspects of PNS biology. In addition to providing trophic support for neuronal development and maintenance through the production of cytokines and growth factors, Schwann cells wrap around axons of motor and sensory neurons to create an insulating myelin sheath that aids impulse conductivity. In recent years, Schwann cells have increasingly known for their support to nerve regeneration while their dysfunction may contribute to neurodegenerative diseases such as ALS.

IFRS1s are spontaneously immortalized cell line established from long term culture of rat dorsal root ganglia and peripheral nerves. IFRS1 cells retain many characteristics of mature Schwann cells and are readily propagated in the supplied media. IFRS1 cells can be used to investigate many aspects of neuronal biology including nerve regeneration, neuroprotective mechanisms, and for the development of novel therapeutic approaches to neurological disorders.

II. Precautions Before Use

Please make sure to review this manual before using the product.

Product should be used under [aseptic operation]. The biosafety level is [Level 1].

Please use the dedicated medium sold separately for culturing this product.

III. Product Warranty

We guarantee against growth failure after the start of culture only if cells have been properly stored in liquid nitrogen, and cultured according to the manual using the dedicated medium and reagents.

The warranty period is [within 6 months of receiving the product].

Please note that the warranty does not apply if there have been changes to the medium or the method of use, if re-frozen cells have been used.

IV. Components

Product Name	Size	Quantity	Storage	Expiration
Rat Schwann Cell Line (IFRS1), Cryopreserved	5 x 10 ⁵ cells/vial	1 vial	Liquid nitrogen storage	6 months

*Store frozen cells in liquid nitrogen if not ready for use after receipt.

*Based on the license agreement of Tokyo Metropolitan Institute of Medical Science, IFRS1 cell are prohibited to provide (distribution, lending, transfer, licensing, etc.) to a third party.

V. General Information

Organism	Fischer344 Rat (adult)
Tissue	Dorsal Root Ganglion and Peripheral Nervous System
Cultural Properties	Adherent
Biosafety Level	1
Morphology	Spindle-shaped
Cell Type	Spontaneous Immortalization
Quality Check	Mycoplasma (-), Immunohistochemical reactivity of cell markers

VI. Dedicated Media (Serum Containing Media) (sold separately)

Product Name	Catalog No.	Volume	Storage	Expiration
Schwann Cell (IFRS1) Medium	SWNMR	250ml	-20°C (frozen) 4°C (after thawing)	Labeled on the bottle (when stored at -20 °C) 3 months after thawing (when stored at 4 °C)

VII. Protocols

*This product is available for passage.

[Thawing Protocol]-for 100mm plastic culture dish

Preparation

- 100mm dish for culturing cells
- Medium for the rat Schwann cell line (IFRS1)

1. Thaw vials of frozen cells by warming them in 37°C warm water for 2 minutes.
2. Transfer the thawed cell solution to a 50 mL centrifuge tube containing 10 mL of medium previously returned to room temperature, mix, and collect the cell solution by 1 mL separating the culture medium in the centrifuge tube and co-washing the vials.
3. Transfer cell suspension to 100mm dish (seeding density: approximately 0.9×10^4 cells/cm²). Incubate the plate at 37°C under 5% CO₂

Note: Centrifugation of cells after thawing are not recommended since these actions are more harmful to the cells than the effect of residual cell freezing medium in the culture and cell attachment.

4. Replace the medium with medium warmed to 37°C once the day after seeding, and then 2 to 3 times a week.

*After 3~5 days, IFRS1 cells reach approximately 80~100% confluent.

[Cell passage Protocol]

Preparation

- 70-90% confluent cells
- HBSS(-) or PBS(-)
- 0.05% Trypsin-EDTA
- Medium for the rat Schwann cell line (IFRS1)

1. Remove 70-90% confluent cells from the CO₂ incubator.
2. Aspirate the supernatant, add 10mL/100mm dish of HBSS(-) or PBS(-), and wash the dish.
3. Aspirate HBSS(-) or PBS(-) and add 1 mL/100 mm dish of 0.05% Trypsin-EDTA.
4. Store in a 37°C CO₂ incubator for about 3 to 4 minutes.
*Since trypsin damages cells, observe the degree of detachment under a microscope, and immediately proceed to the next treatment once almost all the cells have detached.
5. Add 10mL of medium warmed to 37°C
6. Pipette gently and collect the cell suspension in a 50mL centrifuge tube.
7. Centrifuge the cell suspension at 200 ×g for 5 minutes at 4°C, centrifuge, and aspirate the supernatant.
8. Add 10ml of media and gently pipet it to suspend the cell again.
9. Mix the cell suspension and medium in a ratio of 1:3 to 1:4 and seed it in a new culture vessel.
10. 1Incubate in a 37°C incubator in the presence of 5% CO₂. Passage should be performed in the same manner with cells semi-confluent (see Figure 1, C).

[Freezing Cells Protocol]

1. Cell passage protocols 1 to 6 are performed.
2. Centrifuge at 4°C for 5 minutes at 200 ×g, remove the supernatant, add COS banker (KOJ, Cat No. COS-CFM01) to 1×10⁶ cells/mL and suspend.
3. Dispense 1 mL/tube into cryopreservation tubes, and cryopreserve at -80°C by using cell freezing containers (BM-instrument, Cat No. BCS-136 or equivalent products).
4. Store at -80°C for 1 day or more. Transfer to liquid nitrogen.

VIII. Technical information


Figure 1. Immunofluorescence of IFRS1 cells.

A: Anti-p75NTR (green), nuclear staining (blue)

B: Anti-s100 (green), nuclear staining (blue)

C: phase-contrast microscopic image

Magnification : x 100(A,B,C)

IX. References

- 1) Watabe K, Fukuda T, Tanaka J, et al. Spontaneously immortalized adult mouse Schwann cells secrete autocrine and paracrine growth-promoting activities. *J Neurosci Res.* 1995;41:279-290.
- 2) Watabe K, Sakamoto T, Kawazoe Y, Michikawa M, Miyamoto K, Yamamura T, Saya H, Araki N. Tissue culture methods to study neurological disorders: Establishment of immortalized Schwann cells from murine disease models. *Neuropathology* 2003;23:64-74.

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