Rat Visceral Adipocyte H-1

Catalog No.PMC-VACH1C

April 17, 2020

Cell specification

Mesenteric adipocytes, a type of visceral adipocytes, are located along the portal vein that transports nutrients absorbed from the intestinal tract to the liver (Fig. 1). Evidence has shown that excess fat accumulation in the visceral adipose tissue contributes to the pathogenesis of Type II diabetes, hypertension and atherosclerosis. The Rat Visceral Adipocyte Culture Kit contains preadipocytes isolated from rat mesentery and culture medium that induces differentiation of precursor cells into mature adipocytes, finally causes hypertrophy. The kit provides a convenient system for studying the mechanism of adipogenesis as well as for screening drugs that prevent metabolic syndrome such as obesity, diabetes and hypertension by blocking the processes of adipogenesis.


Fig.1

Components

Product Name	Quantity	Amount	Storage Conditions	Stability
Visceral Preadipocytes, Rat	1.5 x 10 ⁶ cells/vial	1 vial	Liquid nitrogen storage	6months

※ Store frozen cells in liquid nitrogen if not ready for use after receipt.

Optimized Culture Medium (selling separately)

Catalog No.	Product Name	Size	Storage Conditions	Expiration date
VACMR	Visceral Adipocyte Culture Medium Ver.1	250 mL	-20°C Freezer	- Written on the bottle
			4°C	3months after thawing
VACM2	Visceral Adipocyte Culture Medium Ver.2	250 mL	-20°C Freezer	- Written on the bottle
			4°C	3months after thawing

Culture Medium components: DMEM/F-12, FBS, antibiotics, etc.

※Medium component comparison

Medium component	Visceral Adipocyte Culture Medium Ver.1 (#VACMR)	Visceral Adipocyte Culture Medium Ver.2 (#VACM2)
Differentiation inducer	×	×
fatty acid	○	○
Serum	○	○
Insulin concentration	10,000 ng/mL	1 ng/mL
Growth factor	×	○

Materials required but not provided

- Variable volume pipettes
- Culture plate, 24-well, flat bottom
- Visceral Adipocyte Culture Medium





Precautions

- Read the instructions carefully before beginning the culture.
- This kit is for research use only, not for human or diagnostic use.
- Always wear gloves and lab coat when handling the cell culture.

Protocols

1. Thaw the Visceral Adipocyte Culture Medium at 4°C
2. Thaw vials of frozen cells by warming them in 37°C warm water for 2 minutes.
3. Transfer thawed cell solution to a 15 mL centrifuge tube containing 10 mL of Visceral Adipocyte Culture Medium. Mix gently and centrifuge at 4°C at 200 x g for 5 minutes.
4. Remove the supernatant and re-suspend cells in 10 mL of Visceral Adipocyte Culture Medium and centrifuge at 4°C at 200 x g for 5 minutes.
5. Remove the supernatant and re-suspend the cell pellet in 6.3 mL of Visceral Adipocyte Culture Medium
6. Dispense 0.5 mL of cell suspension to each well of 24-well plate and incubate the plate at 37°C in a 5% CO₂ humidified incubator.
7. Gently add 0.5 ml of Visceral Adipocyte Culture Medium to each well of 24-well plate on the following day.
8. Carefully change the medium 1 ml per well on the second day after seeding.
9. Change the medium every other day. Be careful to not disturb the cell layer.
 - i. Approximately 3-4 days of culture, preadipocyte culture becomes confluent.
 - ii. Approximately 5 days of culture, cells turn to mature adipocytes.
 - iii. Approximately 8 days of culture, cells become hypertrophic and start detaching from the bottom of the well.

Technical information

COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN

Phone: +81-3-5632-9610

FAX: +81-3-5632-9619

URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA

email: info@cosmobiousa.com

Phone/FAX: (+1) 760-431-4600

URL: www.cosmobiousa.com

For research use only. Not for clinical diagnosis.

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(Rat Visceral Adipocyte, Code No. VACH1C)

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