

TECHNICAL DATA SHEET 576

Page 1 of 1

BioMag Multi-6 Microcentrifuge Tube Separator

Description

The Multi-6 is a durable, plastic magnetic separation unit designed to accommodate six 1.5ml microcentrifuge tubes. The magnets and tube holders are positioned to efficiently magnetically separate microliter volumes common in molecular biology applications.

Use

Magnetic separations can be performed on six 1.5ml microcentrifuge tubes simultaneously. Once the tubes are placed in the tube holders, the magnetic particles are pulled to the side of the unit to complete a separation in less than one minute. Upon completion of the separation, the supernatant may be aspirated, leaving the magnetic particles as a pellet on the side of the tube.

Care

The BioMag Multi-6 unit should be wiped clean with mild detergent.

DO NOT IMMERSE. Harsh chemicals and organic solvents may cause damage to the surface of the unit.

Warning

Strong Magnetic Field. This separation unit should be used and stored at a safe distance from magnetic media such as computers, videotapes and credit cards etc.

Ordering Information:

Cat. #	Description	Size
8MB4111S	BioMag Multi-6 Microcentrifuge Tube Separator	1 ea

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com

