ProSpecT[®] *Cryptosporidium* Microplate Assay

24/96 Tests

INTENDED USE: ProSpecT[®] *Cryptosporidium* Microplate Assay uses monoclonal antibodies for the qualitative detection of *Cryptosporidium* Specific Antigen (CSA) in aqueous extracts of fecal specimens.

SUMMARY AND EXPLANATION OF THE TEST: Cryptosporidiosis has recently been recognized as an important human disease in most areas of the world (6, 9). The causative agent, *Cryptosporidium* spp., has been identified in stool specimens of children and adults in many countries and most states in the United States (6). This parasite has been implicated in severe disease in HIV-infected persons (10), day-care centers (1, 2, 11, 12) and waterborne outbreaks (5, 7, 8) in the United States. Groups at particular risk include immunocompromised persons, especially those with HIV infection, family members and sexual partners of infected patients, children and caretakers in child day care centers, animal handlers and travelers (6). Acute symptoms of cryptosporidiosis may include diarrhea, abdominal pain, nausea and vomiting, fever, malaise, and respiratory problems lasting from several days to more than a month and often leading to persistent infection or death in immunologically deficient patients (6). Infection with *Cryptosporidium* may also be asymptomatic.

Cryptosporidium specific antigens have been found associated with *Cryptosporidium* infections and have been used as the basis of fluorescent and antigen capture immunoassays (3, 4, 13). Remel has identified a *Cryptosporidium* specific antigen (CSA) that is produced by *Cryptosporidium* organisms as they multiply within the host intestinal tract. The antigen is specific to *Cryptosporidium* and has not been found to cross-react with other enteric parasites. The antigen is stable to transport through the host intestinal tract as well as to routine procedures used to collect and transport specimens for microscopic examination.

PRINCIPLES OF THE PROCEDURE: ProSpecT[®] *Cryptosporidium* Microplate Assay is a solid phase immunoassay for the detection of CSA. Diluted stool specimens are added to break-away microplate wells on which anti-CSA antibody is bound. If CSA is present, it is "captured" by the bound antibody. The wells are incubated and then washed to remove unbound material. The enzyme conjugate (monoclonal anti-CSA antibody labeled with horseradish peroxidase enzyme) is added. The wells are incubated and then washed to remove unbound conjugate. In a positive reaction, CSA binds the enzyme to the well. The substrate for the enzyme, TMB, is added. In a positive reaction, the enzyme bound to the well by CSA converts the substrate to a colored reaction product. Color development can be detected visually or spectrophotometrically. In a negative reaction, there is no CSA or an insufficient level of CSA present to bind the enzyme conjugate and no colored reaction product develops.

REAGENTS: The ProSpecT[®] *Cryptosporidium* Microplate Assay includes sufficient reagents to perform 24 or 96 tests.

| | |
|---|----------------------|
| Reagents | 24 / 96 Test |
| Microplate* | 8 wells/strip |
| Coated with rabbit anti-CSA antibodies | 3 strips/12 strips |
| Enzyme Conjugate* | 1 bottle |
| Peroxidase labeled mouse monoclonal anti-CSA with bovine serum and 0.01% thimerosal | 5 ml/25 ml |
| Positive Control | 1 bottle |
| <i>Cryptosporidium</i> Oocyst Extract | 4 ml |
| Negative Control | 1 bottle |
| Human fecal material with 0.02% thimerosal | 4 ml |
| Specimen Dilution Buffer | 1 bottle |
| Buffered solution with rabbit serum and 0.02% thimerosal | 35 ml/110 ml |
| Wash Buffer | 1 bottle |
| 10X concentrated buffer solution with 0.1% thimerosal | 50 ml/110 ml |
| Color Substrate | 1 bottle |
| TMB in buffer | 5 ml/25 ml |
| Stop Solution | 1 bottle |
| 1.0 N Sulfuric acid (corrosive) | 6 ml |

*Note: Do not interchange reagents between kits with different lot numbers.

SYMBOL DEFINITIONS:

IVD For IN VITRO Diagnostic Use **LOT** Lot **REF** Catalog No.

 Store at 2-8°C  Expiration date

WARNINGS AND PRECAUTIONS:

- For IN VITRO Diagnostic Use Only.
- Reagents are provided at the necessary working strength, with the exception of the Wash Buffer concentrate. Do not dilute reagents, except where instructed.
 - Do not use reagents beyond the expiration dates. Expiration dates are printed on each reagent label. Use of reagents beyond the expiration date may affect the accuracy of results.

3. Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes when removing aliquots from reagent bottles.

4. Reagents are prepared from biological materials and should be handled as potentially infectious. Discard using appropriate biohazard procedures.

5. Microplate strips must be stored in the resealable foil pouch to protect microplate wells from moisture.

6. Specimens may contain potentially infectious agents and should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual, "Biosafety in Microbiological and Biomedical Laboratories, 4th Edition".

7. Stool samples must be thoroughly mixed prior to specimen processing to ensure accurate representation of the specimen. DO NOT CONCENTRATE SPECIMENS BEFORE TESTING.

8. Color substrate is sensitive to light exposure. If the reagent is exposed to light and develops color, the reagent must be discarded.

9. Persons who are color blind or visually impaired may not be able to read the test visually and should use spectrophotometric readings to interpret results.

10. Wash Buffer, Specimen Dilution Buffer, Enzyme Conjugate, Negative Control contain thimerosal, which may be irritating to skin, eyes and mucous membranes. In case of contact, flush eyes or skin with copious amounts of water.

11. Discard used Wash Buffer in appropriate biohazard containers.

12. The Wash Buffer (0.1%), Specimen Dilution Buffer (0.02%), Enzyme Conjugate (0.01%) and Negative Control (0.02%) contain thimerosal, which is classified per applicable European Economic Community (EEC) Directives as Dangerous for the environment (N). The following are the appropriate Risk (R) and Safety (S) phrases.

N  R53
S35
May cause long-term adverse effects in the aquatic environment.
This material and its container must be disposed of in a safe way.

13. Stop solution is corrosive and should be handled with care. If this solution comes into contact with skin or clothing, flush area with water for 15 minutes. If eyes are exposed, obtain additional medical attention.

14. The Stop Solution contains 2.8% Sulfuric acid, which is classified per applicable European Economic Community (EEC) Directives as an irritant (Xi). The following are the appropriate Risk (R) and Safety (S) phrases.

Xi  R36
R38
S26
Irritating to eyes
Irritating to skin
In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

PREPARATION OF REAGENTS: Before use, bring all reagents to room temperature (20 - 25°C) and mix gently. Return the unused reagents to the refrigerator after use.

All reagents, except the Wash Buffer, are supplied in ready-to-use bottles. Reagents can be dispensed directly from the dropper bottles or poured out for use with multichannel pipettes. If excess reagent has been poured, the excess should be discarded. Do not pour excess reagent back into the bottle.

Dilute 10x Wash Buffer concentrate to 1x by adding 1 part concentrate to 9 parts distilled or deionized water. Diluted Wash Buffer is stable for 1 month when stored at 2 - 8°C.

STORAGE AND HANDLING: The expiration date of each kit is stated on the package label. Store all components at 2 - 8°C. Unused microplate strips should be stored in the foil pouch containing desiccant to exclude moisture. The Color Substrate should be stored in and used from the light protected bottle in which it is provided. If an aliquot is removed from the original bottle for any reason, do not return unused Color Substrate to the original bottle.

COLLECTION OF STOOL SPECIMENS: Specimens collected for routine ova and parasite examination can be used for the ProSpecT[®] *Cryptosporidium* Microplate Assay. Stool specimens should be collected in clean, leak-proof plastic containers. Fresh, untreated stool specimens should be stored at 2 - 8°C and tested within 48 hours. If fresh specimens cannot be tested within 48 hours, they should be frozen at -20 to -70°C. Stool specimens treated with 10% formalin, MF or SAF fixatives may be refrigerated (2 - 8°C) or stored at room temperature (20 - 25°C) and should be tested within 2 months after collection. Stool specimens collected in Cary Blair Transport Medium (or equivalent) should be refrigerated or frozen and tested within 1 week after collection. Stool specimens that have been concentrated or treated with PVA fixatives are not suitable for use. Stool specimens obtained from rectal swabs and diapers are acceptable for use in the ProSpecT[®] *Cryptosporidium* Microplate Assay.

PROCEDURE NOTES

- Carefully read and follow all instructions in this package insert.
- Allow all reagents and specimens to reach room temperature (20 - 25°C) before use.
- Add reagents to the test wells in the same order throughout the procedure. To avoid contamination do not touch the fluid in the wells with the bottle tips.
- Time each incubation accurately. Start timing after adding reagent to the last well on each plate being tested. To ensure accurate timing, process no more than three 96 well plates at one time. Deviation from the established procedure may alter the performance of the assay.

TEST PROCEDURE :

Required materials provided:

Reagents Transfer pipettes
Microplate Strip Holder & Cover Color Reaction Procedure Card

Required materials not provided:

Stool specimen collection containers Timer that measures minutes
Wash bottle for Wash Buffer Distilled or deionized water

Optional materials not provided:

Microplate reader capable of reading 450nm
Cotton or rayon tipped applicator sticks
Micropipette to deliver volumes to 200µl
Plastic or glass disposable test tubes
Vortex mixer with plate adapter or shaker

Procedure:

- Open the foil pouch, remove the required number of microplate strips and place into a microplate strip holder. Use one well for the Negative Control and one well for the Positive Control. If using less than 8 wells, break off the required number of wells from a strip and return the unused wells to the foil pouch. RESEAL POUCH TIGHTLY TO EXCLUDE MOISTURE AND RETURN TO THE REFRIGERATOR.
- Specimens can be added directly into the wells or pre-diluted in tubes before adding to the wells. Pre-diluted specimens can be held at room temperature (20 - 25°C) for 8 hours or at 2-8°C for 48 hours prior to testing (see below). Choose one of these two methods: See Box **A** for dilution in wells; See Box **B** for pre-dilution in tubes.

A DILUTION IN WELLS

- Unpreserved Solid Specimens:** Label one tube for each specimen. Add **0.4 ml** Specimen Dilution Buffer (SDB) to each tube. Coat **1 swab** with specimen and vigorously mix into SDB. Express as much fluid as possible and discard the swab. Put a transfer pipette into the tube.
- Preserved or Watery Unpreserved Specimens:** Mix by shaking specimen collection containers. No further preparation is necessary.
- Add **4 drops** Negative Control to well A1. Add **4 drops** Positive Control to well B1.
- Add **100 µl** SDB to each specimen well.
- Using transfer pipettes add **2 drops** of each specimen to a well.
Note: Place the opening of the transfer pipettes just inside the wells to avoid splashing into adjacent wells.

PROCEED TO STEP 8

B DILUTION IN TUBES

- Unpreserved Solid Specimens:** Label one tube for each specimen. Add **1 ml** Specimen Dilution Buffer (SDB) to each tube. Coat **1 swab** with specimen and vigorously stir into SDB. Express as much fluid as possible and discard the swab. Put a transfer pipette into each tube.
- Preserved or Watery Unpreserved Specimens:** Label one tube for each specimen. Add **1 ml** SDB to each tube. Mix samples by shaking specimen collection containers. Using transfer pipettes draw up **0.3 ml** (third mark from the tip of the pipette). Expel sample into SDB. Mix by drawing up and down once. Leave transfer pipettes in the tubes.

Diluted specimens may be held for 8 hours at room temperature (20 - 25°C) or 48 hours at 2-8°C.

- Add **4 drops** Negative Control to well A1.
- Add **4 drops** Positive Control to well B1.
- Using transfer pipettes add **0.2 ml** (second mark from the tip of the pipette) of each specimen to a well.

Note: Place the opening of the transfer pipettes just inside the wells to avoid splashing into adjacent wells.

PROCEED TO STEP 8

- Cover the plate and incubate at room temperature (20 - 25°C) for **60 minutes**. Begin timing after the addition of the last specimen.
- Shake out or aspirate the contents of the wells. Wash by completely filling each well with **diluted** Wash Buffer (~350-400 µl/well). Shake out or aspirate all fluid from wells after each wash. Wash a total of **3 times**. After the last wash remove contents and strike plate on clean paper towels or aspirate. Remove as much Wash Buffer as possible but do not allow the wells to dry out at any time.
- Add **4 drops** (200 µl) of Enzyme Conjugate to each well.
- Cover the plate and incubate at room temperature (20 - 25°C) for **30 minutes**.
- Shake out or aspirate and wash each well **5 times** as in step 9.
- Add **4 drops** (200 µl) of Color Substrate to each well.
- Cover the plate and incubate at room temperature (20 - 25°C) for **10 minutes**.
- Add **1 drop** (50 µl) of Stop Solution to each well. Gently tap or vortex the wells until the yellow color is uniform. Read reactions within **10 minutes** after adding stop solution.
- Read visually or spectrophotometrically at 450 nm.

QUALITY CONTROL: Positive and Negative Controls must be included each time the test is performed. The Positive and Negative Controls serve as both reagent and procedural controls. The controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cutoff.

The optical density (O.D.) of the Negative Control should be ≤0.100 at 450nm and should be colorless when read visually. If yellow color equal to 1+ or greater on the Procedure Card is present in the Negative Control, the test should be repeated with careful attention to the wash procedure.

The O.D. of the Positive Control should be ≥0.300 at 450nm, after the O.D. of the Negative Control is subtracted and should be equal to or greater than the 2+ reaction when read visually. If yellow color less than 2+ on the Procedure Card is present in the Positive Control, call for technical assistance.

RESULTS: Refer to enclosed Procedure Card for color interpretations.**Visual:**

1. Read the test results by comparing with the reaction colors on the procedure card.

Positive: yellow color of at least 1+ intensity

Negative: colorless

2. Interpretation of test results:

Positive: If yellow color of at least 1+ intensity develops in the test well, the sample contains CSA and the test is positive.

Note: Tests with faint yellow color (less than 1+) should be repeated.

Negative: A colorless reaction is a negative result and indicates that no CSA or an undetectable level of CSA is present in the sample tested.

Spectrophotometric:

1. Set the spectrophotometer (microplate reader) to read at 450 nm.

2. Read the optical density (O.D.) of the Negative Control.

3. Subtract the O.D. of the Negative Control well from the O.D. readings of the Positive Control well and the test wells before interpreting results.

Note: Readers may be set to blank on the Negative Control well so that the Negative Control well O.D. is automatically subtracted from all of the other readings. If the reader does not have this capability, blank on air and subtract the value of the Negative Control well from the values of the Positive Control and test wells before interpreting results.

4. Read the test results:

Positive: O.D. of ≥ 0.050 blanked value

(i.e. after the O.D. of the Negative Control is subtracted)

Negative: O.D. of < 0.050 blanked value

(i.e. after the O.D. of the Negative Control is subtracted)

5. Interpretation of spectrophotometric results:

Positive: If the blanked O.D. reading is equal to or greater than 0.050 in the test well, the sample contains CSA and the test is positive.

Negative: A blanked O.D. reading less than 0.050 is a negative result and indicates that no CSA or an undetectable level of CSA is present in the sample tested.

***Note:** Any wells that are visually clear but give an O.D. reading that is inconsistent with the visual interpretation should be considered a discrepant reading and examined for the presence of bubbles, small particles in the wells, or an opaque film on the bottom of the well. To remove these, wipe the underside of the wells and read the O.D. again. If the discrepancy between visual and O.D. readings persists, repeat the test.

Limitations of the Procedure: The validity of results with the ProSpecT[®] *Cryptosporidium* Microplate Assay depends on the control reaction performing as expected. See Quality Control section.

A negative test result does not exclude the possibility of the presence of *Cryptosporidium*, and may occur when the antigen level in the sample is below the detection level of the test. Correlation between the amount of antigen in a sample and clinical presentation has not been established.

As with all IN VITRO diagnostic tests, results should be interpreted by the clinician in conjunction with clinical findings and/or other laboratory results.

Proper specimen collection and processing are essential to achieve optimal performance of the assay. Optimal test results are obtained from specimens tested as soon after collection as possible. See Collection of Stool Specimens section.

ProSpecT[®] *Cryptosporidium* Microplate Assay has been classified high complexity.

EXPECTED VALUES: The prevalence of *Cryptosporidium* infection varies in different populations and geographic areas. In the U.S., the incidence of *Cryptosporidium* is approximately 0.5 - 3.0% with higher prevalence rates in children (12) and in homosexual males (5, 6).

PERFORMANCE CHARACTERISTICS:

Sensitivity and Specificity: Studies were conducted to evaluate the performance of the ProSpecT[®] *Cryptosporidium* Microplate Assay with specimens obtained from large reference and hospital laboratories which performed O&P testing. A total of 214 specimens were tested; 81 were positive for *Cryptosporidium* by acid-fast (AF) and 133 were negative. Forty of the *Cryptosporidium* negative specimens contained organisms other than *Cryptosporidium* by AF or O&P. The results of these evaluations are presented below:

| | | Acid Fast | | |
|---|---|-----------|-----|-----|
| | | + | - | |
| ProSpecT [®] <i>Cryptosporidium</i> | + | 81 | 0 | 214 |
| | - | 0 | 133 | |
| | | 81 | 133 | |

Sensitivity 81/81 = 100% (95.5–100%)

Specificity 133/133 = 100% (97.3–100%)

Numbers in parenthesis are 95% confidence intervals.

Clinical studies were conducted to evaluate the performance of the ProSpecT[®] *Cryptosporidium* Microplate Assay. Specimens were obtained from hospital labs and the CDC. Patient populations represented in the specimen pool were symptomatic patients in normal prevalence populations, symptomatic patients in a high prevalence (HIV positive) population and asymptomatic patients from a day care population. Specimens were submitted unpreserved or preserved in 10% formalin or SAF. Samples were tested for *Cryptosporidium* by either acid-fast (AF) or immunofluorescent staining methods (IFA). A total of 212 specimens were tested; 134 were positive for *Cryptosporidium* specific antigen (CSA) and 78 were

negative. The results with the ProSpecT[®] *Cryptosporidium* Microplate Assay are presented below:

| | | Acid Fast | | |
|---|---|-----------|----|-----|
| | | + | - | |
| ProSpecT [®] <i>Cryptosporidium</i> | + | 130 | 0 | 212 |
| | - | 4 | 78 | |
| | | 134 | 78 | |

Sensitivity 130/134 = 97% (92.5–99.2%)

Specificity 78/78 = 100% (95.4–100%)

Numbers in parenthesis are 95% confidence intervals.

A prospective trial was conducted at a large metropolitan hospital. All samples submitted for acid fast staining for *Cryptosporidium* over a period of 4 months were included in the study. Samples were unpreserved and frozen at -20°C prior to testing with the ProSpecT[®] *Cryptosporidium* Microplate Assay. The results of initial testing and the resolved data are presented below. Data was resolved by repeat testing of the 14 AF negative/CSA positive samples. Six of the fourteen were reproducibly positive for CSA. Specific inhibition studies with antibody to CSA showed greater than 50% inhibition in all 6 samples. These 6 samples are considered to be true positives in the resolved data.

| | | Acid Fast | | Resolved | | |
|---|---|-----------|-----|----------|-----|-----|
| | | + | - | + | - | |
| ProSpecT [®] <i>Cryptosporidium</i> | + | 28 | 14 | 34 | 8 | 378 |
| | - | 1 | 335 | 1 | 335 | |
| | | 29 | 349 | 35 | 343 | |

Sensitivity 28/29 = 97% (82.2–99.9%) 34/35 = 97% (85.1–99.9%)

Specificity 335/349 = 96% (93.4–97.8%) 335/343 = 98% (95.5–99.0%)

Numbers in parentheses are 95% confidence intervals.

Analytical Sensitivity: The ProSpecT[®] *Cryptosporidium* Microplate Assay detects approximately 20 nanograms/ml of CSA.

Reproducibility: The inter-assay coefficient of variation of the ProSpecT[®] *Cryptosporidium* Microplate Assay was evaluated by selecting 10 positive specimens with varying optical density readings. Each sample was tested in 10 wells per day for five days. The mean inter-assay CV was 10.6%.

The intra-assay CV was evaluated by testing 24 wells with each of 5 positive specimens. The mean intra-assay CV was 2.52%.

Cross-Reactivity: The ProSpecT[®] *Cryptosporidium* Microplate Assay has been tested with stool specimens found to be O&P positive for a number of fecal parasites. No cross-reactivity was observed with any of the infectious agents listed below.

| | | |
|---------------------------------|----------------------------------|--------------------------------------|
| <i>Ascaris lumbricoides</i> (2) | <i>Entamoeba histolytica</i> (5) | <i>Blastocystis hominis</i> (4) |
| <i>Hymenolepis nana</i> (2) | <i>Endolimax nana</i> (3) | <i>Iodamoeba butschlii</i> (2) |
| <i>Dientamoeba fragilis</i> (4) | <i>Entamoeba hartmanni</i> (2) | <i>Isospora belli</i> (2) |
| <i>Taenia solium</i> (1) | <i>Chilomastix mesnili</i> (1) | <i>Strongyloides stercoralis</i> (2) |
| <i>Entamoeba coli</i> (6) | <i>Giardia lamblia</i> (5) | <i>Trichuris trichiura</i> (1) |

The numbers in parentheses indicate the number of specimens tested.

BIBLIOGRAPHY

- Alpert, G. et al. 1986. "Outbreak of Cryptosporidiosis in a Day-Care Center." *Pediatrics* 77(2):152-157.
- Anon. 1984. "Cryptosporidiosis among children attending day-care centers - Georgia, Pennsylvania, Michigan, California, New Mexico." *Morbidity and Mortality Weekly Report* 33(42):599-601.
- Arrowood, M.J. and C.R. Sterling. 1989. "Comparison of Conventional Staining Methods and Monoclonal Antibody-based Methods for *Cryptosporidium* Oocyst Detection." *J. Clin. Microbiol.* 27(7):1490-1495.
- Chapman, P.A., B.A. Rush and J. McLaughlin. 1990. "An enzyme immunoassay for detecting *Cryptosporidium* in faecal and environmental samples." *J. Med. Microbiol.* 32:233-237.
- D'Antonio, R.G. et al. 1985. "A waterborne outbreak of cryptosporidiosis in normal hosts." *Ann. Intern. Med.* 103:886.
- Dubey, J.P., C.A. Speer and R. Fayer, eds. 1990. "Cryptosporidiosis of Man and Animals." CRC Press.
- Gallaher, M.M. et al. 1989. "Cryptosporidiosis and Surface Water." *Am. J. Public Health* 79(1):39-42.
- Hayes, E.B. et al. 1989. "Large Community Outbreak of Cryptosporidiosis Due to Contamination of a Filtered Public Water Supply." *N. Eng. J. Med.* 320(21):1372-1376.
- Leech, J.H., M.A. Sande and R.K. Root, eds. 1988. "Parasitic Infections." Churchill Livingstone.
- Navin, T.R. and A.M. Hardy. 1984. "Cryptosporidiosis in Patients with AIDS." *J. Infect. Dis.* 155:150.
- Stehr-Green, J.K. et al. 1987. "Shedding of Oocysts in Immunocompetent Individuals Infected with *Cryptosporidium*." *Am. J. Trop. Med. Hyg.* 36(2):338-342.
- Taylor, J.P. et al. 1985. "Cryptosporidiosis Outbreak in a Day-Care Center." *Am. J. Dis. Child.* 139:1023-1025.
- Ungar, B.L.P. 1990. Enzyme-linked immunoassay for detection of *Cryptosporidium* antigens in fecal specimens. *J. Clin. Microbiol.*, 28:2491.

ORDERING

ProSpecT[®] *Cryptosporidium* Microplate Assay 24 Test Kit ...REF 2454024
ProSpecT[®] *Cryptosporidium* Microplate Assay 96 Test Kit ...REF 2454096

For ordering or customer and technical assistance please call: Remel Inc.

Manufactured by:

Remel Inc.

12076 Santa Fe Drive, Lenexa, KS 66215, USA

General Information **800-255-6730**

Technical Services **800-447-3641**

Order Entry **800-447-3635**

Local / International **913-888-0939**

Local / International Fax **913-888-5884**

Website: www.remel.com Email: remel@remel.com

**Authorized Representative in the European Community:**

Remel Europe Ltd.

Clipper Boulevard West, Crossways

Dartford, Kent, DA2 6PT

UK

IFU 24540GB01, Revised January 7, 2005

Printed in the USA

ProSpecT[®] is a registered trademark of Remel Inc.

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- Alpert, G. et al. 1986. "Outbreak of Cryptosporidiosis in a Day-Care Center". *Pediatrics* 77(2):152-157.
- Anon. 1984. "Cryptosporidiosis among children attending day-care centers - Georgia, Pennsylvania, Michigan, California, New Mexico". *Morbidity and Mortality Weekly Report* 33(42):599-601.
- Arrowood, M.J. and C.R. Sterling. 1989. "Comparison of Conventional Staining Methods and Monoclonal Antibody-based Methods for *Cryptosporidium* Oocyst Detection". *J. Clin. Microbiol.* 27(7):1490-1495.
- Chapman, P.A., B.A. Rush and J. McLaughlin. 1990. "An enzyme immunoassay for detecting *Cryptosporidium* in faecal and environmental samples". *J. Med. Microbiol.* 32:233-237.
- D'Antonio, R.G. et al. 1985. "A waterborne outbreak of cryptosporidiosis in normal hosts". *Ann. Intern. Med.* 103:886.
- Dubey, J.P., C.A. Speer and R. Fayer, eds. 1990. "Cryptosporidiosis of Man and Animals". CRC Press.
- Gallagher, M.M. et al. 1989. "Cryptosporidiosis and Surface Water". *Am. J. Public Health* 79(1):39-42.
- Hayes, E.B. et al. 1989. "Large Community Outbreak of Cryptosporidiosis Due to Contamination of a Filtered Public Water Supply". *N. Eng. J. Med.* 320(21):1372-1376.
- Leech, J.H., M.A. Sande and R.K. Root, eds. 1988. "Parasitic Infections". Churchill Livingstone.
- Navin, T.R. and A.M. Hardy. 1984. "Cryptosporidiosis in Patients with AIDS". *J. Infect. Dis.* 155:150.
- Stehr-Green, J.K. et al. 1987. "Shedding of Oocysts in Immunocompetent Individuals Infected with *Cryptosporidium*". *Am. J. Trop. Med. Hyg.* 36(2):338-342.
- Taylor, J.P. et al. 1985. "Cryptosporidiosis Outbreak in a Day-Care Center". *Am. J. Dis. Child.* 139:1023-1025.
- Ungar, B.P. 1990. Enzyme-linked immunoassay for detection of *Cryptosporidium* antigens in fecal specimens. *J. Clin. Microbiol.*, 28:2491.

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BIBLIOGRAPHIE

- Alpert, G. et al. 1986. "Outbreak of Cryptosporidiosis in a Day-Care Center." *Pediatrics* 77(2):152-157.
- Anon. 1984. "Cryptosporidiosis among children attending day-care centers - Georgia, Pennsylvania, Michigan, California, New Mexico." *Morbidity and Mortality Weekly Report* 33(42):599-601.
- Arrowood, M.J. and C.R. Sterling. 1989. "Comparison of Conventional Staining Methods and Monoclonal Antibody-based Methods for *Cryptosporidium* Oocyst Detection." *J. Clin. Microbiol.* 27(7):1490-1495.
- Chapman, P.A., B.A. Rush and J. McLaughlin. 1990. "An enzyme immunoassay for detecting *Cryptosporidium* in faecal and environmental samples." *J. Med. Microbiol.* 32:233-237.
- D'Antonio, R.G. et al. 1985. "A waterborne outbreak of cryptosporidiosis in normal hosts." *Ann. Intern. Med.* 103:886.
- Dubey, J.P., C.A. Speer and R. Fayer, eds. 1990. "Cryptosporidiosis of Man and Animals." CRC Press.
- Gallaher, M.M. et al. 1989. "Cryptosporidiosis and Surface Water." *Am. J. Public Health* 79(1):39-42.
- Hayes, E.B. et al. 1989. "Large Community Outbreak of Cryptosporidiosis Due to Contamination of a Filtered Public Water Supply." *N. Eng. J. Med.* 320(21):1372-1376.
- Leech, J.H., M.A. Sande and R.K. Root, eds. 1988. "Parasitic Infections." Churchill Livingstone.
- Navin, T.R. and A.M. Hardy. 1984. "Cryptosporidiosis in Patients with AIDS." *J. Infect. Dis.* 155:150.
- Stehr-Green, J.K. et al. 1987. "Shedding of Oocysts in Immunocompetent Individuals Infected with *Cryptosporidium*." *Am. J. Trop. Med. Hyg.* 36(2):338-342.
- Taylor, J.P. et al. 1985. "Cryptosporidiosis Outbreak in a Day-Care Center." *Am. J. Dis. Child.* 139:1023-1025.
- Ungar, B.L. 1990. Enzyme-linked immunoassay for detection of *Cryptosporidium* antigens in fecal specimens. *J. Clin. Microbiol.*, 28:2491.

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- Alpert, G. et al. 1986. "Outbreak of Cryptosporidiosis in a Day-Care Center." Pediatrics 77(2):152-157.
- Anon. 1984. "Cryptosporidiosis among children attending day-care centers - Georgia, Pennsylvania, Michigan, California, New Mexico." Morbid. Mortal. Weekly Rep. 33(42):599-601.
- Arrowood, M.J. and C.R. Sterling. 1989. "Comparison of Conventional Staining Methods and Monoclonal Antibody-based Methods for *Cryptosporidium* Oocyst Detection." J. Clin. Microbiol. 27(7):1490-1495.
- Chapman, P.A., B.A. Rush and J. McLaughlin. 1990. "An enzyme immunoassay for detecting *Cryptosporidium* in faecal and environmental samples." J. Med. Microbiol. 32:233-237.
- D'Antonio, R.G. et al. 1985. "A waterborne outbreak of cryptosporidiosis in normal hosts." Ann. Intern. Med. 103:886.
- Dubey, J.P., C.A. Speer and R. Fayer, eds. 1990. "Cryptosporidiosis of Man and Animals." CRC Press.
- Gallaher, M.M. et al. 1989. "Cryptosporidiosis and Surface Water." Am. J. Public Health 79(1):39-42.
- Hayes, E.B. et al. 1989. "Large Community Outbreak of Cryptosporidiosis Due to Contamination of a Filtered Public Water Supply." N. Eng. J. Med. 320(21):1372-1376.
- Leech, J.H., M.A. Sande and R.K. Root, eds. 1988. "Parasitic Infections." Churchill Livingstone.
- Navin, T.R. and A.M. Hardy. 1984. "Cryptosporidiosis in Patients with AIDS." J. Infect. Dis. 155:150.
- Stehr-Green, J.K. et al. 1987. "Shedding of Oocysts in Immunocompetent Individuals Infected with *Cryptosporidium*." Am. J. Trop. Med. 1 Hyg 36(2):338-342.
- Taylor, J.P. et al. 1985. "Cryptosporidiosis Outbreak in a Day-Care Center." Am. J. Dis. Child. 139:1023-1025.
- Ungar, B.P. 1990. Enzyme-linked immunoassay for detection of *Cryptosporidium* antigens in fecal specimens. J. Clin. Microbiol., 28:2491.

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BIBLIOGRAFIA

1. Alpert, G. et al. 1986. "Outbreak of Cryptosporidiosis in a Day-Care Center." *Pediatrics* 77(2):152-157.
2. Anon. 1984. "Cryptosporidiosis among children attending day-care centers - Georgia, Pennsylvania, Michigan, California, New Mexico." *Morbidity and Mortality Weekly Report* 33(42):599-601.
3. Arrowood, M.J. and C.R. Sterling. 1989. "Comparison of Conventional Staining Methods and Monoclonal Antibody-based Methods for *Cryptosporidium* Oocyst Detection." *J. Clin. Microbiol.* 27(7):1490-1495.
4. Chapman, P.A., B.A. Rush and J. McLaughlin. 1990. "An enzyme immunoassay for detecting *Cryptosporidium* in faecal and environmental samples." *J. Med. Microbiol.* 32:233-237.
5. D'Antonio, R.G. et al. 1985. "A waterborne outbreak of cryptosporidiosis in normal hosts." *Ann. Intern. Med.* 103:886.
6. Dubey, J.P., C.A. Speer and R. Fayer, eds. 1990. "Cryptosporidiosis of Man and Animals." CRC Press.
7. Gallaher, M.M. et al. 1989. "Cryptosporidiosis and Surface Water." *Am. J. Public Health* 79(1):39-42.
8. Hayes, E.B. et al. 1989. "Large Community Outbreak of Cryptosporidiosis Due to Contamination of a Filtered Public Water Supply." *N. Eng. J. Med.* 320(21):1372-1376.
9. Leech, J.H., M.A. Sande and R.K. Root, eds. 1988. "Parasitic Infections." Churchill Livingstone.
10. Navin, T.R. and A.M. Hardy. 1984. "Cryptosporidiosis in Patients with AIDS." *J. Infect. Dis.* 155:150.
11. Stehr-Green, J.K. et al. 1987. "Shedding of Oocysts in Immunocompetent Individuals Infected with *Cryptosporidium*." *Am. J. Trop. Med. Hyg.* 36(2):338-342.
12. Taylor, J.P. et al. 1985. "Cryptosporidiosis Outbreak in a Day-Care Center." *Am. J. Dis. Child.* 139:1023-1025.
13. Ungar, B.P. 1990. Enzyme-linked immunoassay for detection of *Cryptosporidium* antigens in fecal specimens. *J. Clin. Microbiol.*, 28:2491.

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Per effettuare ordini o richiedere assistenza tecnica fare riferimento a: Remel Inc.
Remel Inc. 12076 Santa Fe Drive, Lenexa, KS 66215, USA
ProSpec[™] è un marchio registrato di Remel Inc.