Cat. No. SGE-C-452
Anti Collagen IV α 5(IV) Chain, Human (Mono)
Clone: H52

Cat. No. SGE-C452

www.cosmobio.co.jp

Code: SGE-C452
Quantity: 500 μ l
Appearance: Solution. Monoclonal antibody was collected from supernatant by high-density culture and diluted with serum GIT medium.
(Containing 10% FBS).
0.1% NaN₃ was added for preservation.
Clone names: H52 (rat IgG2b/kappa)
Specificity of antibodies:
H52 is specific to NC1 domain of α 5(IV) (Reference 1).
Reactive with NC1 domain of α 5(IV) of rat, rabbit, guinea pig, dog, and bovine.
Preparation of antibodies:
Monoclonal antibody was prepared by the rat lymph node method developed by Shigei Medical Research Institute with synthetic peptide of type IV collagen as immunogen.
Use: Staining of human cryostat sections by indirect immunofluorescence (Acid treatment is necessary.); Western blotting (dilution x10)
Storage: Shipped at 4° C.
Store at 4° C short term (< 1 month). Store under -20 °C, be stable for 2 years.
Avoid freeze / thaw cycle.
Reference: 1) Sado et al. (1995) Establishment by the rat lymph node method of epitope-defined monoclonal antibodies recognizing the six different α chains of human type IV collagen. Histochem. Cell Biol. 104: 267-275.
Producer : Shigei Medical Research Institute

For research use only, Not for diagnostic use.
Please read this manual thoroughly before use.

Manufacturer

Shigei Medical Research Institute

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 COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com

Cat. No. SGE-C-452
Anti Collagen IV α 5(IV) Chain, Human (Mono)
Clone: H52

Cat. No. SGE-C452

www.cosmobio.co.jp

This product can be used for immunostaining of cryostat sections and for western blotting


How to stain human cryostat sections with C-452 (H52)

Materials

1. Acid solution : KCl-HCl buffer, pH 1.5, Ion strength=0.1 (Mix equal volumes of 0.1 M KCl and 0.1 M HCl.)
2. Blocking solution: 10% goat serum (normal goat serum diluted with PBS)
3. First antibody: H52 monoclonal Ab is usually diluted to 1/10 - 1/100 with PBS containing 1 mg/ml of bovine serum albumin before use.
4. Secondary Ab: Affinity purified Ab to rat IgG (Cappel #55760) is usually diluted to 1/50- 1/100 with PBS containing 5% human normal serum before use.

Procedure

1. Cryostat sections: 4 mm
2. Treatment with acetone for 10 min
3. Acid-solution treatment with 0.1 M KCl- HCl (pH1.5) for 10 min
4. Washing with PBS
5. Blocking with 10% normal goat serum for 10 min
6. Washing
7. First Ab for 60 min
8. Washing
9. Secondary Ab for 30 min
10. Washing
11. Mount with Perma Fluor (Tanner #11611)

Human cryostat section stained with monoclonal antibody H52. Glomerular basement membrane (BM), part of tubular BM, and Bowman's capsular BM are stained.

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COSMO BIO CO., LTD.