Cat. No. SGE-C-452
Anti Collagen IV α 5(IV) Chain, Human (Mono)

Clone: H52

Cat. No. SGE-C452

www.cosmobio.co.jp

Code: SGE-C452
Quantity: 500 μ l
Appearance: Solution. Monoclonal antibody was collected from supernatant by high-density culture and diluted with serum GIT medium.
(Containing 10% FBS).
0.1% NaN₃ was added for preservation.
Clone names: H52 (rat IgG2b/kappa)
Specificity of antibodies:
H52 is specific to NC1 domain of α 5(IV) (Reference 1).
Reactive with NC1 domain of α 5(IV) of rat, rabbit, guinea pig, dog, and bovine.
Preparation of antibodies:
Monoclonal antibody was prepared by the rat lymph node method developed by Shigei Medical Research Institute with synthetic peptide of type IV collagen as immunogen.
Use: Staining of human cryostat sections by indirect immunofluorescence (Acid treatment is necessary.); Western blotting (dilution x10)
Storage: Shipped at 4° C.
Store at 4° C short term (< 1 year). Store under -20 °C, be stable for 2 years.
Avoid freeze / thaw cycle.
Reference: 1) Sado et al. (1995) Establishment by the rat lymph node method of epitope-defined monoclonal antibodies recognizing the six different α chains of human type IV collagen. Histochem. Cell Biol. 104: 267-275.
Producer: Shigei Medical Research Institute

For research use only, Not for diagnostic use.
Please read this manual thoroughly before use.

Manufacturer

COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: support@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com

Cat. No. SGE-C-452
Anti Collagen IV α 5(IV) Chain, Human (Mono)
Clone: H52

Cat. No. SGE-C452

www.cosmobio.co.jp

This product can be used for immunostaining of cryostat sections and for western blotting


How to stain human cryostat sections with C-452 (H52)

Materials

1. Acid solution : KCl-HCl buffer, pH 1.5, Ion strength=0.1 (Mix equal volumes of 0.1 M KCl and 0.1 M HCl.)
2. Blocking solution: 10% goat serum (normal goat serum diluted with PBS)
3. First antibody: H52 monoclonal Ab is usually diluted to 1/10 - 1/100 with PBS containing 1 mg/ml of bovine serum albumin before use.
4. Secondary Ab: Affinity purified Ab to rat IgG (Cappel #55760) is usually diluted to 1/50- 1/100 with PBS containing 5% human normal serum before use.

Procedure

1. Cryostat sections: 4 mm
2. Treatment with acetone for 10 min
3. Acid-solution treatment with 0.1 M KCl- HCl (pH1.5) for 10 min
4. Washing with PBS
5. Blocking with 10% normal goat serum for 10 min
6. Washing
7. First Ab for 60 min
8. Washing
9. Secondary Ab for 30 min
10. Washing
11. Mount with Perma Fluor (Tanner #11611)

Human cryostat section stained with monoclonal antibody H52. Glomerular basement membrane (BM), part of tubular BM, and Bowman's capsular BM are stained.

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