

Cat # SL100668 Store at 4 °C

LipoD293™ DNA In Vitro Transfection Reagent

----- A Protocol for Transfections of Mammalian Cell

- ☐ 100 µl
- ☐ 500 µl
- ☐ 1000 µl



15875 Gaither Drive
Gaithersburg, MD 20877
FAX. 301-560-4919
TEL. 301-330-5966
Toll Free. 1-(866)-918-6812
Email: info@signagenlabs.com
Web: www.signagenlabs.com

This product is for laboratory research ONLY and not for diagnostic use

Introduction:

LipoD293™ DNA In Vitro Transfection Reagent is a powerful transfection Reagent that ensures effective and reproducible transfection with invisible cytotoxicity. LipoD293™ is a specially designed cationic lipid formulation which offers extremely high transfection efficiencies of HEK293 cells as well as many mammalian cells. LipoD293™ Reagent, 1.0 ml, is sufficient for 300 to 600 transfections in 24 well plates or 50 to 100 transfections in 6 well plates.

Features:

- Invisible cytotoxicity
- Exceptional transfection efficiency on HEK293 cells
- Efficient transfection with or without serum
- High levels of recombinant protein production
- Simple, robust transfection procedure

Procedures for Transfecting Mammalian Cells:

Cell Seeding (see Table 1):

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal 60-70% confluency at the time of transfection. Serum-free DMEM medium is changed to replace complete serum-containing culture medium 30 minutes before transfection.

Note: High serum levels (>5%) have a moderate inhibitory effect on LipoD293™-mediated transfections. Maximal transfection efficiencies are observed in the absence of serum. Depending upon the cell type, the presence of serum <5% may sometimes improve the overall levels of recombinant protein expression.

Table 1. A Guideline for Seeding Adherent Cells Prior to Transfection in Different Culture Formats

| Culture Dishes | Surface Area (cm ²) | Number of Cells to Seed |
|----------------|---------------------------------|---------------------------|
| T175 Flask | 175 | 0.7–1.4 × 10 ⁷ |
| T75 Flask | 75 | 3.0–6.0 × 10 ⁶ |
| 100 mm Dish | 58 | 2.2–4.4 × 10 ⁶ |
| 60 mm Dish | 21 | 0.9–1.8 × 10 ⁶ |
| 35 mm Dish | 9.6 | 3.5–7.0 × 10 ⁵ |
| 6-well Plate | 9.6 | 4.0–8.0 × 10 ⁵ |
| 12-well Plate | 3.5 | 1.5–3.0 × 10 ⁵ |
| 24-well Plate | 1.9 | 0.8–1.6 × 10 ⁵ |
| 48-well Plate | 1.0 | 4.0–8.0 × 10 ⁴ |
| 96-well Plate | 0.3 | 1.2–2.4 × 10 ⁴ |

Preparation of LipoD293™-DNA Complex and Transfection Procedures

The optimal ratio of LipoD293™/DNA is of 3/1. We recommend using serum-free DMEM with High Glucose to dilute DNA and LipoD293™ Reagent to ensure the optimal size of complex particles.

The following protocol is given for transfection in 24-well plates, refer to **Table 2** for transfection in other culture formats. The optimal transfection conditions for HEK293 cell are given in the standard protocol described below.

- For each well, dilute 1 µg of DNA into 50 µl of serum-free DMEM with High Glucose. Vortex gently and spin down briefly.
- For each well, dilute 3 µl of LipoD293™ solution into 50 µl of serum-free DMEM with High Glucose. Vortex gently and spin down briefly.
- Add the 50 µl LipoD293™ solution to the 50 µl DNA solution all at once. (**Important: do not mix the solutions in the reverse order !**)
- Vortex- mix the solution immediately and spin down briefly to bring drops to the bottom of the tube.
- Incubate for 10 minutes at room temperature.
- Add the 100 µl LipoD293™/ DNA mixture drop-wise onto the medium in each well and homogenize the mixture by gently swirling the plate.
- For maximal transfection efficiency, change the medium to complete serum containing medium 4~5 hours post addition of LipoD293™/DNA complex.
- Check transfection efficiency 24 to 48 hours post transfection.

Table 2. Recommended Amounts for Different Culture Vessel Formats

| Culture Dish | Culture Volume (ml) | Plasmid DNA (µg) | Diluent Volume (mL) | GenJet™ Reagent (µL) |
|--------------|---------------------|------------------|---------------------|----------------------|
| 6-well plate | 1.6 | 2 - 3 | 0.20 | 6 - 9 |
| 35 mm dish | 1.6 | 2 - 3 | 0.20 | 6 - 9 |
| 60 mm dish | 4.5 | 5 | 0.50 | 15 |
| 100 mm dish | 8 | 7 - 8 | 1.0 | 21 - 24 |
| T75 flask | 15 | 18 - 36 | 1.5 | 54 - 108 |
| 250 ml flask | 50 | 50 - 100 | 2.5 | 150 - 300 |

Storage: LipoD293™ DAN In Vitro Transfection Reagent is stable for up to 18 months at 4 °C. This item shipped at ambient temperature