SHIMA Laboratories CO., LTD.

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Citrullinated protein detection antibody

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
1. Senshu, T. et al., Detection of citrulline residues in deiminated proteins on polyvinylidene difluoride membrane. *Anal. Biochem.* 203, 94-100 (1992)
2. Senshu, T. et al., Detection of deiminated proteins in rat skin: probing with a monospecific antibody after modification of citrulline residues. *J. Invest. Dermatol.* 105, 163-169 (1995)
3. Ishigami, A. et al., Abnormal accumulation of citrullinated proteins catalyzed by peptidylarginine deiminase in hippocampal extracts from patients with Alzheimer's disease. *J. Neurosci. Res.* 80, 120-128 (2005)
4. 石神昭人: アルツハイマー病とシトルリン化蛋白質. 医学のあゆみ 255, 889-893 (2015)

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SHIMA Laboratories CO., LTD.

ROI004

**Citrullinated protein detection antibody (for Western blot analysis)**

**Another name: Anti-modified citrulline (AMC) IgG**

Rabbit Polyclonal Antibody (Purified IgG Fraction)

Rabbit AMC IgG

Volume : 0.075 mL

Antigen : Calf thymus histones (Histone H1 depleted) deiminated *in vivo* with recombinant human PAD2 proteins. The resulting deiminated histones were then modified by ureid group specific color reaction.

Host : Rabbit

Supplied As : IgG fraction purified from rabbit serum.

Prepared in 10 mM Tris (pH 7.4), 0.14 M NaCl and 0.05% NaN<sub>3</sub>.

Storage and Stability : 2-8 °C, 2 years.


Tested applications : • Western Blot (1:500-1:1,000 dilution)

To detect deiminated proteins by Western blot using this antibody, it is essential to perform chemical modification (chemical modification of citrulline residues) of PVDF membrane after transfer the electrophoresed proteins. The details refer to experiment protocol.

• We have no data about immunohistochemistry.

Cross Reactivity: Cross reacts with Human, Mouse and Rat deiminated proteins.

Not yet tested in other species.

**[Western blot analysis]**

Lane 1 : Histones (from calf thymus)

Lane 2 : Deiminated histones

Lane 3 : Human skin cornified cell lysate

A : Chemical modification of PVDF membrane

B : Control

Anti-modified citrulline (AMC) IgG (ROI004) at 1:1,000 dilution used.

Arrows indicate modified deiminated histones (Lane 2, 10-15kDa) and deiminated keratins (Lane 3, 55-65kDa).

Signals were not detected in the control membrane (B) without chemical modification of PVDF membrane.

**References :**

1. Senshu, T. et al., Detection of citrulline residues in deiminated proteins on polyvinylidene difluoride membrane. *Anal. Biochem.* 203 94-100 (1992)
2. Senshu, T. et al., Detection of deiminated proteins in rat skin: probing with a monospecific antibody after modification of citrulline residues. *J. Invest. Dermatol.* 105 163-169 (1995)
3. Ishigami, A. et al., Abnormal accumulation of citrullinated proteins catalyzed by peptidylarginine deiminase in hippocampal extracts from patients with Alzheimer's disease. *J. Neurosci. Res.* 80 120-128 (2005)

***For Research Use Only***

---

## **Chemically modification of deiminated proteins - Western blot protocol**

### ■ Reagents (necessary but not included)

- 98%  $\text{H}_2\text{SO}_4$
- 85%  $\text{H}_3\text{PO}_4$
- Acetic acid
- Diacetyl Monoxime [2,3-Butanedione monoxime] (Nacalai Tesque 10923-52, SIGMA-ALDRICH B0753)
- Antipyrine (Wako Pure Chemical 018-11112, SIGMA-ALDRICH, A5882)
- $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$  [Iron(III) Chloride Hexahydrate] (Wako Pure Chemical 091-00872, SIGMA-ALDRICH 236489)

### ■ Stock solutions

#### 1. Reagent I (100 mL) store at room temperature

Water 55 mL

Please add following  $\text{H}_2\text{SO}_4$  and  $\text{H}_3\text{PO}_4$  slowly in small quantities while stirring water in a stirrer because temperature of the solution rises.

$\text{H}_2\text{SO}_4$  20 mL

$\text{H}_3\text{PO}_4$  25 mL

Please add following  $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$  and dissolve it after temperature of the solution reached room temperature.

$\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$  41.6 mg

#### 2. Reagent II (100 mL) store at 4 °C

Diacetyl Monoxime 2 g (final conc.2%)

Antipyrine 1 g (final conc.1%)

Water 80 mL

Please stir it in a stirrer and dissolve it.

Acetic acid (17.4 M) 5.75 mL (final conc.1 M)

Adjust to 100 mL with water

### ■ Experimental protocol

## **Electrophoresis and transfer**

1. Prepare the electrophoresis samples.
2. Perform SDS-polyacrylamide gel electrophoresis (SDS-PAGE) and transfer the proteins to PVDF membrane.

**Optional:** If the expected deiminated proteins are not stably retained on the PVDF membrane, please perform the following procedure to improve protein retention. If it does not need, please proceed to the next step (Chemical modification of deiminated proteins).

1. Wash the transferred PVDF membrane in water.
2. Immerse the PVDF membrane in 0.1% ovalbumin in phosphate-buffered saline (PBS) for 15 min with agitation.
3. Wash the PVDF membrane in water.
4. Immerse the PVDF membrane in 4% paraformaldehyde in PBS for 15 min with agitation.
5. Wash the PVDF membrane several times in water.

#### **Chemical modification of deiminated proteins (citrulline residues)**

6. Prepare chemically modification solution (Regent I : Regent II = 2 : 1).
7. Immerse the PVDF membrane and incubate at 37°C more than 2 hour (or overnight) without agitation.
8. Wash the PVDF membrane several times in water.

~~✕if you check the specificity of reaction, please skip this procedure.~~

#### **Immunoreactions and detection**

1. Blocking  
Block the membrane in 5% Skim milk/TBST (10 mM Tris-HCl, pH 7.4, 0.14 M NaCl, 0.1% Tween-20) for 1 hour at room temperature with constant agitation.
2. Reaction of primary antibody  
Incubate the membrane with citrullinated protein detection antibody [rabbit anti-modified citrulline (AMC) antibody] diluted 1:500-1,000 in 10 ml 5% Skim milk/TBST for 3 hour at room temperature (or overnight at 4 °C).
3. Wash the membrane three times in TBST for 5 minutes each time.
4. Reaction of secondary antibody (you can use a commercial antibody)  
Method of using a Goat Anti-Rabbit IgG, HRP-conjugate for secondary antibody  
Incubate the membrane with secondary antibody diluted in 5% Skim milk/TBST for 1 hour at room temperature.
5. Wash the membrane three times in TBST for 5 minutes each time.
6. Detection  
Use the detection method of your choice. We recommend enhanced chemiluminescence reagent.