



Manual for Topoisomerase II Drug Screening Kit

Catalog No. 1009-1

100 Assays


Catalog No. 1009-2

250 Assays

Lot Number:

Kit Description

The TopoGEN Topoisomerase II Drug Screening Kit contains reagents necessary to identify compounds that affect enzyme activity. This kit will allow detection of two kinds of topoisomerase inhibitors: those that stimulate formation of cleavable complexes and those that antagonize topoisomerase (topo) II action on the DNA. A known topo II poison is included as a control. The DNA substrate (pRYG) included in this assay is ideal for these studies because it contains a single, high affinity topo II cleavage and recognition site. Topo II cuts pRYG primarily at a single site at the 54 bp of alternating purine/pyrimidine DNA (see Spitzner et al., 1990). Given that cleavage analyses require larger amounts of enzyme compared to catalytic assays, it is possible to conserve on enzyme by using a DNA substrate with a high affinity site. The assay system is based upon evaluating the formation of DNA cleavage products which may be nicked, open circular DNA or linear DNA. The products are then resolved by ethidium bromide gel electrophoresis as described in the kit (refer to diagram below).



Storage of Reagents

The DNAs should be stored at 4° C and the VP-16 and buffers stored at -20° C upon receipt.

Product Application and Disclaimer

This product is not licensed or approved for administration to humans or animals. It may be used with experimental animals only. The product is for in vitro research diagnostic studies only. The product is non-infectious and non-hazardous to human health. This information is based on present knowledge and does not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. TopoGEN, Inc. shall not be held liable for product failure due to mishandling and incorrect storage by end user. TopoGEN's liability is limited to credit or product replacement.

TopoGEN, Inc.

For orders: www.topogen.com

page 1

For Technical Questions: 614-451-5810

or support@topogen.com ©TopoGEN, Inc. 2007

Kit Contents

(Volumes are given for the 100 assay kit size; for 250 assay kits, multiply volumes by 2.5.)

1. pRYG DNA, supercoiled (25 ug in 100 ul) = 0.25 ug/ul
2. Marker DNA, Linear pRYG; 0.05 ug/ul in 1x gel loading buffer (25 ul, load 2 ul as marker).
- *4. 10x Topo II Incomplete Assay Buffer A contains 0.5 M Tris-HCl (pH 8), 1.50 M NaCl, 100 mM MgCl₂, 5 mM Dithiothreitol, 300 ug BSA/ml
- *5. 10x ATP Buffer B contains 20mM ATP in water.

**You must mix Buffers A and B together prior to make a 5x Complete Assay Buffer. To prepare a fresh stock of the 5x Assay buffer: Add equal volumes Buffer A and B to give the Complete 5x Assay Buffer (example, if you need 50 ul of 5x Complete Buffer for a single experiment, mix 25 ul of Buffer B with 25 ul of Buffer A). The Complete 5x Buffer MUST be made fresh for each experiment. Prepare only the amount needed fresh each day. DO NOT SAVE THE 5X COMPLETE ASSAY BUFFER, IT IS NOT STABLE.*

6. VP-16 (etoposide) control drug (Cat# 4140) in lyophilized state. Add 250 ul of methanol or DMSO to give a 10 mM Stock solution. Use at 1/10 volume or 1mM final concentration to detect cleavages.
7. 10% sodium dodecyl sulfate (SDS) (300 ul): to terminate reactions, use 1/10 volume.
8. 10x gel loading buffer (300 ul): 0.25% bromophenol blue, 50% glycerol: Use 0.1 vol.

Purified Topoisomerase II is not included in this kit but is available on line at www.topogen.com (Catalog number #2000H1

TOPOISOMERASE II DRUG SCREENING KIT PROTOCOL

Reaction volumes should be 20-30 ul final volume (limited by the volume that can be loaded into the wells of the agarose gel). The reactions should be assembled on ice in microcentrifuge tubes (water, buffer, and DNA, test compound and enzyme last). After adding enzyme, the tubes should be transferred to a heating block to initiate the reaction. Reactions should be incubated 30 min (37°C), terminated by rapid addition of 1/10 volume of 10% SDS followed by digestion with proteinase K (not supplied) prior to loading the gel. SDS is added to reactions at 37° to facilitate trapping the enzyme in a cleavage complex. If the reactions are heated, cooled or treated with high salt prior to SDS, the topo II breakage and resealing equilibrium may be altered and breaks can reseal. A sample reaction is shown below. The investigator may elect to use assay buffer if relaxation (catalytic) activity is being tested or cleavage buffer if cleavage complex formation is being tested. The latter is more difficult to detect and as noted below requires higher levels of enzyme. Note also that cleavage buffer is not optimal but will support relaxation activity.

For Reaction Mixture of 20 ul:

H ₂ O	Variable , make up to volume (20 ul in this case)
5x Complete Buffer	4 ul (prepared fresh as described above)
DNA	1 ul (200 to 250 ng is sufficient)
Test Compound	Variable. (Limited by solvent effects, see controls below) (Inhibitor control should be used at a final concentration of 0.1 mM or use 2 ul of a 1 mM stock).
Topoisomerase II	Variable. Enzyme is supplied at 2 to 20 units/ul (depends on lot#). Usually, 2 ul of stock enzyme is enough to detect cleavages; however, up to 4 ul can be tested in a 20 ul reaction but not more (due to salts in the enzyme preparation). Note that the amount of cleavage will be low, but still detectable. You should not expect to see 100% conversion of substrate to linear DNA (see comments below).

- Incubate 30 minutes at 37°C and stop by addition of 2 ul 10% SDS.

- Add proteinase K to 50 ug/ml, (incubate 37°C for 15 min.); add 0.1 vol. loading buffer. Samples may be loaded directly onto the agarose gel at this point. Alternatively, the samples can be cleaned up by extraction and then loaded: Add equal volume (20 ul) of Chloroform: isoamyl Alcohol or CIA (24:1), vortex briefly; spin in a microfuge for 5 sec. Withdraw blue colored upper aqueous phase and load onto agarose gel. Note that CIA extraction is optional but usually can improve the cosmetic quality of the agarose gel separation.

- Run 1% agarose until the dye front is at the bottom of the gel. Stain with 0.5 ug/ml ethidium bromide and photodocument results. This is a “non-ethidium bromide” gel separation which is optimal for resolving relaxed and supercoiled DNAs (see gel data, page 5). We recommend that each researcher use this method to verify enzyme activity, PRIOR to testing new compounds. Once enzyme activity is demonstrated, the researcher may wish to run gels containing 0.5 ug/ml ethidium bromide. The buffer and gel should contain ethidium bromide DURING the electrophoretic separation (followed by extensive destaining with water prior to photodocumentation). The ethidium bromide gel separation method allows one to clearly resolve cleavage products (open circular, linear DNA) from the closed DNA forms (relaxed circular and supercoiled circular DNA).



TOPOISOMERASE II DRUG SCREENING KIT PROTOCOL

Marker DNAs and controls (see Fig. 1) are extremely important. Any nuclease-free agarose of reasonable quality can be used. Agarose gels (1 %) and running buffers can be any standard nondenaturing electrophoresis buffer (example, 1 x TAE buffer, 50x stock is 242 g Tris base, 57.1 ml glacial acetic acid and 100 ml of 0.5 M EDTA) and ethidium bromide to a final concentration of 0.5 µg/ml (5 µl of a 10 mg/ml ethidium bromide (EB) stock per 100 ml of running buffer). Include EB in the gel and running buffer. Run gels at 1.5-2 V/cm (measured between electrodes) until the dye front has traveled about 80% of total gel length. The DNA can be monitored during electrophoresis with a handheld ultraviolet light. After running, gels should be destained (30 min) in water or buffer and photographed using a standard photodyne system. As noted above, we recommend EB gels to allow you to clearly resolve linear DNA (see below) from the topoisomers. If, however, you wish to determine whether your test compound is a catalytic inhibitor, you should run the gels without ethidium to resolve the topoisomers. The non-EB gels are also ideal for documenting relaxation activity of topo II.

Helpful hints and interpretation

Solvent controls are especially important to ensure that the drug solvent (DMSO or methanol for example) are not interfering with topo II activity. Also, a quick CIA extraction prior to loading the gels is a good idea since the gels will be much more cosmetic; however, the samples must first be treated with proteinase K to digest the bound protein. CIA extractions may also be required if the test compound affects the mobility of the DNA or fluorescent detection (intercalators, strong DNA binding agents, etc.). Note also, to confirm the formation of covalent complexes, you can run the following control (refer to diagram on the cover page above). Two reactions are carried out with the inhibitor; both are stopped with SDS. One is digested with proteinase K and the other is undigested. The covalent topo II/DNA complexes in the reaction that was not digested with proteinase K will be lost by CIA extraction (they will partition to the interphase) and you should lose the linears or possibly open circular DNA complexes. In contrast, the proteinase K digested sample should contain expected cleavage products (linear and nicked DNAs). Fig 1 shows a typical Gel (EB) result with the markers and controls:




Figure 1


Lane 1: Control, shows that incubating plasmid DNA with topo II gives relaxed DNA product. Note that resolution between relaxed and supercoiled DNA is not always ideal (proper conditions must be worked out empirically for each lab's gel unit system). Since a positive control of enzyme activity is critical, it may be a good idea to run both EB and non-EB gels by splitting the sample in half and loading the two different gels with equal amounts of reaction.

Lane 2: A marker of pRYG DNA. Note that a small amount of nicked OC (open circular) DNA is always present in plasmid preps.

Lane 3 shows the solvent control that demonstrates that the solvent alone is not affecting topo II activity.

Lane 4: Linear DNA marker (essential to mark position of cleavage product!).

Lane 5 shows cleavage in presence of a control poison, like m-AMSA or VP16. In the case of m-AMSA, a slight increase in nicked OC DNA is seen along with linear DNA cleavage products. Note that not all of the supercoiled DNA is converted to cleavage product. As discussed below, the cleavage event is a low yield process in most assay screens of the type.



The following controls and considerations are important (see Fig. 1)

1. Reaction containing enzyme without test drug or solvent: To detect enzymatic activity in the absence of test compound. This reaction should be carried out in assay buffer to detect relaxation.
2. Reaction containing solvent alone without enzyme: To ensure that solvent (that the test compound is dissolved in) does not affect mobility of DNA species in the assay. Also to make sure that DNA nicking independent of enzyme is not occurring. (Often DMSO is used; test a concentration of solvent at the highest concentration of compound being tested, see example in next point.)
3. Reaction containing enzyme plus solvent: To ensure that solvent (ex. DMSO or methanol) is not affecting the activity of purified enzyme. This control should be carried out with assay buffer to allow you to detect relaxation activity (run a non-EB gel). Test the concentration of solvent corresponding to the highest concentration of test compound (example if DMSO is solvent and the test drug is dissolved in 50% DMSO and a 1:100 dilution of stock drug is tested, one must show that 0.5% DMSO is not causing inhibition of enzyme activity). DMSO concentrations in this range can (depending on the DMSO source) affect topo II activity; this must be pre-determined to ensure that a range of drug concentration is used that does not cause a solvent inhibition artifact. As noted above, solvents may also cause nicking of the DNA substrate. This can create serious complications, obviously.
4. Marker DNAs: Linear DNA should be included as a reference marker. Also, we strongly recommend that supercoiled DNA be incubated in cleavage buffer without enzyme at 37° for 30 min to make sure that nicking is not occurring in the absence of topo II enzyme.
5. Control reaction with Topo II and VP-16 (etoposide) or another known topo II poison (not required with every gel, but should be performed initially to reconstruct the inhibition one should expect). As noted above, a drug control will yield linear pRYG cleavage complexes provided sufficient enzyme is used.
6. The cleavage reaction is a low yield event. That is, the amount of cleavage product will generally be only weakly detected, even with well established poisons like VP-16 (etoposide). This is because the cleavage event actively consumes the enzyme leading to stoichiometric complexes between topo II/DNA; thus, the cleavage conversion uses up the enzyme. You should expect to see some detectable level of linear DNA but it is rare to see complete conversion of supercoiled DNA into linear DNA cleavage products. It is more likely that you will see 5-25% conversion of substrate to linears. If higher cleavage conversion is desired, one can convert to a radioactive assay format using 32P-end labeled DNA and/or use high levels of topo II in order to improve sensitivity.

References

Muller et al., Biochemistry 27: 8369 (1988)
Spitzner et al., Nucleic Acids Research 18: 1-11 (1990)

To see a current list of recent papers that have used our products, go to:

http://www.topogen.com/html/new_papers.html



TopoGEN, Inc.

For orders: www.topogen.com

page 4

For Technical Questions: 614-451-5810

or support@topogen.com ©TopoGEN, Inc. 2007

Topo II Drug Screening Results

Reactions indicated above each lane were carried out with topo II (4 units, based on relaxation activity or 20 units based on decatenation assay) and 250 ng plasmid DNA. Reactions were incubated for 30 min at 37 deg C, stopped with 0.1 vol of 10% SDS, digested with proteinase K and CIA extracted as described above. A 1% agarose gel containing 0.5ug/ml EB was run until the dye front was at the bottom of the gel (2 v/cm). Formation of linear DNA is detected but the cleavages are relatively weak with VM26 under these conditions. This a typical result one sees when screenin for topo II poisons.



**Ethidium bromide
included in gel and
buffer**

Analysis of topo II relaxation activity using agarose gels without ethidium bromide (EB):

Notice that the above gel (contains EB in gel and running buffer during electrophoresis) is optimal for resolving cleavage products (linear, nicked or OC DNA); however it is not ideal for resolving form I supercoiled DNA (substrate) from relaxed DNA (product). Thus, it is not ideal to detect topo II relaxation activity. Non-ethidium bromide gels (ie, no EB in gel or buffer DURING electrophoresis) are ideal for this purpose. Gels are run and stained with EB after electrophoresis is completed. The gel to the left shows a typical result. It is easy to see that Supercoiled DNA is relaxed to form a series of topoisomers (“Relaxed” bands in the gel). Since it is essential to have good enzyme activity in order to detect cleavage products, this analysis is important to verify that topo II is active. In general, it makes sense to run both EB and non-EB gels for each experiment. Simply divide the reactions in half and run one EB and one non-EB gel.

