



Handbook for TA-cloning



pLUG[®] TA-cloning Vector Kit

pLUG[®]-Multi TA-cloning Vector Kit

Mar 2008



iNtRON
Biotechnology

pLUG[®] and pLUG[®]-Multi TA-cloning vectors were covered by Korean patent or patent pending.

Patent : 0418015, Patent Application : 0109265

pLUG[®] is a registered trademark of iNtRON Biotechnology, Inc.

© 2008 iNtRON Biotechnology, Inc. 08-0301-01



Table

Kit Contents and Storage	1
Notice to Purchaser	1
Related Products	2
Introduction	3
Overview of TA-cloning	3
Characteristics	3
How to works	4
Flow-chart of TA-cloning	4
Features of the pLUG® and pLUG®-Multi TA-cloning Vectors	5
Circle Map and Multiple Cloning Sites of the pLUG® and pLUG®-Multi TA-cloning Vectors	6
TA-cloning Protocol	8
Producing PCR products	9
Cloning into the pLUG® or pLUG®-Multi TA-cloning vector	10
Transforming the competent cells	10
Screening of a recombinant plasmid	11
Troubleshooting Guide	12
Appendix	14
Sequences of the pLUG® and pLUG®-Multi TA-cloning vectors	14
Restriction sites in the pLUG® and pLUG®-Multi TA-cloning vectors	16



Kit Contents and Storage

pLUG® and pLUG®-Multi TA-cloning vector kits

pLUG® TA-cloning Vector Kit

Cat. No.	11041
pLUG® TA-cloning vector	20 µl
5X Ligation buffer	100 µl
T4 DNA ligase	20 µl
Distilled water	1 ml

pLUG®-Multi TA-cloning Vector Kit

Cat. No.	11051
pLUG®-Multi TA-cloning vector	20 µl
5X Ligation buffer	100 µl
T4 DNA ligase	20 µl
Distilled water	1 ml

Shipping conditions


pLUG® and pLUG®-Multi TA-cloning vector kits are shipped in a refrigerated condition.

Storage recommendations

It is recommended the products to be stored at -20°C. Always avoid multiple freeze-thaw cycles or exposure to frequent temperature changes. These fluctuations can greatly alter product stability.

Notice to Purchaser

These products are for laboratory research use only.



Related Products

Additional materials listed are available from iNtRON Biotechnology. More information is available on the Internet at www.intronbio.com.

PCR-related products			
	Product	Size	Cat. No.
Dried PreMix	<i>Maxime</i> [™] PCR PreMix (<i>i-Taq</i> [™])	96 tubes	25025
		480 tubes	25026
	<i>Maxime</i> [™] PCR PreMix (<i>i-StarTaq</i> [™])	96 tubes	25165
		480 tubes	25167
DNA polymerase	<i>i-Taq</i> [™] DNA polymerase	250U	25021
		500U	25022
	<i>i-StarTaq</i> [™] DNA polymerase	250U	25161
		500U	25162
	<i>i-Max</i> [™] II DNA polymerase	250U	25261
	<i>i-Max</i> [™] DNA polymerase	250U	25041
<i>i-StarMax</i> [™] DNA polymerase	250U	25171	
Master mix solution	2X PCR Master mix Solution (<i>i-Taq</i> [™])	1 ml	25027
		5 ml	25028
	2X PCR Master mix Solution (<i>i-StarTaq</i> [™])	1 ml	25166
	2X PCR Master mix Solution (<i>i-Max</i> [™] II)	1 ml	25266
dNTPs	dNTP set	1 ml/each dNTP	32111
	dNTP mixture	2 ml	32292

Transformation-related products			
	Product	Size	Cat. No.
Competent cells	DH5 α (10^7)	1 ml	15045
	DH5 α (10^9)	1 ml	15046
	JM109 (10^7)	1 ml	15047
	JM109 (10^9)	1 ml	15048
	TOP10 (10^7)	1 ml	15049
	TOP10 (10^9)	1 ml	15050
	Blue/white colony selection	α -Complementation Solution	100 rxn.



Introduction

Overview of TA-cloning

TA-cloning technology exploits the terminal transferase activity of some DNA polymerases such as *Taq* DNA polymerase and other non-proofreading DNA polymerase. These enzymes preferentially add a 3'-end A-overhang to PCR products. This allows the direct insertion of such PCR products into the prelinearized cloning vector, which has a T-overhang on each 3'-end. This eliminates the need for restriction enzyme digestion of the vector or insert, primers with built-in restriction sites, or specially designed adapters, resulting in a much more efficient and robust cloning procedure. This technique is especially useful when compatible restriction sites are not available for the subcloning of DNA fragments.

Characteristics

pLUG® and pLUG®-Multi TA-cloning vectors were designed for the convenient and direct PCR cloning and serving a credible blue/white colony selection. The credible blue/white colony selection could be realized by iNtRON's proprietary DCS™ technology. DCS™ technology gives high percentage of true white colonies with an anticipated recombinant plasmid. DCS™ technology excludes the possibility of selection of colonies having a parental plasmid used for the preparation of TA-cloning vector or a re-circulated plasmid after loss of T-overhang of TA-cloning vector.

Note: False white colony without cloned DNA may be occurred by the undesirable mutation in the lacZ alpha gene or its regulatory elements during the plasmid propagation and the insertion of non-specific PCR products or primer-dimers into cloning vector, which may shift reading-frame.

pLUG® and pLUG®-Multi TA-cloning vectors contain several engineered restriction-enzyme recognition sites around the TA-cloning site allowing easy restriction analysis of recombinant plasmids or re-cloning to another vector. Especially, the restriction-enzyme recognition sites of pLUG®-Multi TA-cloning vector around the TA-cloning site are located in mirror-repeat pattern, which is useful in re-cloning to another vector.

pLUG® and pLUG®-Multi TA-cloning vector kits offer:

- High cloning efficiency
- High percentage of true white colony
- Credible blue / white colony selection



Introduction

How to works

The PCR product and TA-cloning vector with a T-overhang on each 3'-end can be ligated under the action of T4 DNA ligase. The mechanism of T4 DNA ligase is to form covalent phosphodiester bonds between 3'-hydroxyl ends of A-overhang in PCR product with 5'-phosphate ends of T-overhang in TA-cloning vector. ATP is required for the ligase reaction.




Figure 1. Schematic diagram for TA-cloning

Flow-chart of TA-cloning




Figure 2. Overview of TA-cloning using the pLUG® or pLUG®-Multi TA-cloning vector.



Features of the pLUG® and pLUG®-Multi TA-cloning Vectors

M13 forward/ reverse priming sites

The M13 forward/ reverse priming sites in the pLUG® and pLUG®-Multi TA-cloning vectors allow convenient sequencing.

LacZ alpha sequence

The fragment of lacZ alpha sequence in the pLUG® and pLUG®-Multi TA-cloning vectors is able to complement beta-galactosidase activity. The lacZ alpha sequence reduces the time to screen for positive clones.

Multiple cloning region

The multiple cloning region is located around the TA-cloning site in the pLUG® and pLUG®-Multi TA-cloning vectors. The restriction-enzyme recognition sites of pLUG®-Multi TA-cloning vector around the TA-cloning site are located in mirror-repeat pattern.

Origin of replication


The pLUG® TA-cloning vector has *f1* and *ColE1* origins of replication. The pLUG®-Multi TA-cloning vector has *ColE1* origin of replication which is responsible for the replication of plasmid.

Ampicillin resistance gene

The pLUG® and pLUG®-Multi TA-cloning vectors have *bla* gene encoding for beta-lactamase that confer resistance to ampicillin.

Circle Map and Multiple Cloning Sites of the pLUG® and pLUG®-Multi TA-cloning Vectors

Map and features of pLUG® TA-cloning vector



Element	Position (bp)
LacZα fragment	216 ~ 643
LacZα start codon	216
M13 reverse priming site	204 ~ 220
M13 forward priming site	439 ~ 454
Multiple cloning site	272 ~ 404
Ampicillin resistance gene	1,188 ~ 2,308
<i>f1</i> origin of replication	594 ~ 1,050
<i>ColE1</i> origin of replication	2,042 ~ 2,082

Figure 3. Map of pLUG® TA-cloning vector

Multiple cloning sequence of pLUG® TA-cloning vector

LacZα Start

CAGGAAACAGCTATGACCATGATTACGCCAAGCTCGAAATTAACCCCTCACTAAAGGGAACAAAAGC
 GTCCTTTGTCGATACTGGTACTAATGCGGTTTCGAGCTTTAATTGGGAGTGATTTCCCTTGTTTTCG

M13 reverse priming site

Sac I Not I Xba I BamHI Sma I Pst I EcoRI

TGGAGCTCCACCGGGTGCGGCCGCTCTAGAACTAGTGGATCCCCGGGCTGCAGGAATTCGATA
 ACCTCGAGGTGGCGCCACCGCCGGCGAGATCTTGATCACCTAGGGGGCCCGACGTCCTTAAGCTAT

Hind III Hind III Sal I

TCAAGCTTCAGAGCTA PCR Product AGACGTGGCAAGCTTATCGATAACCGTCGACCTTCAGGGGGGG
 AGTTCGAAAGGTCTCGA TCTGCACCGTTCGAAATAGCTATGGCAGCTGGAAGTCCCCCCC

Kpn I

CCCAGTACCAATTGCGCCCTATAGTGAGTCGTATTACAATTCACTGGCCGTCGTTTTTAC AACGTCGT
 GGGCCATGGGTTAAGCGGGATATCACTCAGCATAATGTTAAGTGACCGGCAGCAAAATGTTGCAGCA

M13 forward priming site

Circle Map and Multiple Cloning Sites of the pLUG® and pLUG®-Multi TA-cloning Vectors

Map and features of pLUG®-Multi TA-cloning vector





Figure 4. Map of pLUG®-Multi TA-cloning vector

Multiple cloning sequence of pLUG®-Multi TA-cloning vector





TA-cloning Protocol

<Rapid ligation method>

1. Thaw 5X ligation buffer, TA-cloning vector, and distilled water.
2. Prepare a ligation-reaction mixture according to the following scheme:

Component	Volume (μl)
TA-cloning vector	1
PCR product	1 ~ 4
Distilled water	Adjust to 7
5X Ligation buffer	2
T4 DNA ligase	1
Total volume	10

3. Briefly mix the ligation-reaction mixture by pipetting and incubate the tube at 20°C for 10~15 minutes.
4. Proceed transformation according to the standard protocol or store ligation-reaction mixture at -20°C until use.
5. For blue/white colony selection, the standard protocol can be used.

<Standard overnight-ligation method: Optional>

1. Thaw 10X ligation buffer, TA-cloning vector, and distilled water.
2. Prepare a ligation-reaction mixture according to the following scheme:

Component	Volume (μl)
TA-cloning vector	1
PCR product	1 ~ 4
Distilled water	Adjust to 8
10X Ligation buffer	1
T4 DNA ligase	1
Total volume	10

3. Briefly mix the ligation-reaction mixture by pipetting and incubate the tube overnight at 4°C or 16°C (e.g., in a water bath or heating block).
4. Proceed transformation according to the standard protocol or store ligation-reaction mixture at -20°C until use.
5. For blue/white colony selection, the standard protocol can be used.



TA-cloning Protocol

Producing PCR products

• PCR primer

Though a conventional PCR primer can be routinely used in TA-cloning, the use of primer with high purity is recommended. According to recent reports^{**} the different bases in the 5'-end of the forward and reverse primers could effect on the TA-cloning efficiency of PCR products. The addition of base A on the 5'-end of reverse and/or forward primers improved cloning efficiency. Moreover, when two bases (AA, AT, AC, AG) were added to the 5'-end of pair primers, TA-cloning efficiency was improved markedly. To facilitate the A-addition of PCR product, add a 10 minute extension at 72°C to the end of the PCR cycling.

* Adenosine added on the primer 5' end improved TA cloning efficiency of polymerase chain reaction products. *Analytical biochemistry* 363: 163-165, 2007

** Change of 5' terminal nucleotides of PCR primers causing variable TA-cloning efficiency. *Journal of integrative plant biology* 49: 382-385, 2007

• PCR product

The use of fresh PCR products in TA-cloning is strongly recommended. The potential presence of exonucleases may cause degradation of the nucleotide overhangs of PCR products, reducing the cloning efficiency. While it is not recommended, some PCR products can be successfully cloned after 1 week of storage at +4°C.

For efficient cloning, PCR products need to be purified to eliminate residual dNTP, enzyme, primers, and salt as well as undesired non-specific bands. Purification of PCR products prior to ligation will generally result in higher transformation efficiencies. There are several choices to clean up the PCR product. When you get a single band on a gel, the PCR product can be purified with *PCRquick-spin*™ PCR product purification kit (Cat. No.17202). If you do not obtain a single discrete band from your PCR, you may purify the PCR product with *MEGA-spin*™ agarose gel extraction kit (Cat. No.17183). During gel-purification of PCR products, avoid long exposure to UV-light. To reduce nicking of the DNA, crystal violet staining is an alternative for band visualization.

• Blunt-ended PCR products

PCR products generated using proofreading DNA polymerases such as *Pfu* and *Vent* can be used in TA-cloning procedures after the addition of a 3'-end A-overhang.

The addition of a 3'-end A-overhang can be performed as follows:

1. Add 1 unit of *Taq* DNA polymerase to the blunt-ended PCR product sample (~100 µl) and mix well. It is not necessary to change the buffer.
2. Incubate at 72°C for 10 minutes.
3. Purify the PCR product with a silica-based spin column or phenol extraction following ethanol precipitation.



TA-cloning Protocol

Cloning into the pLUG® or pLUG®-Multi TA-cloning vector

- *Optimizing insert to vector ratio*

For the ligation reaction, the optimal molar ratio of insert (i.e., PCR product) to vector has to be optimized. We recommend using a 5~10 fold molar excess of PCR product over TA-cloning vector. In some cases, a lower ratio of PCR product to vector may be sufficient for efficient ligation. The pLUG® and pLUG®-Multi TA-cloning vectors are about 3 kbp and the kit suggests adding 50 ng (1 µl) of the vector.

Table 1. Guide for the amount of PCR product to use in the ligation reaction.

PCR product size (bp)	pLUG® TA-cloning vector		pLUG®-Multi TA-cloning vector	
	5-times molar excess (ng)*	10-times molar excess (ng)*	5-times molar excess (ng)*	10-times molar excess (ng)*
500	41.9	83.8	45	90
1,000	83.8	167.7	90	180
1,500	125.8	251.5	135	270

* Calculated for 50 ng vector using the following equation:

$$\begin{aligned} & \text{ng PCR product required} \\ & = [50 \text{ ng} \times \text{PCR product size (bp)} \times \text{molar ratio}] / [\text{vector size (bp)}] \end{aligned}$$

- *Ligation-reaction*

Set up a 10 µl ligation-reaction. The 5X ligation buffer was optimized for rapid ligation which requires short ligation times at room temperature.

If you see precipitate in the ligation buffer, warm it briefly at 37°C to dissolve the salts. It is important to mix the solution completely before use to avoid localized concentrations of salts. If using non-purified PCR product, do not add more than 2 µl PCR product.

Transforming the competent cells

The pLUG® and pLUG®-Multi TA-cloning vectors are compatible with a wide range of chemically competent cells which are available from iNtRON Biotechnology. In particular, DH5α, JM109 and TOP10 competent cells are appropriate to blue/white colony selection (i.e., blue/white screening capability) and exhibit the high transformation efficiency. A 3~10 µl aliquot of ligation mixture may be used to transform competent cells. Transforming the competent cells can be performed according to a standard transformation method. Transformants should be plated onto LB plates with ampicillin, X-gal and IPTG.



TA-cloning Protocol

Background colonies may appear following transformation if the PCR template was plasmid DNA containing a resistance gene for the antibiotic used for colony selection (i.e., an ampicillin-resistance gene). In this case the PCR product should be gel-purified prior to ligation to remove template plasmid DNA.

Screening of a recombinant plasmid

Colonies with a recombinant plasmid will be white. Screening a recombinant colony with cloned DNA can be performed with white colony by colony PCR or enzymatic analysis of the plasmid isolated from its culture.

- *Colony PCR can be performed as follows:*

1. Pick a white colony with an autoclaved toothpick or pipette tip and swirl it into 20~30 μ l of dH₂O in a microcentrifuge tube.
2. Heat in a boiling water bath or hot (90~100°C) heating block for 2 minutes (can be used directly for PCR).
3. Spin for 2 minutes at top speed in microcentrifuge and transfer the supernatant into a new microcentrifuge tube.
4. Use 1~2 μ l of supernatant as a template in 20 μ l PCR reaction and perform the PCR.
5. Subject the PCR product to agarose gel electrophoresis.

- *Enzymatic analysis can be performed as follows:*

1. Pick a white colony with an autoclaved toothpick or pipette tip and swirl it into 1~3 ml LB media containing ampicillin.
2. After cultivation at 37°C, perform the isolation of plasmid DNA using our *DNA-spin*™ plasmid DNA extraction kit.
3. Treat the isolated plasmid with appropriate restriction enzyme located around the TA-cloning site.
4. Subject enzyme digest to agarose gel electrophoresis.



Troubleshooting Guide

If you do not obtain the results you expect, refer to following table. Please contact iNtRON Technical Support Department if you have questions on use of the products.

☛ **E-mail: consult@intronbio.com. Tel: +82-31-778-7830.**

Problem	Possible cause	Suggestion
No colony was obtained from transformation	Bacteria were not competent	Check the transformation efficiency of competent cells used. We recommend use highly competent cells ($\geq 1.0 \times 10^8$ cfu/ μ g circular plasmid DNA)
	Incorrect transformation procedure	Make sure that the appropriate transformation procedure was used
	Use of incorrect antibiotic on plates	Agar plates should include ampicillin; Check the ampicillin concentration added to the plates; Use fresh antibiotic plates
Low number of white colonies	A-tailing of the PCR product was not efficient	Verify that PCR amplification was performed using <i>Taq</i> DNA polymerase
		Use fresh PCR products. Efficiency may be reduced after as little as 1 day of storage because of loss of A-overhang of PCR products
		Addition of extension time at 72°C to the end of the PCR cycling will be helpful for increasing of A-tailing
	Insert to vector ratio was too low	Use of a 5~10 fold molar excess of PCR product over TA-cloning vector is recommended for ligation
	Ligation mixture used in transformation was insufficient	Increase the amount of the ligation mixture added to the transformation reaction or increase the amount of the transformants plated
	Ligation reaction time was too short	Extend the ligation reaction time. Increasing the ligation time can yield higher numbers of recombinant colonies
	Presence of inhibitor in PCR mixture	The high salt content of PCR reactions can inhibit ligation and transformation. Though your PCR product is shown as a single band, purify the PCR product prior to ligation, e.g., using the <i>PCRquick-spin</i> ™ PCR product purification kit
	Overexposure of PCR product to UV-light	The PCR product should be exposed to long-wave UV-light for as short a time as possible. Overexposure of PCR product to UV-light will lead to the formation of pyrimidine dimers that cannot be ligated efficiently



Troubleshooting Guide

Problem	Possible cause	Suggestion
White colonies did not have insert	Cells transformed with residual plasmid DNA from PCR reaction	If the PCR template was plasmid DNA containing a resistance gene for the antibiotic used for colony selection, colonies may appear following transformation. In that case, PCR product should be gel-purified prior to ligation reaction
	Non-specific PCR products or primer-dimers cloned into TA-cloning vector	Gel-purify the PCR product prior to performing the cloning reaction. Use only specified band to the ligation reaction Verify that the PCR primer design and quality
Only white colonies were obtained	No IPTG or X-gal in plates	X-gal has to be included in the plate for blue/white colony selection. If you use competent cell which genotype shows <i>lacI^s</i> , it requires IPTG to induce expression from the <i>lac</i> promoter
	Inappropriate bacterial strain was used for blue/white colony selection	Ensure that the bacterial strain used for transformation has <i>lacIqZΔM15</i> genotype
Most colonies were blue or light blue	The insert did not interrupt the reading-frame of the <i>lacZ</i> gene	If the insert is small (< 300bp) and the number of its bases including the 3'-end A-overhang is multiple of 3, the recombinant colony may be light blue colony.
	Inappropriate bacterial strain was used for blue/white colony selection	Ensure that the bacterial strain used for transformation carries a <i>lacZ</i> mutation. Phenotype of colony is always blue which is <i>lacZ⁻</i> strains
	Self-ligation of TA-cloning vector	Frequent freezing/thawing may induce loss of 3'-end T-overhang in the TA-cloning vector. Loss of the T-overhangs results in blunt-end ligation of vector itself and the colony with this plasmid is blue Nuclease may degrade the T-overhang in TA-cloning vector. Use only the provided distilled water and ligation buffer in the ligation reaction
White colonies do not grow in liquid culture	Ampicillin-sensitive satellite colonies	Pick large white colonies. Be sure that ampicillin is fresh



Appendix

Sequences of the pLUG® and pLUG®-Multi TA-cloning vectors

• pLUG® TA-cloning vector (2,982 bp)

```

10      20      30      40      50      60      70      80      90      100
GCGCCCAATA CGCAAACCGC CTCTCCCCGC GCGTTGGCCG ATTCATTAAT GCAGCTGGCA CGACAGGTTT CCCGACTGGA AAGCGGGCAG TGAGCGCAAC
110     120     130     140     150     160     170     180     190     200
GCAATTAATG TGAGTTAGCT CACTCATTAG GCACCCAGG CTTTACACTT TATGCTTCCG GCTCGTATGT TGTGTGGAAT TGTGAGCGGA TAACAATTTT
210     220     230     240     250     260     270     280     290     300
ACACAGGAAA CAGCTATGAC CATGATTACG CCAAGCTCGA AATTAACCCT CACTAAAGGG AACAAAAGCT GGAGCTCCAC CGCGGTGGCG GCCGCTCTAG
310     320     330     340     350     360     370     380     390     400
AACTAGTGGA TCCCCCGGGC TGCAGGAATT CGATATCAAG CTTCAGAGC *TA*GACGTGGC AAGCTTATCG ATACCGTCGA CCTTCAGGGG GGGCCCCGGTA
410     420     430     440     450     460     470     480     490     500
CCCAATTTCG CCTATAGTGA GTCGTATTAC AATTCAGTGG CCGTCGTTTT ACAACGTCGT GACTGGGAAA ACCCTGGCGT TACCCAACCT AATCGCCTTG
510     520     530     540     550     560     570     580     590     600
CAGCACATCC CCCTTTCGCC AGCTGGCGTA ATAGCGAAGA GGCCCCCACC GATCGCCCTT CCCAACAGTT GCGCAGCCTG AATGGCGAAT GGCAAATGTT
610     620     630     640     650     660     670     680     690     700
AAGCGTTAAT ATTTTGTAA AATTCGCGTT AAATTTTGT TAAATCAGCT CATTTTAA CCAATAGGCC GAAATCGGCA AAATCCCTTA TAAATCAAAA
710     720     730     740     750     760     770     780     790     800
GAATAGACCG AGATAGGGTT GAGTGTGTT CCAGTTTGGG ACAAGAGTCC ACTATTAAAG AACGTGGACT CCAACGTCAA AGGGCGAAAA ACCGTCTATC
810     820     830     840     850     860     870     880     890     900
AGGGCGATGG CCCACTACGT GAACCATCAC CCTAATCAAG TTTTTGGGG TCGAGTGCC GTAAAGCACT AAATCGGAAC CCTAAAGGGA GCCCCCGATT
910     920     930     940     950     960     970     980     990     1000
TAGAGCTTGA CGGGGAAAGC CGGCGAACGT GGCGAGAAAG GAAGGGAAGA AAGCGAAAGG AGCGGGCGCT AGGGCGCTGG CAAGTGTAGC GGTCACGCTG
1010    1020    1030    1040    1050    1060    1070    1080    1090    1100
CGCGTAACCA CCACACCCGC CGCGCTTAAT GCGCCGCTAC AGGGCGCGTC AGGTGGCACT TTTCCGGGAA ATGTGCGCGG AACCCCTATT TGTTTATTTT
1110    1120    1130    1140    1150    1160    1170    1180    1190    1200
TCTAAATACA TTCAAATATG TATCCGCTCA TGAGACAATA ACCCTGATAA ATGCTTCAAT AATATTGAAA AAGGAAGAGT ATGAGTATTC AACATTTCCG
1210    1220    1230    1240    1250    1260    1270    1280    1290    1300
TGTCGCCCTT ATTCCCTTTT TTGCGGCATT TTGCCTTCCT GTTTTTGCTC ACCCAGAAAC GCTGGTGAAA GTAAAAGATG CTGAAGATCA GTTGGGTGCA
1310    1320    1330    1340    1350    1360    1370    1380    1390    1400
CGAGTGGGTT ACATCGAAGT GGATCTCAAC AGCGGTAAGA TCCTTGAGAG TTTTCGCCCC GAAGAAGCTT TTCCAATGAT GAGCACTTTT AAAGTCTGTC
1410    1420    1430    1440    1450    1460    1470    1480    1490    1500
TATGTGGCGC GGTATTATCC CGTATTGACG CCGGGCAAGA GCAACTCGGT CGCCGCATAC ACTATTCTCA GAATGACTTG GTTGAGTACT CACCAGTAC
1510    1520    1530    1540    1550    1560    1570    1580    1590    1600
AGAAAAGCAT CTTACGGATG GCATGACAGT AAGAGAATTA TGCAGTGTCT CCATAACCAT GAGTGATAAC ACTGCGGCCA ACTTACTTCT GACAACGATC
1610    1620    1630    1640    1650    1660    1670    1680    1690    1700
GGAGGACCGA AGGAGCTAAC CGCTTTTTTG CACAACATGG GGGATCATGT AACTCGCCTT GATCGTTGGG AACCGGAGCT GAATGAAGCC ATACCAAACG
1710    1720    1730    1740    1750    1760    1770    1780    1790    1800
ACGAGCGTGA CACCACGATG CCTGTAGCAA TGGCAACAAC GTTGCGCAAA CTATTAACCTG GCGAACTACT TACTTAGCT TCCCGGCAAC AATTAATAGA
1810    1820    1830    1840    1850    1860    1870    1880    1890    1900
CTGGATGGAG GCGGATAAAG TTGCAAGACC ACTTCTGCGC TCGGCCCTTC CGGCTGGCTG GTTTATGCTG GATAAATCTG GAGCCGTTGA GCGTGGGTCT
1910    1920    1930    1940    1950    1960    1970    1980    1990    2000
CGCGGTATCA TTGCAGCACT GGGGCCAGAT GGTAAGCCCT CCCGTATCGT AGTTATCTAC ACGACGGGGA GTCAGGCAAC TATGGATGAA CGAAATAGAC
2010    2020    2030    2040    2050    2060    2070    2080    2090    2100
AGATCGCTGA GATAGGTGCC TCACTGATTA AGCATTGGTA ACTGTCAGAC CAAGTTTACT CATATATACT TTAGATTGAT TTAAAACCTC ATTTTTAATT
2110    2120    2130    2140    2150    2160    2170    2180    2190    2200
TAAAAGGATC TAGGTGAAGA TCCTTTTTGA TAATCTCATG ACCAAAATCC CTTAACGTGA GTTTTCGTTT CACTGAGCGT CAGACCCCGT AGAAAAGATC
2210    2220    2230    2240    2250    2260    2270    2280    2290    2300
AAAGGATCTT CTTGAGATCC TTTTTTCTG CGCGTAATCT GCTGCTTGCA AACAAAAAAA CCACCGCTAC CAGCGGTGGT TTGTTTGGCG GATCAAGAGC
2310    2320    2330    2340    2350    2360    2370    2380    2390    2400
TACCAACTCT TTTTCCGAG GTAAGTGGCT TCAGCAGAGC GCAGATACCA AATACTGTCC TTCTAGTGTA GCCGTAGTTA GGCCACCACT TCAAGAACTC
2410    2420    2430    2440    2450    2460    2470    2480    2490    2500
TGTAGCACC CTTACATACC TCGTCTGCT AATCCTGTTA CCAGTGGCTG CTGCCAGTGG CGATAAGTCG TGTCTTACC GTTTGGACTC AAGACGATAG
2510    2520    2530    2540    2550    2560    2570    2580    2590    2600
TTACCGGATA AGGCGCAGCG GTCGGGCTGA ACGGGGGT CGTGCACACA GCCCAGCTTG GAGCGAACGA CCTACACCGA ACTGAGATAC CTACAGCGTG
2610    2620    2630    2640    2650    2660    2670    2680    2690    2700
AGCTATGAGA AAGCGCCACG CTTCCCGAAG GGAGAAAGGC GGACAGGTAT CCGGTAAGCG GCAGGGTCGG AACAGGAGAG CGCACGAGGG AGCTTCCAGG
2710    2720    2730    2740    2750    2760    2770    2780    2790    2800
GGGAAACGCC TGGTATCTTT ATAGTCCTGT CGGGTTTCGC CACCTCTGAC TTGAGCGTGC ATTTTTGTGA TGCTCGTCAG GGGGGCGGAG CCTATGGAAA
2810    2820    2830    2840    2850    2860    2870    2880    2890    2900
AACGCCAGCA ACGCGCCCTT TTTACGGTTC CTGGCCTTTT GCTGGCCTTT TGCTCACATG TTCTTTCTG CGTTATCCC TGATTCTGTG GATAACCGTA
2910    2920    2930    2940    2950    2960    2970    2980
TTACCGCCTT TGAGTGAGCT GATACCGCTC GCCGCGCCG AACGACCGAG CGCAGCGAGT CAGTGAGCGA GGAAGCGGAA GA

```




Appendix

Sequences of the pLUG® and pLUG®-Multi TA-cloning vectors

• pLUG®-Multi TA-cloning vector (2,766 bp)

```

10      20      30      40      50      60      70      80      90      100
GCGCCCAATA CGCAAACCGC CTCTCCCCGC GCGTTGGCCG ATTCATTAAT GCAGCTGGCA CGACAGGTTT CCCGACTGGA AAGCGGGCAG TGAGCGCAAC
110     120     130     140     150     160     170     180     190     200
GCAATTAATG TGAGTTAGCT CACTCATTAG GCACCCAGG CTTTACACTT TATGCTTCCG GCTCGTATGT TGTGTGGAAT TGTGAGCGGA TAACAATTTT
210     220     230     240     250     260     270     280     290     300
ACACAGGAAA CAGCTATGAC CATGATTACG AATTCGAGCT CGGTACCCGG GGATCCTCTA GAGTCGACCT GCAGGCATGC AAGCTTCCAG AGC*TA*GACGT
310     320     330     340     350     360     370     380     390     400
GGCAAGCTTG CATGCCTGCA GGTGACTCTT AGAGGATCCC CGGGTACCGA GCTCGAATTC GAGCTTGGCA CTGGCCGTCG TTTTACAACG TCGTACTGG
410     420     430     440     450     460     470     480     490     500
GAAAACCCTG GCGTTACCCA ACTTAATCGC CTTGCAGCAC ATCCCCCTTT CGCCAGCTGG CGTAATAGCG AAGAGGCCCG CACCGATCGC CCTTCCCAAC
510     520     530     540     550     560     570     580     590     600
AGTTGCGCAG CCTGAATGGC GAATGGCGCC TGATGCGGTA TTTTCTCCTT ACGCATCTGT GCGGTATTTT ACACCGCATA TGGTGCCTC TCAGTACAAT
610     620     630     640     650     660     670     680     690     700
CTGCTCTGAT GCCGCATAGT TAAGCCAGCG CCGACACCCG CCAACACCCG CTGACGCGCC CTGACGGGCT TGTCTGCTCC CGGCATCCGC TTACAGACAA
710     720     730     740     750     760     770     780     790     800
GCTGTGACCG TCTCCGGGAG CTGCATGTGT CAGAGTTTTT CACCGTCATC ACCGAAACGC GCGAGACGAA AGGGCCTCGT GATACGCCTA TTTTATAGG
810     820     830     840     850     860     870     880     890     900
TTAATGTCAT GATAATAATG GTTCTTTAGA CGTCAGGTGG CACTTTTCGG GGAAATGTGC GCGGAACCCC TATTTGTTTA TTTTCTAAA TACATTCAAA
910     920     930     940     950     960     970     980     990     1000
TATGTATCCG CTCATGAGAC AATAACCCCTG ATAAATGCTT CAATAATATT GAAAAAGGAA GAGTATGAGT ATTCAACATT TCCGTGTGCG CCTTATTCCC
1010    1020    1030    1040    1050    1060    1070    1080    1090    1100
TTTTTTGCGG CATTTTGCCT TCCTGTTTTT GCTCACCAG AAACGCTGGT GAAAGTAAA GATGCTGAAG ATCAGTTGGG TGCACGAGTG GGTACATCG
1110    1120    1130    1140    1150    1160    1170    1180    1190    1200
AACTGGATCT CAACAGCGGT AAGATCCTTG AGAGTTTTCG CCCCGAAGAA CGTTTTCCAA TGATGAGCAC TTTTAAAGTT CTGCTATGTG GCGCGGTATT
1210    1220    1230    1240    1250    1260    1270    1280    1290    1300
ATCCCGTATT GACGCCGGGC AAGAGCAACT CCGTCCCGCG ATACACTATT CTCAGAATGA CTTGGTTGAG TACTCACCAG TCACAGAAAA GCATCTTACG
1310    1320    1330    1340    1350    1360    1370    1380    1390    1400
GATGGCATGA CAGTAAGAGA ATTATGCAGT GCTGCCATAA CCATGAGTGA TAACACTGCG GCCAACTTAC TTCTGACAAC GATCGGAGGA CCGAAGGAGC
1410    1420    1430    1440    1450    1460    1470    1480    1490    1500
TAACCGCTTT TTTGCACAAC ATGGGGGATC ATGTAACTCG CTTGATCGT TGGGAACCGG AGCTGAATGA AGCCATACCA AACGACGAGC GTGACACCAC
1510    1520    1530    1540    1550    1560    1570    1580    1590    1600
GATGCCTGTA GCAATGGCAA CAACGTTGCG CAAACTATTA ACTGGCGAAC TACTTACTCT AGCTTCCCGG CAACAATTA TAGACTGGAT GGAGCGGGAT
1610    1620    1630    1640    1650    1660    1670    1680    1690    1700
AAAGTTGCAG GACCACTTCT GCGCTCGGCC CTTCCGGCTG GCTGGTTTAT TGCTGATAAA TCTGGAGCCG GTGAGCGTGG GTCTCGCGGT ATCATTGCAG
1710    1720    1730    1740    1750    1760    1770    1780    1790    1800
CACTGGGGCC AGATGGTAAG CCCTCCCGTA TCGTAGTTAT CTACACGACG GGGAGTCAGG CAACTATGGA TGAACGAAAT AGACAGATCG CTGAGATAGG
1810    1820    1830    1840    1850    1860    1870    1880    1890    1900
TGCTCACTG ATTAAGCATT GGTAAGTGTG AGACCAAGTT TACTCATATA TACTTTAGAT TGATTTAAAA CTTCATTTTT AATTTAAAAG GATCTAGGTG
1910    1920    1930    1940    1950    1960    1970    1980    1990    2000
AAGATCCTTT TTGATAATCT CATGACCAAA ATCCCTTAAC GTGAGTTTTT GTTCCACTGA GCGTCAGACC CCGTAGAAAA GATCAAAGGA TCTTCTTGAG
2010    2020    2030    2040    2050    2060    2070    2080    2090    2100
ATCCTTTTTT TCTGCGGTA ATCTGCTGCT TGCAAACAAA AAAACCACCG CTACCAGCGG TGGTTTGTGT GCCGGATCAA GAGCTACCAA CTCTTTTTCC
2110    2120    2130    2140    2150    2160    2170    2180    2190    2200
GAAGGTAACT GGCTTCAGCA GAGCGCAGAT ACCAAATACT GTCCTTCTAG TGTAGCCGTA GTTAGGCCAC CACTTCAAGA ACTCTGTAGC ACCGCCTACA
2210    2220    2230    2240    2250    2260    2270    2280    2290    2300
TACCTCGCTC TGCTAATCCT GTTACCAGTG GCTGCTGCCA GTGGCGATAA GTCGTGTCTT ACCGGGTTGG ACTCAAGACG ATAGTTACCG GATAAGGCGC
2310    2320    2330    2340    2350    2360    2370    2380    2390    2400
AGCGGTCGGG CTGAACGGGG GGTTCTGTGA CACAGCCCAG CTTGGAGCGA ACGACCTACA CCGAACTGAG ATACCTACAG CGTGAGCTAT GAGAAAGCGC
2410    2420    2430    2440    2450    2460    2470    2480    2490    2500
CACGCTTCCC GAAGGGAGAA AGGCGGACAG GTATCCGGTA AGCGGCAGGG TCGGAACAGG AGAGCGCACG AGGGAGCTTC CAGGGGAAA CGCTTGGTAT
2510    2520    2530    2540    2550    2560    2570    2580    2590    2600
CTTTATAGTC CTGTCCGGTT TCGCCACCTC TGACTTGAGC GTCGATTTTT GTGATGCTCG TCAGGGGGGC GGAGCCTATG GAAAAACGCC AGCAACGCGG
2610    2620    2630    2640    2650    2660    2670    2680    2690    2700
CCTTTTTTACG GTTCTTGCC TTTTGTGTGC CTTTTGTCTA CATGTTCTTT CTGCGTTTAT CCCCTGATTC TGTGGATAAC CGTATTACCG CCTTTGAGTG
2710    2720    2730    2740    2750    2760
AGCTGATACC GCTCGCCGCA GCCGAACGAC CGAGCGCAGC GAGTCAGTGA GCGAGGAAGC GGAAGA

```



Appendix

Restriction sites in the pLUG® and pLUG®-Multi TA-cloning vectors

• *pLUG® TA-cloning vector*

Table 2. Restriction enzymes that cut the pLUG® TA-cloning vector.

Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence	Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence
Acc113I	1	1487	agt/act	Bsa29I	1	368	at/cgat
Acc16I	2	572, 1745	tgc/gca	BsaAI	1	818	yac/gtr
Acc65I	1	397	g/gtacc	BsaHI	1	1428	gr/cgyc
AccB1I	4	130, 397, 855, 2015	g/gyrcc	BsaI	1	1901	ggtctc
AccBSI	5	189, 297, 965, 1129, 2930	gagcgg	BsaOI	6	292, 553, 1450, 1599, 2522, 2946	cgry/cg
AccI	1	377	gt/mkac	BsaWI	3	1672, 2503, 2650	w/ccggw
AclNI	1	302	a/ctagt	BscI	1	368	at/cgat
AcsI	3	326, 620, 631	r/aatty	Bse118I	2	919, 1883	r/ccggy
AcyI	1	1428	gr/cgyc	BseCI	1	368	at/cgat
AflIII	1	2856	a/crygt	Bsh1285I	6	292, 553, 1450, 1599, 2522, 2946	cgry/cg
AhdI	1	1968	gacnnn/nngtc	BshNI	4	130, 397, 855, 2015	g/gyrcc
Alw21I	4	276, 1300, 1385, 2546	gwgwcw/c	BsiEI	6	292, 553, 1450, 1599, 2522, 2946	cgry/cg
Alw44I	2	1296, 2542	g/tgcac	BsiHKAI	4	276, 1300, 1385, 2546	gwgwcw/c
AlwNI	1	2447	cagnnn/ctg	BsiI	2	1304, 2688	ctcgtg
Ama87I	1	314	c/ycgrg	BsoBI	1	314	c/ycgrg
ApaI	1	395	gggcc/c	Bsp106I	1	368	at/cgat
ApaLI	2	1296, 2542	g/tgcac	Bsp120I	1	391	g/ggccc
ApoI	3	326, 620, 631	r/aatty	Bsp143II	3	969, 977, 2616	rgcgc/y
AseI	3	46, 105, 1793	at/taat	BspCI	2	553, 1599	cgat/cg
AsnI	3	46, 105, 1793	at/taat	BspDI	1	368	at/cgat
Asp700I	1	1368	gaann/nnttc	BspHI	2	1128, 2136	t/catga
Asp718I	1	397	g/gtacc	BspLU11I	1	2856	a/catgt
AspEI	1	1968	gacnnn/nngtc	BspXI	1	368	at/cgat
AspHI	4	276, 1300, 1385, 2546	gwgwcw/c	BsrBI	5	189, 297, 965, 1129, 2930	gagcgg
AvaI	1	314	c/ycgrg	BsrDI	2	1732, 1914	gcaatg
AviII	2	572, 1745	tgc/gca	BsrFI	2	919, 1883	r/ccggy
BamHI	1	308	g/gatcc	BssAI	2	919, 1883	r/ccggy
BanI	4	130, 397, 855, 2015	g/gyrcc	BssSI	2	1304, 2688	ctcgtg
BanII	3	276, 395, 893	grgcy/c	BstD102I	5	189, 297, 965, 1129, 2930	gagcgg
BanIII	1	368	at/cgat	BstDSI	1	280	c/crygg
BbiII	1	1428	gr/cgyc	BstH2I	3	969, 977, 2616	rgcgc/y
Bbv12I	4	276, 1300, 1385, 2546	gwgwcw/c	BstI	1	308	g/gatcc
BcgI	2	370, 1452	cgannnnntgc				
BcoI	1	314	c/ycgrg				
BglI	2	582, 1850	gccnnn/nggc				
BpmI	2	274, 1883	ctggag				



Appendix

Restriction sites in the pLUG® and pLUG®-Multi TA-cloning vectors

• *pLUG® TA-cloning vector (Continued)*

Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence
BstMCI	6	292, 553, 1450, 1599, 2522, 2946	cgry/cg
BstSFI	6	320, 412, 1037, 1722, 2400, 2591	c/tryag
BstX2I	7	308, 1321, 1338, 2106, 2118, 2204, 2215	r/gatcy
BstXI	1	284	ccannnnn/ntgg
BstYI	7	308, 1321, 1338, 2106, 2118, 2204, 2215	r/gatcy
BstZI	1	289	c/ggccg
Bsu15I	1	368	at/cgat
CciNI	1	289	gc/ggccgc
Cfr10I	2	919, 1883	r/ccggy
Cfr42I	1	283	ccgc/gg
Cfr9I	1	314	c/ccggg
CfrI	4	35, 289, 438, 1575	y/ggccr
ClaI	1	368	at/cgat
DraI	3	1390, 2082, 2101	ttt/aaa
DraII	1	391	rg/gnccy
DraIII	1	818	cacnnn/gtg
DrdI	2	773, 2754	gacnnnn/nngtc
DsaI	1	280	c/crygg
EaeI	4	35, 289, 438, 1575	y/ggccr
EagI	1	289	c/ggccg
Eam1104I	2	541, 1179	ctcttc
Eam1105I	1	1968	gacnnn/nngtc
EarI	2	541, 1179	ctcttc
Ecl136II	1	274	gag/ctc
EclHKI	1	1968	gacnnn/nngtc
EclXI	1	289	c/ggccg
Eco24I	3	276, 395, 893	grgcy/c
Eco255I	1	1487	agt/act
Eco31I	1	1901	ggtctc
Eco32I	1	334	gat/atc
Eco52I	1	289	c/ggccg
Eco57I	3	387, 1286, 2334	ctgaag
Eco64I	4	130, 397, 855, 2015	g/gyrcc
Eco88I	1	314	c/ycgrg

Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence
EcoICRI	1	274	gag/ctc
EcoO109I	1	391	rg/gnccy
EcoRI	1	326	g/aattc
EcoRV	1	334	gat/atc
FriOI	3	276, 395, 893	grgcy/c
FspI	2	572, 1745	tgc/gca
GsuI	2	274, 1883	ctggag
HaeII	3	969, 977, 2616	rgcgc/y
Hin1I	1	1428	gr/cgyc
HincII	1	378	gty/rac
HindII	1	378	gty/rac
HindIII	2	338, 361	a/agctt
Hsp92I	1	1428	gr/cgyc
KpnI	1	401	ggtac/c
Ksp632I	2	541, 1179	ctcttc
KspI	1	283	ccgc/gg
MfiI	7	308, 1321, 1338, 2106, 2118, 2204, 2215	r/gatcy
MroNI	1	919	g/ccggc
MslI	4	282, 1197, 1556, 1715	caynn/nrrtg
Msp17I	1	1428	gr/cgyc
MspA1I	6	54, 282, 522, 1332, 2273, 2518	cmg/ckg
NaeI	1	921	gcc/ggc
NgoAIV	1	919	g/ccggc
NgoMI	1	919	g/ccggc
NotI	1	289	gc/ggccgc
NspBII	6	54, 282, 522, 1332, 2273, 2518	cmg/ckg
NspI	1	2860	rcatg/y
Ple19I	2	553, 1599	cgat/cg
PshBI	3	46, 105, 1793	at/taat
Psp124BI	1	276	gagct/c
Psp1406I	2	1366, 1739	aa/cgtt
PspAI	1	314	c/ccggg
PspALI	1	316	ccc/ggg
PspOMI	1	391	g/ggcc
PstI	1	324	ctgca/g



Appendix

Restriction sites in the pLUG® and pLUG®-Multi TA-cloning vectors

• *pLUG® TA-cloning vector (Continued)*

Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence
PvuI	2	553, 1599	cgat/cg
PvuII	2	54, 522	cag/ctg
RcaI	2	1128, 2136	t/catga
SacI	1	276	gagct/c
SacII	1	283	ccgc/gg
Sall	1	376	g/tcgac
ScaI	1	1487	agt/act
SfcI	6	320, 412, 1037, 1722, 2400, 2591	c/tryag
Sfr303I	1	283	ccgc/gg
SmaI	1	316	ccc/ggg
SpeI	1	302	a/ctagt
SspI	2	610, 1163	aat/att
SstI	1	276	gagct/c
SstII	1	283	ccgc/gg
VneI	2	1296, 2542	g/tgcac
VspI	3	46, 105, 1793	at/taat
XbaI	1	296	t/ctaga
XhoII	7	308, 1321, 1338, 2106, 2118, 2204, 2215	r/gatcy
XmaI	1	314	c/ccggg
XmaIII	1	289	c/ggccg
XmnI	1	1368	gaann/nnttc

Table 3. Restriction enzymes that do not cut the pLUG® TA-cloning vector.

AatI, AatII, AccB7I, AccIII, AfeI, AflII, AgeI, AocI, Aor51HI, AscI, AspI, AtsI, AvrII, Ball, BbeI, BbrPI, BbsI, BbuI, Bbv16II, BclI, BfrI, BglII, BlnI, BlpI, BpiI, Bpu1102I, Bpu14I, BpuAI, BsaBI, BsaMI, Bse21I, Bse8I, BseAI, BsePI, BseRI, BsgI, Bsh1365I, BsiMI, BsiWI, BsmBI, BsmI, Bsp119I, Bsp13I, Bsp1407I, Bsp1720I, Bsp19I, Bsp68I, BspEI, BspMI, BspTI, BsrBRI, BsrGI, BssHII, BssT1I, Bst1107I, Bst98I, BstBI, BstEII, BstPI, BstSNI, Bsu36I, CellI, CpoI, Csp45I, CspI, CvnI, Eco105I, Eco130I, Eco147I, Eco47III, Eco72I, Eco81I, Eco91I, EcoNI, EcoO65I, EcoT14I, EcoT22I, EheI, ErhI, Esp1396I, Esp3I, FauNDI, FbaI, FseI, HpaI, KasI, Kpn2I, Ksp22I, LspI, MamI, MfeI, MluI, MluNI, Mph1103I, MroI, MscI, MspCI, MunI, Mva1269I, NarI, NcoI, NdeI, NheI, NruI, NsiI, NspV, PacI, PaeI, PaeR7I, Pfl23II, PflMI, PinAI, PmaCI, Pme55I, PmeI, PmlI, Ppu10I, PpuMI, PshAI, Psp5II, PspEI, PspLI, PstNHI, RsrII, SapI, SbfI, SexAI, SfiI, Sfr274I, SfuI, Sgfi, SgrAI, SmiI, SnaBI, SphI, SphI, SrfI, Sse8387I, SseBI, SspBI, StuI, StyI, SunI, SwaI, Tth111I, Van91I, Vha464I, XcmI, XhoI, Zsp2I



Appendix

Restriction sites in the pLUG® and pLUG®-Multi TA-cloning vectors

• *pLUG®-Multi TA-cloning vector*

Table 4. Restriction enzymes that cut the pLUG®-Multi TA-cloning vector.

Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence	Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence
AatII	1	833	gacgt/c	BcoI	2	246, 339	c/ycgrg
Acc113I	1	1271	agt/act	BglI	2	516, 1634	gccnnnn/nggc
Acc16I	2	506, 1529	tgc/gca	BpmI	1	1667	ctggag
Acc65I	2	242, 343	g/gtacc	BsaII	3	526, 830, 1212	gr/cgyc
AccBII	5	130 242 343 525 1799	g/gyrcc	BsaI	1	1685	ggcttc
AccBSI	3	189, 913, 2714	gagcgg	BsaOI	5	487, 1234, 1383, 2306, 2730	cgry/cg
AccI	2	264, 323	gt/mkac	BsaWI	3	1456, 2287, 2434	w/ccggw
AcsI	2	230, 355	r/aatty	Bse118I	1	1667	r/ccggy
AcyI	3	526, 830, 1212	gr/cgyc	Bsh1285I	5	487, 1234, 1383, 2306, 2730	cgry/cg
AflIII	1	2640	a/crygt	BshNI	5	130, 242, 343, 525, 1799	g/gyrcc
AhdI	1	1752	gacnnn/nngtc	BsiEI	5	487, 1234, 1383, 2306, 2730	cgry/cg
Alw21I	6	240, 353, 587, 1084, 1169, 2330	gwgcw/c	BsiHKAI	6	240, 353, 587, 1084, 1169, 2330	gwgcw/c
Alw44I	3	583 1080 2326	g/tgcac	Bsil	3	781, 1088, 2472	ctcgtg
AlwNI	1	2231	cagnnn/ctg	BsmBI	2	714, 768	cgcttc
Ama87I	2	246, 339	c/ycgrg	BsoBI	2	246, 339	c/ycgrg
ApaLI	3	583, 1080, 2326	g/tgcac	Bsp143II	2	529, 2400	rgcgc/y
ApoI	2	230, 355	r/aatty	BspCI	2	487, 1383	cgat/cg
AseI	3	46, 105, 1577	at/taat	BspHI	3	807, 912, 1920	t/catga
AsnI	3	46, 105, 1577	at/taat	BspLU11I	1	2640	a/catgt
Asp700I	1	1152	gaann/nnttc	BspMI	2	272, 323	acctgc
Asp718I	2	242, 343	g/gtacc	BsrBI	3	189, 913, 2714	gagcgg
AspEI	1	1752	gacnnn/nngtc	BsrDI	2	1516, 1698	gcaatg
AspHI	6	240, 353, 587, 1084, 1169, 2330	gwgcw/c	BsrFI	1	1667	r/ccggy
AvaI	2	246, 339	c/ycgrg	BssAI	1	1667	r/ccggy
AviII	2	506, 1529	tgc/gca	BssSI	3	781, 1088, 2472	ctcgtg
BamHI	2	251, 334	g/gatcc	BstD102I	3	189, 913, 2714	gagcgg
BanI	5	130, 242, 343, 525, 1799	g/gyrcc	BstH2I	2	529, 2400	rgcgc/y
BanII	2	240, 353	grgcy/c	BstI	2	251, 334	g/gatcc
BbeI	1	529	ggcgc/c	BstMCI	5	487, 1234, 1383, 2306, 2730	cgry/cg
BbiII	3	526, 830, 1212	gr/cgyc	BstSFI	5	269, 316, 1506, 2184, 2375	c/tryag
BbuI	2	279, 314	gcatg/c	BstX2I	8	251, 334, 1105, 1122, 1890, 1902, 1988, 1999	r/gatcy
Bbv12I	6	240, 353, 587, 1084, 1169, 2330	gwgcw/c	BstYI	8	251, 334, 1105, 1122, 1890, 1902, 1988, 1999	r/gatcy
BcgI	1	1236	cgannnnntgc	Cfr10I	1	1667	r/ccggy



Appendix

Restriction sites in the pLUG® and pLUG®-Multi TA-cloning vectors

• *pLUG®-Multi TA-cloning vector (Continued)*

Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence
Cfr9I	2	246, 339	c/ccggg
CfrI	3	35, 372, 1359	y/ggccr
DraI	3	1174, 1866, 1885	ttt/aaa
DraII	1	772	rg/gnccy
DrdI	2	669, 2538	gacnnnn/nngtc
EaeI	3	35, 372, 1359	y/ggccr
Eam1104I	2	475, 963	ctcttc
Eam1105I	1	1752	gacnnn/nngtc
EarI	2	475, 963	ctcttc
Ecl136II	2	238, 351	gag/ctc
EclHKI	1	1752	gacnnn/nngtc
Eco24I	2	240, 353	grgcy/c
Eco255I	1	1271	agt/act
Eco31I	1	1685	ggtctc
Eco57I	2	1070, 2118	ctgaag
Eco64I	5	130, 242, 343, 525, 1799	g/gyrcc
Eco88I	2	246, 339	c/ycgrg
EcoCRI	2	238, 351	gag/ctc
EcoO109I	1	772	rg/gnccy
EcoRI	2	230, 355	g/aattc
EheI	1	527	ggc/gcc
Esp3I	2	714, 768	cgtctc
FauNDI	1	578	ca/tatg
FriOI	2	240, 353	grgcy/c
FspI	2	506, 1529	tgc/gca
GsuI	1	1667	ctggag
HaeII	2	529, 2400	rgcgc/y
HinII	3	526, 830, 1212	gr/cgyc
HincII	2	265, 324	gty/rac
HindII	2	265, 324	gty/rac
HindIII	2	281, 304	a/agctt
Hsp92I	3	526, 830, 1212	gr/cgyc
KasI	1	525	g/gcgcc
KpnI	2	246, 347	ggtac/c
Ksp632I	2	475, 963	ctcttc
MfiI	8	251, 334, 1105, 1122, 1890, 1902, 1988, 1999	r/gatcy
MslI	3	981, 1340, 1499	caynn/nnrtg

Enzyme Name	# of Site	Positions of Sites (bp)	Recognition Sequence
Msp17I	3	526, 830, 1212	gr/cgyc
MspAII	6	54, 456, 650, 1116, 2057, 2302	cmg/ckg
NarI	1	526	gg/cgcc
NdeI	1	578	ca/tatg
NspBII	6	54, 456, 650, 1116, 2057, 2302	cmg/ckg
NspI	4	279, 314, 727, 2644	rcatg/y
PaeI	2	279, 314	gcatg/c
Ple19I	2	487, 1383	cgat/cg
PshBI	3	46, 105, 1577	at/taat
Psp124BI	2	240, 353	gagct/c
Psp1406I	2	1150, 1523	aa/cgtt
PspAI	2	246, 339	c/ccggg
PspALI	2	248, 341	ccc/ggg
PstI	2	273, 320	ctgca/g
PvuI	2	487, 1383	cgat/cg
PvuII	2	54, 456	cag/ctg
RcaI	3	807, 912, 1920	t/catga
SacI	2	240, 353	gagct/c
Sall	2	263, 322	g/tcgac
SbfI	2	273, 320	cctgca/gg
ScaI	1	1271	agt/act
SfcI	5	269, 316, 1506, 2184, 2375	c/tryag
SmaI	2	248, 341	ccc/ggg
SphI	2	279, 314	gcatg/c
Sse8387I	2	273, 320	cctgca/gg
SspI	1	947	aat/att
SstI	2	240, 353	gagct/c
VneI	3	583, 1080, 2326	g/tgcac
VspI	3	46, 105, 1577	at/taat
XbaI	2	257, 328	t/ctaga
XhoII	8	251, 334, 1105, 1122, 1890, 1902, 1988, 1999	r/gatcy
XmaI	2	246, 339	c/ccggg
XmnI	1	1152	gaann/nnttc



Appendix

Restriction sites in the pLUG® and pLUG®-Multi TA-cloning vectors


- *pLUG®-Multi TA-cloning vector (Continued)*

Table 5. Restriction enzymes that do not cut the pLUG®-Multi TA-cloning vector.

AatI, AccB7I, AccIII, AclNI, AfeI, AflIII, AgeI, AocI, Aor51HI, ApaI, AscI, AspI, AtsI, AvrII, Ball, BanIII, BbrPI, BbsI, Bbv16II, BclI, BfrI, BglII, BlnI, BlpI, BpiI, Bpu1102I, Bpu14I, BpuAI, Bsa29I, BsaAI, BsaBI, BsaMI, BscI, Bse21I, Bse8I, BseAI, BseCI, BsePI, BseRI, BsgI, Bsh1365I, BsiMI, BsiWI, BsmI, Bsp106I, Bsp119I, Bsp120I, Bsp13I, Bsp1407I, Bsp1720I, Bsp19I, Bsp68I, BspDI, BspEI, BspTI, BspXI, BsrBRI, BsrGI, BssHII, BssT1I, Bst1107I, Bst98I, BstBI, BstDSI, BstEII, BstPI, BstSNI, BstXI, BstZI, Bsu15I, Bsu36I, CciNI, CellI, Cfr42I, ClaI, CpoI, Csp45I, CspI, CvnI, DraIII, DsaI, EagI, EclXI, Eco105I, Eco130I, Eco147I, Eco32I, Eco47III, Eco52I, Eco72I, Eco81I, Eco91I, EcoNI, EcoO65I, EcoRV, EcoT14I, EcoT22I, ErhI, Esp1396I, FbaI, FseI, HpaI, Kpn2I, Ksp22I, KspI, LspI, MamI, MfeI, MluI, MluNI, Mph1103I, MroI, MroNI, MscI, MspCI, MunI, Mva1269I, NaeI, NcoI, NgoAIV, NgoMI, NheI, NotI, NruI, NsiI, NspV, PacI, PaeR7I, Pfl23II, PflMI, PinAI, PmaCI, Pme55I, PmeI, PmlI, Ppu10I, PpuMI, PshAI, Psp5II, PspEI, PspLI, PspOMI, PstNHI, RsrII, SacII, SapI, SexAI, SfiI, Sfr274I, Sfr303I, SfuI, SgfI, SgrAI, SmiI, SnaBI, SpeI, SpsII, SrfI, SseBI, SspBI, SstII, StuI, StyI, SunI, Swal, Tth111I, Van91I, Vha464I, XcmI, XhoI, XmaIII, Zsp2I



pLUG® TA-cloning Vector Kit
pLUG®-Multi TA-cloning Vector Kit



iNtRON
Biotechnology

Room703, Jungang Induspia V, Sangdaewon-Dong,
Joongwon-Gu, Seongnam, Gyeonggi-Do, 462-120 KOREA
TEL : 82-505-550-5600 | FAX : 82-505-550-5660
www.intronbio.com