

TREVIGEN® Instructions

For Research Use Only. Not For Use In Diagnostic Procedures

Superoxide Dismutase Assay Kit

Reagent kit for the analysis of Superoxide Dismutase in
cell extracts

Sufficient reagents for 100 Reactions

Catalog# 7500-100-K

Superoxide Dismutase Assay Kit

Catalog# 7500-100-K

| <u>Table of Contents</u> | <u>Page</u> |
|---|-------------|
| I. Background | 1 |
| II. Precautions and Limitations | 2 |
| III. Materials Supplied | 2 |
| IV. Reagents/Equipment Required But Not Supplied | 2 |
| V. Reagent Preparation | 2 |
| VI. Assay Protocol | 2 |
| VII. Data Interpretation and SOD Activity Determination | 4 |
| VIII. References | 4 |
| IX. Related Products Available from Trevigen | 5 |

I. Background

The production of superoxide radicals by mitochondria is a substantial contributor to, if not the primary cause of pathology associated with neurodegenerative diseases, ischemia reperfusion injury, atherosclerosis and aging.¹ Superoxide Dismutases (SOD) catalyze the dismutation of the superoxide radical ($O_2^{\bullet -}$) into hydrogen peroxide (H_2O_2) and elemental oxygen (O_2) into the intermembrane space or mitochondrial matrix (Fig. 1), and thus provides an important defense against the toxicity of superoxide radicals.²




Figure 1. Hydrogen peroxide production by SODs.

Overexpression of SOD protects murine fibrosarcoma cells from apoptosis and promotes cell differentiation.³ SOD also inhibits adriamycin-induced apoptosis in murine peritoneal macrophages.⁴ Superoxide ions, generated from the conversion of xanthine to uric acid and hydrogen peroxide by xanthine oxidase (XOD), converts NBT to NBT-diformazan. NBT-diformazan absorbs light at 550 nm. SODs reduce superoxide ion concentrations and thereby lower the rate of NBT-diformazan formation. The extent of reduction in the appearance of NBT-diformazan is a measure of SOD activity present in your experimental sample (Fig. 2).




Figure 2. XOD and SOD cooperation in the inhibition of NBT-diformazan formation.

Trevigen's SOD assay is free of interference by other catalytic activities, and is ideal for assaying SOD in mammalian cell lysates. The kit contains the proper lysis buffer and the reagents needed for 100 experimental tests, 50 positive controls, and 50 negative controls. Unlike some other assay kits for SOD, this system is not greatly disturbed by trace metals. Each assay requires only about five minutes, and after a simple calculation, the percent inhibition of the formation of NBT-diformazan by SOD is converted to the relative activity of the sample.

II. Precautions and Limitations

1. For Research Use Only. Not for use in diagnostic procedures.
2. The physical, chemical and toxicological properties of the provided products may not yet have been fully investigated, therefore, Trevigen recommends the use of gloves, lab coats, and eye protection while using these chemical reagents. Trevigen assumes no liability for damage resulting from handling or contact with these products. MSDS are available on request.

III. Materials Supplied

| Component | Quantity | Storage | Catalog # |
|-------------------------|-------------|---------|-------------|
| SOD (1unit/ μ l) | 200 μ l | 4 °C | 7500-100-01 |
| 25X SOD Reaction Buffer | 12.0 ml | 4 °C | 7500-100-02 |
| Xanthine Solution | 1.5 ml | 4 °C | 7500-100-03 |
| NBT Solution | 6.0 ml | 4 °C | 7500-100-04 |
| XOD Solution | 2.0 ml | 4 °C | 7500-100-05 |
| 20X Cell Lysis Solution | 12.0 ml | 4 °C | 7500-100-06 |

IV. Reagents/Equipment Required But Not Supplied

Equipment

1. Visible spectrophotometer to read absorbance at 550 nm
2. Cuvettes (disposable or quartz, with at least a 1.5 ml volume)
3. Pipettor
4. Pipette tips
5. Pipette aid
6. Pasteur pipette and bulb
7. Centrifuge (for cell lysis)
8. Timer

Reagents

1. High quality, double-distilled H_2O

V. Reagent Preparation

Prior to each experiment, prepare the necessary amount of 1X Cell Lysis Solution by diluting the 20X Cell Lysis Solution with dH_2O . All other reagents are ready for use. Store all components at 4 °C until needed and avoid contamination.

VI. Assay Protocol

A. Cell Lysate Preparation

1. Detach adherent cells by gentle trypsinization. Count the cells and centrifuge at 250 x g for 10 minutes at 4 °C. Wash the cells once with cold 1 x PBS.
2. Suspend the pellet with 500 μ l of Lysis Solution per $1-5 \times 10^6$ cells. Mix thoroughly by repeated pipetting. We recommend a sample volume of 400 μ l.

- Transfer the suspension to a 1.5 ml tube and centrifuge at 12-14,000 x g for 5 minutes at 4 °C. Place the supernatant into a clean 1.5 ml tube. Store on ice if you intend to assay for SOD immediately, or freeze at -80 °C for future use.

B. Tissue Lysate Preparation

- Liver and other tissues may be lysed and processed in isotonic buffer (10mM Tris-Cl (pH 7.4), 200 mM mannitol, 50 mM sucrose, 1 mM EDTA) as described (7).

C. SOD Assay Procedure

- The assay is performed at room temperature. All components, except cell lysates, should be brought to room temperature before use. The total reaction volume is 1.5 ml. The volume of the reagent components is 107.5 µl. Therefore, the volume of deionized water required is: 1500 µl – 107.5 µl – vol. sample. Briefly vortex each reagent immediately before use.
- To a disposable cuvette add the following components in order:

| | |
|----------------------|---------------|
| dH ₂ O | From Step C-1 |
| 25 x Reaction Buffer | 60 µl |
| Xanthine Solution | 7.5 µl |
- Mix thoroughly by repeated pipetting with a clean Pasteur pipette.
- Add 30 µl of NBT Solution and repeat Step C-3.
- Add your cell lysate and repeat Step C-3.
- Place the cuvette into a spectrophotometer, read absorbance at 550 nm or set the absorbance reading to zero.
- Just before use, briefly vortex the Xanthine Oxidase (XOD) Solution and add 10 µl to the cuvette. Quickly repeat Step C-3.
- Immediately place the cuvette in the spectrophotometer, start a timer or stopwatch and record the absorbance reading every 30 seconds for 5 minutes. The first time point will be at 30 seconds and the final time point will be at 5 minutes 30 seconds.

D. Controls

1. Negative Control


A negative control must be performed. It includes all components except SOD or cell lysate. In this case, the increase in absorbance due to generation of superoxide radical proceeds maximally.

2. Positive Control

The kit contains sufficient SOD for generating 10 standard curves or 50 individual positive controls. The SOD has an activity of 1 unit/µl (1 unit is the amount of SOD which inhibits the rate of increase in absorbance due to NBTdiformazan formation by 50%). A typical standard curve would include the following SOD concentrations: 0.01 unit, 0.1 unit, 1 unit, and 10 units. For the 0.01 and 0.1 unit points, dilute 1 µl of SOD to 100 µl and 10 µl, respectively,

with 1 x Reaction Buffer and use 1 µl of each. Add 1 µl and 10 µl of undiluted SOD for the 1 unit and 10 unit activity points, respectively. Briefly vortex the SOD immediately before use.

Figure 3: Plot of SOD concentration vs. % inhibition of the rate of increase of absorbance at 550nm due to the reduction of NBT to NBT-diformazan by the superoxide radical (O₂^{•-}).



VII. Data Interpretation and SOD Activity Determination

- Determine the rate of increase in absorbance units (A) per minute for the negative control and for the test sample(s).

$$\frac{A_{550} @ 5\text{min.}30 \text{ sec.} - A_{550} @ 30 \text{ sec.}}{5 \text{ min}} = \Delta A_{550n}/\text{minute}$$

- Determine the % inhibition for the test sample(s):

$$\frac{[(\Delta A_{550n}/\text{minute})_{\text{negative control}} - (\Delta A_{550n}/\text{minute})_{\text{test}}]}{(\Delta A_{550n}/\text{minute})_{\text{negative control}}} \times 100 = \% \text{ Inhibition}$$

VIII. References

- Kussmaul, L., and Hirst J. 2006. The mechanism of superoxide production by NADH:ubiquinone oxidoreductase (complex I) from bovine heart mitochondria. PNAS **103**, 7607-12.
- Szeto, H.H. 2006. Mitochondria-targeted peptide antioxidants: novel neuroprotective agents. AAPS J. **8**, E521-31.
- Zhao, Y., K.K. Kiningham, S.M. Lin, and D.K. St. Clair. 2001. Overexpression of MnSOD protects murine fibrosarcoma cells from apoptosis and promotes a differentiation program upon treatment with 5-azacytidine. Antioxid. Redox Signal. **3**:375-386.

4. Dominguez-Rodriguez, J.R., et al. 2001. *In vivo* inhibition by antioxidants of adriamycin-induced apoptosis in murine peritoneal macrophages. *Anticancer Res.* **21**:1869-1872.
5. Beauchamp, C. and I. Fridovich. 1971. Superoxide Dismutase: Improved assays and an assay applicable to acrylamide gels. *Anal. Biochem.* **44**:276-287.
6. Fridovich, I. 1989. Superoxide dismutase: An adaptation to a paramagnetic gas. *J. Biol. Chem.* **264**. 7761-7764.
7. Beyer, W.F. and I. Fridovich. 1987. Assaying for superoxide dismutase: Some consequences of minor changes in conditions. *Anal. Biochem.* **161**: 559-566.
8. Sutherland, M.W. and B.A. Learmonth. 1997. The tetrazolium dyes MTS and XTT provide new quantitative assays for superoxide and superoxide dismutase. *Free Rad. Res.* **27**. 283-289.
9. Okado-Matsumoto, A. and I. Fridovich. 2001. Subcellular distribution of superoxide dismutases (SOD) in rat liver. *J. Biol. Chem.* **276**:38388-38393.

IX. Related Products Available From Trevigen

Contact Trevigen for details of our unique product line for studying DNA damage and repair. All of Trevigen's kits include highly qualified enzymes, substrates, buffers, full instructions for use, and a synopsis specific for your kit.

PARP Assay Kits:

| Catalog # | Description | Size |
|-------------|--|------------|
| 4667-50-K | PARP Activity Assay Kit | 50 tests |
| 4677-096-K | HT Universal Colorimetric PARP Assay w/ Histone Coated Strip Wells | 96 samples |
| 4676-096-K | Universal Chemiluminescent PARP Assay w/Histone Coated Strip Wells | 96 samples |
| 4667-250-01 | Recombinant Human PARP Enzyme | 250 µl |
| 4668-100-1 | Recombinant Human PARP (High Specific Activity) | 1000 Units |

DNA Damage Antibodies:

| Catalog # | Description | Size |
|-------------|---------------|--------|
| 4410-PC-100 | Fen-1 | 100 µl |
| 4411-PC-100 | γ-H2AX | 100 µl |
| 2372-PC-050 | p53 Ack317 | 50 µl |
| 2370-PC-050 | p53 Ack379 | 50 µl |
| 2371-PC-050 | p53 Ack382 | 50 µl |
| 2381-PC-100 | p53 total | 100 µl |
| 4350-MC-100 | UVssDNA | 100 µg |
| 4431-MC-100 | XPF | 100 µg |
| 4421-MC-100 | XRCC1 | 100 µg |
| 4354-MC-050 | anti-8-oxo-dG | 50 µl |

CometAssay™ Kits:

| Catalog # | Description | Size |
|------------|-----------------------------------|------------|
| 4250-050-K | CometAssay™ Kit | 50 samples |
| 4251-050-K | CometAssay™ Silver Staining Kit | 50 samples |
| 4252-040-K | CometAssay™ Higher Throughput Kit | 40 samples |
| 4253-096-K | CometAssay™ Kit 96 Wells | 96 samples |

FLARE™ Assay Kits:

| Catalog # | Description | Damage Recognized | Size |
|-------------|----------------------|---|-------------|
| 4040-100-FK | Fpg Kit | 8-oxoguanine, DNA containing formamidopyrimidine moieties | 75 samples |
| 4040-100-FM | | | 100 samples |
| 4045-01K-FK | Endonuclease III Kit | Thymine Glycol, 5,6-dihydro-thymine, urea, 5-hydroxy-6-hydrothymine, 5,6-dihydro-uracil, alloxan, 5-hydroxy-6-hydrouracil, uracil glycol, 5-hydroxy-5-methylhydantoin, 5-hydroxycytosine, 5-hydroxy-uracil, methyl-tartronylurea, thymine ring saturated or fragmentation product | 75 samples |
| 4045-01K-FM | | | 100 samples |
| 4130-100-FK | hOGG1 Kit | 8-oxoguanine, DNA containing formamidopyrimidine moieties | 75 samples |
| 4130-100-FM | | | 100 samples |
| 4055-100-FK | T4-PDG Kit | Cis-syn isomers of cyclobutane pyrimidine dimers | 75 samples |
| 4055-100-FM | | | 100 samples |
| 4065-100-FK | cv-PDG Kit | Cis-syn and trans-syn isomers of cyclobutane pyrimidine dimers | 75 samples |
| 4065-100-FM | | | 100 samples |
| 4100-100-FK | UVDE Kit | Cyclobutane pyrimidine dimers, (6-4) photoproducts | 75 samples |
| 4100-100-FM | | | 100 samples |

Oxidative Damage Kits

| Catalog # | Description | Size |
|------------|-------------------------------------|-----------|
| 7511-100-K | HT Glutathione Assay Kit | 384 tests |
| 7512-100-K | HT Glutathione Peroxidase Assay Kit | 480 tests |
| 7513-500-K | HT Glutathione Reductase Assay Kit | 500 tests |
| 7501-500-K | HT Superoxide Dismutase Assay Kit | 500 tests |

The product accompanying this document is intended for research use only and is not intended for diagnostic purposes or for use in humans.

Trevigen, Inc.

8405 Helgerman Ct. Gaithersburg, MD 20877

Tel: 1-800-873-8443 • 301-216-2800

Fax: 301-560-4973

e-mail: info@trevigen.com

www.trevigen.com



Please
Recycle