1. Introduction

CellEase®Bacteria II is the reagent kit for extraction of genomic DNA from bacterial cells rapidly and efficiently. You can prepare PCR (Polymerase Chain Reaction) grade template DNA with simple protocol (mix with reagent and test samples and then incubate without any purification steps).

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3. Principle

The reagent A disrupts the cells and stabilizes the genomic DNA from the samples. The reagent B degrades the cell extracts quickly. You can prepare the DNA from the cells efficiently in short time. These reagents don't have any inhibitor for PCR and also inactivate the PCR inhibitors from the cell debris. Thus, you can use the DNA samples with CellEase® Bacteria II directly to PCR. These CellEase® Bacteria II reagents don't include any toxic or harmful chemical compounds.

4. Use

Various bacteria, Yeast, etc. (Pretreatment (freeze and thawing) may necessary in case of yeast cells.)

5. Preparation

Prepare the CellEase® mixture before use. Mix well the reagent A, B (1:1). The mixture is not able to store for long time.

Table 1. Preparation of CellEase® mixture

Number of samples	Reagent A (μl)	Reagent B (μl)
1	2	2
5	10	10
10	20	20
25	50	50
50	100	100

6. Standard protocol (DNA template for PCR)

- 1) If sample is bacterial colony, a little bacteria

picked up is dissolved in 100μl distilled water, add 5μl dissolved sample into micro tube. If sample is liquid culture, add 5μl dissolved sample into micro-tube. (The optimal sample amount will vary with the type of samples).

- 2) Add 4μl of CellEase® mixture into the test samples.
- 3) Incubate the samples at 72℃ for 6 minutes.
- 4) Then, continually incubate the samples at 94℃ for 3 minutes.
- 5) Take appropriate amount of samples and use as a template DNA of PCR.


Fig. 1. The Schematic diagram for the extraction of DNA

5. Caution

- 1) You can adjust the volume of CellEase® mixture as amount of sample volume.
- 2) If it isn't enough DNA extracts with standard protocol. You can use a little longer incubation time (up to 60min) on the step of 72℃.
- 3) The rate of the mixture of reagent A, B and distilled water should be keep correctly.
- 4) The DNA extracts should be added less than 20% of the total volume of PCR.
- 5) The DNA extracts can be kept at 4℃ for a few days. Alternatively it would be stored at -20℃ until use.
- 6) If necessary, remove debris in DNA extraction fraction by using centrifugation.
- 7) CellEase®Bacteria II is for research use only. Please don't use CellEase®Bacteria II to the other purpose.


Culture (*Staphylococcus. aureus*)

Temp. 30°C

Medium LB

Time 18hr

M 1 2 3 4 5 6 7 8 9 10

M Marker(500bp ladder)

1 Undiluted Sample Add 5µl of DNA extract to PCR

2 $\times 10$ dilution

3 $\times 10^2$ dilution

4 $\times 10^3$ dilution

5 $\times 10^4$ dilution

6 Undiluted Sample Add 6µl of DNA extract to PCR

7 $\times 10$ dilution

8 $\times 10^2$ dilution

9 $\times 10^3$ dilution

10 $\times 10^4$ dilution

Protocol of CellEase Bacteria II

- 1) Mix the Reagent A and B (2ul Reagent A, 2ul Reagent B)
- 2) Directly or stepwise diluted Bacterial cells (5ul) were transferred to the tube (usually use 0.2ml or 0.5ml tubes for PCR).
- 3) Add 4ul of the mixture to the samples (5µl).
- 4) Incubate at 72°C for 6 minutes. Then incubate at 94°C for 3 minutes.
- 5) Transfer 5 or 6ul of extracts to PCR reaction mixture and amplify the target DNA fragment.

PCR reagent

5 or 6ul	Test sample
5.0 ul	$\times 10$ buffer (+Mg ²⁺)
5.0 ul	dNTPs
1.0 ul	Forward Primer (10pmol/ul)
1.0 ul	Reverse Primer (10pmol/ul)
0.5 ul	Ex Taq DNA polymerase
(5 U/ul, Takara, Japan)	

Fill up to 50ul by distilled water

PCR Cycle

94°C	1min	35 Cycles
94°C	30sec	
55°C	30sec	
72°C	60sec	
72°C	4min	

※ We can provide application data for the other kind of samples. Please don't hesitate to contact us.

※ CellEase® is the registered trademark of Biocosm Inc.

For research use only. Not for clinical diagnosis.

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