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RML  
ANAEROBIC CDC QUAD #3

# 07-401 (1/2)  
TI NO. 2410

**USE:** This product is used in qualitative procedures to identify anaerobes commonly encountered in the clinical laboratory.

**HISTORY:** The "quadrant plate" system used to identify anaerobes was developed by Lombard and Dowell at the Communicable Disease Center in Atlanta, Georgia.

**PRINCIPLE:** The fermentation of mannitol, lactose, or rhamnose is detected by acid production with a change in the Brom thymol blue indicator to yellow. Gelatin hydrolysis is determined by adding mercuric chloride to precipitate the undigested gelatin in the medium.

**FORMULA:**

Lombard Dowell Agar Base		Vitamin K.....	10.0mg
Trypticase.....	5.0 g	Agar.....	20.0 g
Yeast extract.....	5.0 g	Gelatin 1.....	4.0 g
Sodium chloride.....	2.5 g	Mannitol 2.....	6.0 g
Sodium sulfite.....	0.1 g	Lactose 3.....	6.0 g
L-tryptophan.....	0.2 g	Rhamnose 4.....	6.0 g
L-cystine.....	0.4 g	Glucose.....	1.0 g
Hemin.....	10.0mg	Deionized water.....	1000ml
pH 7.5 ± @ 25C			

**PRECAUTIONS:** This product is "For In Vitro Diagnostic Use" and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers and media after their use. Directions should be read and followed carefully.

**STORAGE:** This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 2-8C until used. Do not freeze or overheat. Allow product to come to room temperature before use. Do not incubate prior to use.

**PRODUCT DETERIORATION:** This product should not be used if (a) there is evidence of dehydration, (b) the product is contaminated, (c) the color has changed, (d) the expiration date has passed, or (e) there are other signs of deterioration.

**SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION:** The laboratory should provide suitable collection containers and instructions for their proper use. Specimens should be transported to the laboratory without delay and protected from excessive heat or cold. If there is any delay in processing a clinical specimen, the specimen should be placed in a suitable transport medium such as Amies, Stuarts, Cary-Blair, a pre-reduced media or gassed-out tube, Transgrow or Jembec System. General rules applicable to all clinical specimens:

1. Quantity of specimen should be sufficient to permit thorough examination.
2. The specimen should be collected properly and should be representative of the infected area.
3. Care should be taken to prevent contamination of specimen.
4. Specimens should be taken to the laboratory promptly.
5. Specimens should be obtained before any antibiotics are administered to the patient. If therapy was initiated prior to collection, this should be noted on the forms sent with the specimen.

For further directions on the collection and transportation of specimens, consult the chapter "Collection, Handling and Processing of Specimens" in the ASM Manual, 3rd edition, 1980.

**PROCEDURE:** The usual clinical microbiological equipment is required for procedures involving this product. The media and equipment required will depend on the identification scheme employed by the microbiologist.

**INSTRUCTIONS:** Streak each quadrant of the plate with a well isolated colony of the organism you wish to identify. Incubate at 35C in anaerobic conditions for 48 hours. Flood the Gelatin quadrant with mercuric chloride.

**INTERPRETATION:** Fermentation of any of the 3 carbohydrates is detected by a change in the indicator to yellow. Gelatin hydrolysis is indicated by a clearing around the colony growth, with the background remaining opaque.

**MATERIALS REQUIRED BUT NOT SUPPLIED:** (1) Loop sterilization device, (2) Inoculating loop, swab, collection containers, (3) Incubators, anaerobic jars or candle jars, (4) Supplemental media, (5) Quality control organisms, (6) Mercuric chloride solution-Remel #21-226.

**USER QUALITY CONTROL:** Use "Performance Test" procedures to inoculate the media with the following organisms to demonstrate positive and negative reactions for each media. Incubate and treat as described in the "Instructions" section.

remel

REGIONAL MEDIA LABORATORIES INC.

12076 SANTA FE DRIVE-LENEXA, KANSAS 66215—Phone (913)888-0939 Wats (800)255-6730 Ka. Wats (800)332-4377

<u>ORGANISM</u>	<u>MEDIA</u>	<u>RESULT</u>
Clos. perfringens	Gelatin	Positive
Clos. ramosum	Gelatin	Negative
Clos. ramosum	Mannitol	Positive
Clos. perfringens	Mannitol	Negative
Clos. perfringens	Rhamnose	Positive
Clos. ramosum	Rhamnose	Negative
Clos. perfringens	Lactose	Positive
Clos. difficile	Lactose	Negative

**STERILITY TEST:** A representative sample of each lot of media should be incubated at 35C or other appropriate temperatures at which it will be used. As a general rule, for a lot of 100 or less, a 3-5% sample should be tested. For a larger lot, 10 random plates or tubes is recommended. Media containing blood should be incubated at 35C for 48 hours and at 25C for up to 7 days as a check for psychrophiles. Report any irregularities to the manufacturer. For additional information on sterility test procedures, consult "Practical Quality Control Procedures for the Clinical Microbiology Laboratory", Cumitech 3, ASM, 1976.

**LIMITATIONS:** This product is only part of the overall scheme for identification. Procedures for biochemical and serological tests for identification may be found in appropriate references.

**BIBLIOGRAPHY:**

1. Cumitech 3, Practical Quality Control Procedures for the Clinical Microbiology Laboratory, ASM, 1976.
2. Dowell, O. R., G. L. Lombard, and A. Y. Armfield, Procedures for Use of Differential Agar Media in the Identification of Anaerobic Bacteria, U. S. Dept. HEW, PHS, CDC, Aug. 1981.
3. Finegold, S. M., W. J. Martin and E. G. Scott, Bailey and Scott's Diagnostic Microbiology, 5th edition, The C. V. Mosby Co., St. Louis, 1978.
4. Lennette, E. H., A. Balows, W. J. Hausler, Jr. and J. P. Truant, Manual of Clinical Microbiology, 3rd edition, 1980.

**CATALOG NUMBER:**

02-410 15X100mm Quad Plate 10/package

April 1981