
PORPHYRIN TEST AGAR

INTENDED USE

REMEL's Porphyrin Test Agar is a solid medium recommended for use in qualitative procedures to speciate the genus *Haemophilus* based on their ability to convert delta-aminolevulinic acid (ALA) to porphyrins.

SUMMARY AND EXPLANATION

Members of the genus *Haemophilus* require hemin (X factor) and/or nicotinamide adenine dinucleotide (NAD - V factor). The need for either one or both factors provides the main means of differentiation of these organisms. Biberstein et al. studied the action of *Haemophilus* strains on ALA and found a perfect correlation between absence of hemin requirement and the ability to convert ALA to porphyrins.¹ White and Granick found that strains of the species *H. influenzae*, *H. aegyptius* and *H. haemoglobinophilus* all lack the enzyme capacities to convert ALA to protoporphyrin.² This accounts for their dependence on hemin for growth. In 1974, Kilian described two rapid tests used to detect the ability of non-hemin requiring *Haemophilus* species to synthesize heme precursors from ALA.³ Lund and Blazevic determined the porphyrin test to be more rapid and accurate than the satellite test used previously to speciate *Haemophilus* organisms.⁴

PRINCIPLE

Porphyrin Test Agar contains hemoglobin as a source of hemin (X factor). GCHI Enrichment is a defined supplement which provides NAD (V factor), vitamins, amino acids, coenzymes, dextrose, and other nutrients for *Haemophilus* species. This medium provides an accurate means of determining X factor requirement based on the observation that X-independent *Haemophilus* strains excrete porphobilinogen and porphyrins, all of which are intermediates in the hemin biosynthetic pathway when supplied with ALA. After incubation, exposure to longwave ultraviolet light will cause an orange-red fluorescence to appear from the bacterial cells that contain porphyrins.

REAGENTS (CLASSICAL FORMULAE)*

Base Medium

Casein Peptone.....	7.5 g	Corn Starch	1.0 g
Meat Peptone.....	7.5 g	Delta-aminolevulinic Acid	Q.S.
Monopotassium Phosphate.....	1.0 g	Hemoglobin Solution	333.0 ml
Dipotassium Phosphate	4.0 g	Agar	10.0 g
Sodium Chloride.....	5.0 g	Deminerlized Water.....	667.0 ml
*GCHI Enrichment.....	10.0 ml		

pH 6.9 +/- 0.2 @ 25°C

*GCHI Enrichment:

Vitamin B-12.....	0.01 g	Glucose	100.0 g
L-Glutamine.....	10.0 g	NAD.....	0.25 g
Adenine	1.0 g	Coccarboxylase.....	0.1 g
Guanine Hydrochloride	0.03 g	Ferric Nitrate.....	0.02 g
P-Aminobenzoic Acid.....	13.0 mg	Thiamine Hydrochloride	3.0 mg
L-Cystine	1.1 g	Cysteine Hydrochloride	25.9 g
Deminerlized Water.....	1000.0 ml		

*Adjusted as required to meet performance standards.

PROCEDURE

1. Inoculate the medium by making a pencil streak across the plate with the suspected *Haemophilus* isolate.
2. Incubate in 5-10% CO₂ at 35-37°C for 18-24 hours.

INTERPRETATION OF THE TEST:

Positive Test - orange red fluorescence

Negative Test - no fluorescence

QUALITY CONTROL

All lot numbers of Porphyrin Test Agar have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL ORGANISM

Haemophilus aphrophilus ATCC® 19415
Haemophilus parahaemolyticus ATCC® 10014
Haemophilus parainfluenzae ATCC® 7901
Haemophilus influenzae ATCC® 10211

INCUBATION

CO₂, 18-24 h @ 35°C
CO₂, 18-24 h @ 35°C
CO₂, 18-24 h @ 35°C
CO₂, 18-24 h @ 35°C

RESULTS

Positive
Positive
Positive
Negative

BIBLIOGRAPHY

1. Biberstein, E.L., P.D. Mini, and M.G. Gills, 1963. J. Bacteriol. 86:814-819.
2. White, D.C. and S. Granick. 1963. J. Bacteriol. 85:842-850.
3. Kilian, M. 1974. Acta. Pathol. Microbiol. Scand. Sect. B. 82:835-842.
4. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C.

Refer to the front of the manual for **General Information** regarding precautions, product storage and deterioration, specimen collection, storage and transportation, materials required, quality control, and limitations.

ATCC® is a registered trademark of American Type Culture Collection.

IFU 1688, Revised January 16, 2003

Printed in U.S.A.

remel Apogent.

12076 Santa Fe Drive, Lenexa, KS 66215, USA

General Information: (800) 255-6730 Technical Service: (800) 447-3641 Order Entry: (800) 447-3635

Local/International Phone: (913) 888-0939 International Fax: (913) 895-4128

Website: www.remel.com Email: remel@remel.com