

研究用



Bio-X Diagnostics

BIO-FLUO SVC (BIO K 196)

**IMMUNOFLUORESCENT ASSAY FOR THE DETECTION OF
SPRING VIRAEemia OF CARP (SVC) VIRUS
IN CELL CULTURES**

I - INTRODUCTION

Spring viraemia of carp is a contagious viral disease of the Cyprinidae. Other species, such as the sheatfish (*Silurus glanis*), are also sensitive to this virus (Pasco et al., 1987). The cause of the disease is a rhabdovirus named Rhabdovirus carpio (Fijan et al., 1971), which has been described in detail by de Kinkelin and Le Berre (1974). This viral infection is found in all parts of Europe where Cyprinidae are farmed, but the disease also exists in other areas. It affects fish of all ages, the most typical victims being the "one summer" fry at the time of the spring warming of the water (Baudouy et al., 1980). The disease has a high death rate. The clinical signs of contamination are petechial haemorrhages of the skin and gills, dark coloring of the tegument, exophthalmia and a distended abdomen. Loss of balance is also seen in diseased fish. The internal lesions are characterised by petechial haemorrhages of the viscera, fibrinous peritonitis, and catarrhal or necrotic enteritis. Whilst the serological traces of an SVC viral infection indicate that serology may be a valid alternative for studying the health status of a carp population (Hattenberger-Baudouy et al., 1989), laboratory diagnosis of the disease usually involves identification of the virus in cell cultures.

II - PRINCIPLE OF THE TEST

The infected specimens are ground up in a mortar with the help of sand, then put in solution in an antibiotic-supplemented culture medium. The preparation is centrifuged and a 24-well cell culture plate inoculated with a serial dilution of the supernatant. After 1 hour's incubation at optimal temperature culture medium is added to each well and the plate is incubated until a cytopathogenic effect is observed. At this point, the cell preparation is fixed, then rinsed. SVC-specific monoclonal antibody is then added and the plate returned to the incubator. After this first incubation with monoclonal antibody the plate is rinsed, then the conjugate, goat anti-mouse FITC (fluorescein-coupled mouse-immunoglobulin-specific rabbit immunoglobulin), is added to each well and the plate incubated once more. The plate is then rinsed, mounting medium added to each well, and the cell layer observed under an inverted microscope equipped for fluorescence. If Rhabdovirus carpio is present, a green color will be seen at the sites of viral replication.

III - COMPOSITION OF THE KIT

- 1 X 100-ml bottle of washing solution: 10-fold concentrate of PBS to be diluted in 900 ml of distilled water: washing solution.

- 1 X 25-ml bottle of fixative (ready-to-use acetone solution): fixation solution.
- 1 bottle containing lyophilised monoclonal antibody specific for the virus causing spring viraemia of carp: SVC-specific monoclonal antibody.
- 1 bottle of lyophilised fluorescein-coupled anti-mouse immunoglobulin conjugate: Rabbit anti Mouse Ig FITC Conjugate.
- 1 X 25-ml bottle of a buffer solution containing Evan's Blue for diluting the conjugate: dilution solution.
- 1 X 25-ml bottle of mounting medium: mounting medium

IV - PRECAUTIONS FOR USE

- This test may be used for in vitro diagnosis only.
- The reagents must be stored at between 4 and 8° C; they may be used until the shelf-life date on the package.
- Do not use reagents from other kits.
- Discard solutions contaminated by bacteria or fungi.
- Some bottles contain merthiolate or thimerosal. This product is toxic if inhaled or if it comes in contact with the skin. Take the usual precautions when handling these bottles.
- Avoid all risks of environmental contamination by inactivating all solutions likely to contain viruses that are pathogenic for fish with a 2% sodium hypochlorite solution.

V – PROCEDURE

1. *Extracting the virus*

1.1. Preparing the specimens

Take from moribund fish or fresh corpses approximately 1-gramme fragments of spleen, kidney and brain tissue. Mix these fragments with oven-sterilised sand and grind the mixture in a mortar. After complete homogenisation is achieved add 2 ml of culture medium containing 2% foetal calf serum and antibiotics (inoculation medium). For example, one may use a mixture of 200 IU of penicillin, 200 µg of streptomycin, and 200 µg of kanamycin per ml of culture medium. A ready-to-use antibiotic + antifungal mixture that has been optimised for this purpose is available from Bio-X. This mixture avoids the problems encountered when cell cultures are inoculated with heavily-contaminated specimens. For small fish, the entire corpse may be homogenised in the mortar, ideally after the intestines have been resected.

1.2. Centrifugation of specimens

The homogenised preparation is centrifuged at between 2,000 and 4,000 g for 15 minutes at 4° C. The supernatant is collected for the subsequent steps.

1.3. Dilution of specimens

1:10, 1:100, and 1:1,000 dilutions of the supernatant are made using the inoculation medium.

2. *Isolating the virus*

2.1. Cell line selection

Cell lines FHM, EPC and BF.2 are susceptible to SVC virus. These cells may be grown in Eagle's modified MEM or with better results in Glasgow's MEM supplemented with 10% foetal calf serum, 10% phosphate tryptose and a mixture of antibiotics at the standard concentration. If a CO2 incubator is not available, the medium may be buffered at pH 7.4 with 0.16 M Tris-HCl. The optimal temperature for growth is 30° C for

stages, divide the solution into portions of the appropriate volumes and freeze at -20° C. Just before starting the test, dilute the mother solution twentyfold with a 1:10 dilution of the concentrated washing solution. Add 200 µl of the dilute monoclonal antibody solution to each well and incubate the plate at room temperature for 1 hour.

3.4. Washing the plate

Proceed as described above in point 3.2.

3.5. Adding the conjugate

Just before starting the test, dilute the mother solution twentyfold with the conjugate dilution solution. Add 200 µl of the dilute conjugate solution to each well and incubate the plate at room temperature for 1 hour.

3.6. Washing the plate

Proceed as described above in point 3.2.

3.7. Adding the mounting medium

Add 500 µl of mounting medium to each well.

3.8. Reading the results

Examine the plate under an inverted microscope equipped for fluorescence using an excitation filter designed for fluorescein. If such equipment is not available the plate may be examined by means of a normal light microscope fitted for fluorescence, provided that the plate is turned upside down and low-power lenses are used. For higher magnification (40X), use lenses with long focal lengths. This will enable you to focus on the cells despite the thickness of the plastic.

3.9. Interpreting the results

If the specimen contains the virus responsible for spring viraemia of carp, fluorescent green plaques will be visible in the cell layer. The higher the initial titre of virus in the inoculated material, the greater the number of these plaques. The dimensions of these plaques will depend on the amount of time that the virus is allowed to replicate. These plaques may exist even in the absence of a cytopathogenic effect. The fluorescence must occur in the same focal plane as the cell layer. If "crystals" of fluorescence occur, the diluted solution may be passed through a 0.22 µm filter to eliminate the aggregates. At low dilutions the specimen itself may be toxic to the cell layer. In this case or in any other situation leading to degeneration of the cell layer it may become difficult to read the results. In such an event, it is preferable to repeat the entire test in order to obtain an unambiguous diagnosis.

VI – REFERENCES

Baudouy A.-M., Danton M. et Merle G. (1980) Virémie printanière de la carpe : étude expérimentale de l'infection évoluant à différentes températures. Ann. Virol. (Inst. Pasteur) 131 E, 479-488.

Fijan N., Petrinc Z., Sulimanovic D. & Zwillenberg L.O., (1971)
Isolation of the viral causative agent from the acute form of infectious dropsy of carp. Vet. Arh 41, 125-135.

Hattenberger-Baudouy A.-M., Danton M. et Merle G., (1989) Spring viraemia of carp (SVC) : persistence of

neutralising antibodies following experimental infection. Diseases of fish and shellfish, IV EAFP International Conference, Santiago de Compostela, September 24-28, p 63.

de Kinkelin P. et Le Berre M., (1974) Rhabdovirus des poissons : II.-Propriétés in vitro du virus de la virémie printanière de la carpe. Ann. Microbiol. (Inst. Pasteur) 125 A, 113-124.

de Kinkelin P., Michel Ch. et Ghittino P. (1986) Précis de pathologie des poissons. Office International des Epizooties et Institut National de la Recherche Agronomique Editeurs, Paris, 348 p.

Pasco L., Torchy C. et de Kinkelin P. (1987) Infection expérimentale de l'alevin de silure (*Silurus glanis* L.) par le virus de la virémie printanière de la carpe (VPC). Bull. Fr. Pêche Piscic., 307, 84-88.

the FHM and EPC cells and 25° C for the BF.2 cells (de Kinkelin et al., 1986).

2.2. Preparation of the cellular substrate

The cells are kept in a Roux flask at their optimal growth temperature. One to two days before use the cells are treated with trypsin to separate them from their backing, then seeded on a 24-well plate. As susceptibility to the virus depends on the cells' age, it is advisable to use them 24-48 hours after their transfer to the plate. To guarantee the quality of the diagnosis, the cell layer must be in perfect condition at the time of inoculation with the specimens.

2.3. Inoculation

The culture medium is eliminated by turning the plate upside down over a receptacle. Use a sharp movement so as to avoid adsorption of the culture medium on the outer surface of the well. In carrying out this step, hold the plate at a reasonable distance from the receptacle to avoid all risks of contamination from splashes. After emptying the plate, quickly deposit the different dilutions of specimens, for the cell layer must be kept moist at all times. The specimens must be deposited in the wells very delicately so as not to damage the cell layer. If automatic microtip pipettes or Pasteur pipettes are used, place the tip of the pipette against the side wall of the well and release the sample material slowly. 200-µl aliquots of the different dilutions are placed in each well. Incubate the plate at 15° C for 1 hour.

2.4. Adding the inoculation medium

At the end of the viral adsorption period add gently to each well 1 ml of the 2% foetal calf serum culture medium (inoculation medium).

2.5. Incubating the plate

The plate is kept in an incubator (under 5% CO₂) at the optimal temperature for viral growth (15 °C). It is inspected daily until a cytopathogenic effect is observed. This consists of the occurrence of dense, spherical cells whose destruction results in the formation of plaques.

3. Identifying the virus

3.1. Fixation

When the cytopathogenic effect is clearly visible, which happens, for example, in the case of spring viraemia of carp, in 48-72 hours, the cell layer is fixed with the acetone solution. Eliminate the culture medium by turning the plate upside down, then add, gently, 500 µl of fixative per well. Special care must be taken to avoid adding liquid to the wells violently or letting the preparation dry out during the various steps for revealing the virus, for the cell layer could separate from the well. Incubate the plate at 4° C for 20 minutes, then eliminate the fixative by turning the plate upside down and placing it, face down, on absorbent paper for a few seconds.

3.2. Washing the plate

Add delicately 1 ml of the 1 X concentrated washing solution to each well. Wait a few minutes, then remove the washing solution by turning the plate over. Repeat this operation, then place the plate face down on absorbent paper.

3.3. Adding the monoclonal antibody

Dissolve the lyophilised monoclonal antibody in 500 µl of distilled water. If the kit is to be used in several