Labiase

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品番: OZ-30EX

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Cat#OZ-30EX

(from *Streptomyces fulvissimus* TU-6)

Labiase, produced by a submerged culture of *Streptomyces fulvissimus* TU-6, is a new enzyme preparation that lyses effectively cell walls of numerous lactic acid bacterium. This preparation is used as a tool for studies of cell walls structure of lactic acid bacterium and the preparation of plasmid DNA, intracellular enzyme and protoplast from lactic acid bacterium and so on.

SPECIFICATIONS :

Appearance.....Lyophilized powder
(containing lactose)
Activity..... ≥ 10 units/g
(N-Acetyl- β -D-glucosaminidase)
Optimum pH and temperature.....See Fig. 1,2
pH and thermal stability.....See Fig. 1,2

ASSAY FOR N-ACETYL- β -D-GLUCOSAMINIDASE ACTIVITY :

Unit Definition

One unit of enzyme catalyzes the release of 1 μ mole of 4-nitrophenol from 4-nitrophenyl N-acetyl- β -D-glucosaminide per minute at 37°C, pH 4.3.

Method

Reaction mixture		
Substrate	5 mM 4-nitrophenyl N-acetyl- β -D-glucosaminide in 10-fold diluted McIlvaine buffer, pH 4.3	0.95 mL
Enzyme	suitably diluted enzyme	0.05 mL
		Total volume 1.00 mL

Procedure

After incubation for 10 minutes at 37°C, terminate by 2 mL of 1 M Na₂CO₃. Read A₄₀₅.

Calculation

$$\text{Enzyme unit} = \frac{A_{405} \times 3}{17.7 \times 10^3 \times 10} \times 10^3$$

The molar extinction coefficient of free 4-nitrophenol under these conditions is 17,700.

ASSAY FOR LYTIC ACTIVITY TOWARD BACTERIA:

Method

Reaction mixture		
Buffer	McIlvaine buffer, pH 6.0	0.1 mL
Substrate	Bacterial suspension(OD ₆₆₀ \approx 10)	0.1 mL
Enzyme	5% solution	0.1 mL
Distilled water		0.7 mL
		Total volume 1.0 mL

Procedure

(Before adding the enzyme solution, preincubation is carried out for 5 minutes at 37°C.)
After incubation for 2 hours at 37°C with gentle shaking, OD₆₆₀ of the mixture is determined.

Calculation

$$\text{Percentage decrease in OD}_{660} = \frac{d0 - dt}{d0} \times 100$$

d(0 or t) : OD₆₆₀ of reaction mixture after 0 or t hours

Result

See Table 1.

The extent of lysis of bacterial cells by **Labiase** varies with bacterial strain, growth stage, or cultural condition.

STORAGE :

Lyophilized preparation is stable for at least 1 year when stored at 4°C.


Fig. 1 Thermal-stability and activity toward *Levilactobacillus brevis* IFO3345


Fig. 2 pH-stability and activity toward *Levilactobacillus brevis* IFO3345

- Fig. 1 Stability: The enzyme solution was adjusted to pH 6.4 and stored at each temperature for 1 hour. The lytic activity was evaluated after reacting at 37°C for 1 hour.
- Activity: The lytic activity was evaluated after reacting at each temperature for 10 minutes at pH 6.4.
- Fig. 2 Stability: The enzyme solution was adjusted to each pH and stored at 25°C for 20 hours. The lytic activity was evaluated after reacting at 37°C for 30 minutes.
- Activity: The lytic activity was evaluated after reacting at 37°C for 30 minutes at each pH.

Table 1 Lysis of bacterial cells by Labiase

Strains	Lysis
<i>Leuconostoc mesenteroides</i> subsp. <i>mesenteroides</i> IFO 3832	+
<i>Leuconostoc lactis</i> IFO 12455	+
<i>Enterococcus faecalis</i> IFO3971	±
<i>Lactobacillus acidophilus</i> IFO 13951	+
<i>Lactiplantibacillus plantarum</i> subsp. <i>plantarum</i> IAM 1041	+
<i>Lactocaseibacillus paracasei</i> subsp. <i>paracasei</i> IFO 3533	+
<i>Lactocaseibacillus rhamnosus</i> IAM 1118	+
<i>Levilactobacillus brevis</i> IFO 3345	+
<i>Limosilactobacillus fermentum</i> IFO 3071	+
<i>Limosilactobacillus fermentum</i> IFO 3961	+
<i>Fructilactobacillus fructivorans</i> IFO 13954	+
<i>Lactobacillus delbrueckii</i> subsp. <i>bulgaricus</i> JCM 1002	±
<i>Lactobacillus delbrueckii</i> subsp. <i>delbrueckii</i> JCM 1012	±
<i>Lactobacillus helveticus</i> JCM 1120	+
<i>Lactobacillus helveticus</i> JCM 1005	+
<i>Lentilactobacillus kefir</i> JCM 5818	+
<i>Fructilactobacillus sanfranciscensis</i> JCM 5668	+
<i>Streptococcus thermophilus</i> IFO 13957	+
<i>Streptococcus mutans</i> JCM 5705	±
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> IFO12732	±
<i>Bacillus subtilis</i> ISW1214	+

Lysis rate: (+) ≥ 80%, 80% > (±) ≥ 20%

Please note that the above data is an example. The values vary widely depending on the condition.

To be used for research only. DO NOT use for human gene therapy or clinical diagnosis.

Distributor

2792 Loker Avenue West, Suite 101
 Carlsbad, CA 92010 USA
 URL : <https://www.cosmobiousa.com/>
 e-mail : info@cosmobiousa.com
 Phone : +1-760-431-4600
 Fax : +1-760-431-4604

Manufacture

4-9, Imazu Dezaikae-Cho, Nishinomiya-Shi,
 Hyogo 663-8227 JAPAN
 TEL : +81-798-32-2169
 FAX : +81-798-34-7475
 URL : <https://www.ozeki.co.jp>