Labiase

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Labiase

(from *Streptomyces fulvissimus* TU-6)

Labiase, produced by a submerged culture of *Streptomyces fulvissimus* TU-6, is a new enzyme preparation that lyses effectively cell walls of numerous lactic acid bacterium. This preparation is used as a tool for studies of cell walls structure of lactic acid bacterium and the preparation of plasmid DNA, intracellular enzyme and protoplast from lactic acid bacterium and so on.

SPECIFICATIONS :

Appearance.....	Lyophilized powder (containing lactose)
Activity.....	≥10 units/g (β-N-Acetyl-D-glucosaminidase)
Optimum pH and temperature.....	See Fig. 1,2
pH and thermal stability.....	See Fig. 1,2

ASSAY FOR β-N-ACETYL-D-GLUCOSAMINIDASE ACTIVITY :

Unit Definition

One unit of enzyme catalyzes the release of 1 μmole of *p*-nitrophenol from *p*-nitrophenyl N-acetyl-β-D-glucosaminide per minute at 37°C, pH 4.3.

Method

Reaction mixture

Substrate	5 mM <i>p</i> -nitrophenyl N-acetyl-β-D-glucosaminide in 10-fold diluted McIlvaine buffer, pH 4.3	0.95 mL
Enzyme	suitably diluted enzyme	0.05 mL
		Total volume 1.00 mL

Procedure

After incubation for 10 minutes at 37°C, terminate by 2 mL of 1 M Na₂CO₃. Read A₄₀₅.

Calculation

$$\text{Enzyme unit} = \frac{A_{405} \times 3}{17.7 \times 10^3 \times 10} \times 10^3$$

The molar extinction coefficient of free *p*-nitrophenol under these conditions is 17,700.

ASSAY FOR LYtic ACTIVITY TOWARD BACTERIA:

Method

Reaction mixture

Buffer	McIlvaine buffer, pH 6.0	0.1 mL
Substrate	Bacterial suspension(OD ₆₆₀ =10)	0.1 mL
Enzyme	5% solution	0.1 mL
Distilled water		0.7 mL
Total volume		1.0 mL

Procedure

(Before adding the enzyme solution, preincubation is carried out for 5 minutes at 37°C.) After incubation for 2 hours at 37°C with gentle shaking, OD₆₆₀ of the mixture is determined.

Calculation

$$\text{Percentage decrease in OD}_{660} = \frac{d_0 - d_t}{d_0} \times 100$$

d(0 or t) : OD₆₆₀ of reaction mixture after 0 or t hours

Result

See Table 1.

The extent of lysis of bacterial cells by **Labiase** varies with bacterial strain, growth stage, or cultural condition.

STORAGE :

Lyophilized preparation is stable for at least 1 year when stored at 4°C.


Fig. 1 pH-stability and activity toward
Levilactobacillus brevis IFO3345


Fig. 2 Thermal-stability and activity toward
Levilactobacillus brevis IFO3345

Fig. 1 Stability: The enzyme solution was adjusted to each pH and stored at 25°C for 20 hours. The lytic activity was evaluated after reacting at 37°C for 30 minutes.

Activity: The lytic activity was evaluated after reacting at 37°C for 30 minutes at each pH.

Fig. 2 Stability: The enzyme solution was adjusted to pH 6.4 and stored at each temperature for 1 hour. The lytic activity was evaluated after reacting at 37°C for 1 hour.

Activity: The lytic activity was evaluated after reacting at each temperature for 10 minutes at pH 6.4.


Table 1 Lysis of bacterial cells by Labiase

Strains	Lysis
<i>Leuconostoc mesenteroides</i> subsp. <i>mesenteroides</i> IFO 3832	+
<i>Leuconostoc lactis</i> IFO 12455	+
<i>Enterococcus faecalis</i> IFO3971	±
<i>Lactobacillus acidophilus</i> IFO 13951	+
<i>Lactiplantibacillus plantarum</i> subsp. <i>plantarum</i> IAM 1041	+
<i>Lacticaseibacillus paracasei</i> subsp. <i>paracasei</i> IFO 3533	+
<i>Lacticaseibacillus rhamnosus</i> IAM 1118	+
<i>Levilactobacillus brevis</i> IFO 3345	+
<i>Limosilactobacillus fermentum</i> IFO 3071	+
<i>Limosilactobacillus fermentum</i> IFO 3961	+
<i>Fructilactobacillus fructivorans</i> IFO 13954	+
<i>Lactobacillus delbrueckii</i> subsp. <i>bulgaricus</i> JCM 1002	±
<i>Lactobacillus delbrueckii</i> subsp. <i>delbrueckii</i> JCM 1012	±
<i>Lactobacillus helveticus</i> JCM 1120	+
<i>Lactobacillus helveticus</i> JCM 1005	+
<i>Lentilactobacillus kefiri</i> JCM 5818	+
<i>Fructilactobacillus sanfranciscensis</i> JCM 5668	+
<i>Streptococcus thermophilus</i> IFO 13957	+
<i>Streptococcus mutans</i> JCM 5705	±
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> IFO12732	±
<i>Bacillus subtilis</i> ISW1214	+

Lysis rate: (+) $\geq 80\%$, 80% > (±) $\geq 20\%$

Please note that the above data is an example. The values vary widely depending on the condition.

Distributor

COSMO BIO Co., LTD.
Inspiration for Life Science

2792 Loker Avenue West, Suite 101

Carlsbad, CA 92010 USA

<https://www.cosmobiousa.com/>

e-mail : info@cosmobiousa.com

Phone / Fax : +1-760-431-4600

Manufacturer

ÖZEKI OZEKI CORPORATION