BETA LYSIN DISK TI No. 21120

INTENDED USE

REMEL's Beta Lysin Disk is a reagent-impregnated disk recommended for use in qualitative procedures to detect the production of CAMP factor by group B streptococci.

SUMMARY AND EXPLANATION

Christie, Atkins, and Munch-Petersen observed a lytic phenomenon when hemolytic and nonhemolytic group B streptococci were grown in a zone of staphylococcal beta lysin activity.¹ The CAMP test, named after the originators, detects the production of a proteinlike compound which group B streptococci produce. This substance acts synergistically with the beta toxin produced by some strains of *Staphylococcus aureus* to produce potent hemolysis. A filter paper disk impregnated with beta lysin may be used as a substitute for the beta lysin producing staphylococcus strain.

PRINCIPLE

The CAMP test is based on the formation of a zone of synergistic hemolysis at the intersect point of the Beta Lysin Disk and the group B streptococci on sheep blood agar. The reaction is caused by the beta toxin reacting with sphingomyelin in the sheep membrane to produce ceramide and leads to disorganization of the lipid bilayer of the sheep erythrocyte membrane, resulting in an area of complete lysis that is shaped like an arrowhead or crescent.²

REAGENTS

Beta Lysin Disk (25 Disks/Vial) (REMEL #21-120)
Reactive Ingredient: Beta Lysin obtained from *Staphylococcus aureus*

PRECAUTIONS

This product is *For In Vitro Diagnostic Use* and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after their use. Directions should be read and followed carefully.

STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 2-8°C until used. Allow product to come to room temperature before use. Do not incubate prior to use.

PRODUCT DETERIORATION

This product should not be used if (1) the color has changed from light beige, (2) the expiration date has passed, (3) the desiccant has changed from blue to pink, or (4) there are other signs of deterioration. Protect disks from moisture by removing from the vial only those disks necessary for testing. Promptly replace the cap and return the vial to 2-8°C.

SPECIMEN COLLECTION, STORAGE, AND TRANSPORTATION
Specimens should be collected and handled following recommended guidelines.³

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swab, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Forceps.

PROCEDURE

1. Allow a Blood Agar (TSA w/ 5% Sheep Blood) plate to equilibrate to room temperature.
2. Using forceps, place a Beta Lysin Disk on the agar surface.
3. Streak the streptococcal isolate approximately 2-3mm from the edge of the disk.
4. Incubate aerobically or in 5-10% CO₂ for 18-24 hours at 35-37°C.
5. Observe for a crescent-shaped zone of hemolysis at the intersection of the disk and the organism to indicate a positive test.

INTERPRETATION OF THE TEST

Positive Test - crescent-shaped zone of hemolysis

Negative Test - no hemolytic zone or a small helmet-shaped hemolytic zone

QUALITY CONTROL

All lot numbers of the Beta Lysin Disk have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL	INCUBATION	RESULTS
<i>Streptococcus agalactiae</i> ATCC® 12386	Aerobic, 18-24h @ 35-37°C	Positive
<i>Streptococcus pyogenes</i> ATCC® 12344	Aerobic, 18-24h @ 35-37°C	Negative

LIMITATIONS

1. Isolate tested should be a streptococcus.
2. Sheep blood agar plates containing unwashed sheep erythrocytes that have not been pretested for anti-camp factor may give a false negative CAMP test.
3. Some group A streptococci may demonstrate a helmet-shaped zone of hemolysis, especially if incubated in CO₂.
4. An agar depth of 1.5mm will result in the most distinct zone of hemolysis.
5. The Beta Lysin Disk for CAMP factor is for the presumptive identification of group B streptococci. Further testing is recommended for confirmation.

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=====

PRODUCT: 21120 LOT: 204646 EXPIRATION: 2005/10/04

PRODUCT DESCRIPTION: BETA LYSIN DISK

This product has been manufactured, processed, and packaged in accordance with Quality Systems Regulation, 21 CFR Part 820. Representative samples were tested per Remel Quality Control specifications and were found to meet performance criteria for this product.

CHALLENGE ORGANISMS	SOURCE	RESULTS
Streptococcus agalactiae	ATCC 12386	POSITIVE
Streptococcus pyogenes	ATCC 19615	NEGATIVE

Microbial Load:

Product is evaluated for performance specifications only.

Appearance: LIGHT BEIGE DISKS WITH RESIDUE

pH: NA

Unit: VIAL

Volume: 25 DISK/VIAL

Storage: 2-8C

Reference IFU:

A. Waddle
Quality Control Supervisor

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Remel Inc. 913-888-0939