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Cumitech 3, Practical Quality Control Procedures for the Clinical Microbiology Laboratory, ASH, 1976.

Paller, A., Schleifler, K.H., 1981, Modified oxidase and benzidine tests for separation of Staphylococci from Micrococci. J. Clinical Micro., 13: 1031-1035.

Pinegold, S.M., V.J. Martin, and E.G. Scott, Bailey and Scott's Diagnostic Microbiology, 5th edition, The C. V. Mosby Co., St. Louis, 1978.

Kovacs, N., 1956, Identification of *Pseudomonas pyocyanea* by the oxidase reaction. Nature (London). 178:703.

Lennette, E.H., A. Balows, V.J. Hausler, Jr., and J.P. Truant, Manual of Clinical Microbiology, 3rd edition, 1980.

Tarrand, F & Groschel, D.H., 1982, Rapid, Modified Oxidase Test for Oxidase-variable Bacterial Isolates. J. Clinical Micro., Oct. pp. 72-774.

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F A X 03(798)3239

(8 6 0 6 0 1)

USE: This product is recommended for use in qualitative procedures as a method of rapid and accurate differentiation of *Staphylococcus* from *Micrococcus* by the detection of the oxidase enzyme.

HISTORY: The oxidase method was originally described by Kovacs in 1956 as a method of differentiating gram-negative bacilli.¹ In 1981, Faller and Schleifer tested a modification of this oxidase reagent to differentiate *Staphylococcus* from *Micrococcus*.² Kovacs used a 1% aqueous solution of tetramethyl-p-phenylenediamine hydrochloride with sodium ascorbate added to prevent autoxidation. Faller and Schleifer modified Kovacs formula by using tetramethyl-p-phenylenediamine (TMPD) in dimethyl sulfoxide (DMSO). Sodium ascorbate was found to have a strong reducing power which resulted in negative results. DMSO provides solubility for the reagent, stability against autoxidation and enhances the permeation of the TMPD into the bacterial cell.⁴ Stanier et al. suggested that the oxidase test is an indirect proof for the presence of a c-type cytochrome. All micrococci contain cytochrome c, whereas all staphylococci with the exception of *Staphylococcus sciuri* lack this type of cytochrome.³

PRINCIPLE: In the presence of atmospheric oxygen, the oxidase enzyme reacts with the oxidase reagent and cytochrome c to form the colored compound, indophenol.

CLASSICAL FORMULA:*

Tetramethyl-p-phenylenediamine.....Q.S. DMSO.....Q.S.

*Adjusted as required to meet performance standards.

PRECAUTIONS: This product is "For In Vitro Diagnostic Use" and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers and media after their use. Directions should be read and followed carefully.

STORAGE: Microdase Disks should be protected from light prior to use. The product should be stored in its original container at 2-8C until used. Do not freeze.

PRODUCT DETERIORATION: This product should not be used if the expiration date has passed, or if the color of the disk has changed to a dark blue color.

PROCEDURE: The usual clinical microbiological equipment is required for procedures involving this product. Other media and equipment required will depend on the identification scheme employed by the microbiologist. Because of the extreme sensitivity of this reagent, the instructions must be followed very carefully.

INSTRUCTIONS: Testing should be performed on aerobic catalase-positive, gram-positive cocci. All colonies should be taken from a Blood Agar plate and growth must be less than 24 hours old, 15-18 hours being optimal. Place the disk in an empty petri dish or on a glass slide. **DO NOT REHYDRATE THE DISK.** Using an applicator stick, inoculate the disk with several colonies (a visible inoculum), being sure not to use a mixed culture. Apply the colonies to the TOP of the disk. Examine for up to 2 minutes and interpret a blue or purple-blue color as positive and no change (white to tan) as being negative.

MATERIALS REQUIRED BUT NOT SUPPLIED: (1) Loop sterilization device, (2) Inoculating loop, swab, collection containers, (3) Incubators, anaerobic jars or candle jars, (4) Supplemental media, (5) Quality control organisms, (6) Applicator sticks, slides or petri dishes.

QUALITY CONTROL: All lot numbers of Microdase Disks have been tested using the following quality control organisms and have been found to be acceptable.

21-132 (7/2)

CONTROL ORGANISM	INCUBATION	RESULTS
<i>Micrococcus luteus</i> ATCC 4698	Aerobic, 2 minutes @ 25-30C	(+) Blue color within 30 seconds
<i>Micrococcus kristinae</i> ATCC 27570	Aerobic, 2 minutes @ 25-30C	(+) Blue color within 2 minutes
<i>Staphylococcus aureus</i> ATCC 25923	Aerobic, 2 minutes @ 25-30C	(-) No color change at 2 minutes

LIMITATIONS:

1. Test results on organisms older than 24 hours or those taken from media other than Blood Agar are not always accurate. It is therefore recommended that the colonies tested be from a fresh culture from a Blood Agar plate.
2. Colonial morphology is often misleading, it is recommended that a gram stain and a catalase test be performed prior to the performance of a Microdase test, and that the Microdase test be performed prior to or in conjunction with coagulase, mannitol and other biochemical testing for *Staphylococcus* and *Micrococcus* speciation.
3. *Staphylococcus sciuri* is the only *Staphylococcus* species now recognized to give a positive Microdase reaction, it is recommended that if *Staphylococcus sciuri* is suspected, a Lysostaphin test (REMEL #21-130) be performed as a differentiating test. *Staphylococcus sciuri* is generally considered an animal pathogen so when working with human infections it is not a routine isolate and in most cases is not considered a potential pathogen.
4. Microdase is not designed for routine testing for oxidase activity in organisms other than *Staphylococcus* and *Micrococcus*.

BIBLIOGRAPHY:

1. Kovacs, N., *Nature* (London), 178:703, 1956.
2. Faller, A., and K.H. Schleifer, *J. Clin. Microbiol.*, 13:1031-1035, 1981.
3. Stanier, R.Y., M. Doudoroff and E.A. Adelberg, *General Microbiology*, 3rd edition, The Macmillan Press, Ltd., London, 1972.
4. Baker, J.S., *J. Clin. Microbiol.*, 19:875-879, 1984.
5. Cumitech 3, *Practical Quality Control Procedures for the Clinical Microbiology Laboratory*, ASM, 1976.
6. Finegold, S.M., W.J. Martin and E.G. Scott, *Bailey and Scott's Diagnostic Microbiology*, 5th edition, The C.V. Mosby Co., St. Louis, 1978.
7. Lennette, E.H., A. Balows, W.J. Hausler, Jr. and H.J. Shadomy, *Manual of Clinical Microbiology*, 4th edition, 1985.

ATCC is a registered trademark of American Type Culture Collection

Revised 12/22/89