



For research use only. Not for clinical diagnosis.

Catalog No. CK027

dbcAMP concentrated solution-100X
(dbcAMP-100X)

BACKGROUND

The addition of this component improves the normal oocyte maturation rate in basic medium for porcine oocyte maturation (POM) and/or high performance basic medium for porcine oocyte maturation (HP-POM). This is sterilized ready to use.

Storage	Stored at 4°C in a dark room. Do not freeze.
Package	0.5 ml
Expiration date	The expiration date is shown on the label of the medium bottle (4 months after the medium preparation)
Note	dbcAMP-100X is for research use only. Do not use in animals, humans, or diagnostic procedures.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619





For research use only. Not for clinical diagnosis.

Catalog No. CK027

dbcAMP 濃縮液 (100 倍液)
(dbcAMP-100X)

卵子成熟培養に用いるブタ卵子成熟用基本培地(POM)、又は高性能ブタ卵子成熟用基本培地(HP-POM)に添加することにより、ブタ卵子の正常成熟に効果を示し、結果として媒精における多精子受精率を低下させます。この濃縮液は滅菌されておりますので、成熟用基本培地に直接添加してご利用になれます。

保存法

本濃縮液は、冷暗所 (2~5℃) に保存して下さい。凍結厳禁。

包装

0.5 ml

有効期間

製造後 4 ヶ月

使用上の注意

- 本品は、ブタ卵子の正常成熟に利用できます。他動物等の卵子の培養については試験を行っておりません。
- ヒト卵子や胚には絶対に使用しないで下さい。
- 生体への注入は行わないで下さい。



人と科学のステキな未来へ

コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
URL : <http://www.cosmobio.co.jp/>

● 営業部 (お問い合わせ)

TEL : (03) 5632-9610 FAX : (03) 5632-9619

TEL : (03) 5632-9620



コスモ・バイオ株式会社