

Medaka vitellogenin standard for ELISA assay

Background:	Vitellogenin, the egg yolk precursor protein, is induced by estrogens. On the other hand, environmental estrogens are known to disrupt endocrine system in animals. Recently, vitellogenin has proven to be an ideal marker for environmental estrogens.
Code:	MV-STD
Lot:	
lot size:	0.5 mL
Concentration:	100ng/ml
Presentation:	Supplied as a liquid in 0.01M Tris , 0.15M NaCl, 1% BSA, 0.1% Sodium azide, pH 7.5
Source:	Ascites of 17 β -estradiol treated female medaka
Appearance:	Clear liquid
Application:	ELISA
Storage:	below -20° C for at least 6 months (recommendation: at -80° C for long term storage) Avoid freeze thaw cycles
Note:	After preparation of working standards, these solutions should be stored at -20° C for up to two weeks.
Warning:	For research use only Contains sodium azide as preservative

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp