

Medaka vitellogenin standard for ELISA assay

Code : EBT-MV-STD-EX

Lot :

Lot size : 0.5 mL

Concentration : 100 ng /mL

Presentation : Supplied as a liquid in 0.01M Tris , 0.15M NaCl,
1% BSA, 0.1% Sodium azide, pH 7.5

Background: Vitellogenin, the egg yolk precursor protein, is induced by estrogens.
On the other hand, environmental estrogens are known to disrupt endocrine system in animals. Recently, vitellogenin has proven to be an ideal marker for environmental estrogens.

Source : Ascites of 17β -estradiol treated female medaka

Appearance : Clear liquid

Applications : ELISA

Storage : below -20°C for at least 6 months
(recommendation: at -80°C for long term storage)
Avoid freeze thaw cycles

Note : After preparation of working standards, these solutions should be stored at -20°C for up to two weeks.

Warning : For research use only
Contains sodium azide as preservative

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

e-mail : export@cosmobio.co.jp

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618