

POLYCLONAL ANTIBODY

For research use only, Not for diagnostic use

Catalog No. 106

Anti Androstenedione-11 α

Product type	Primary antibody
Immunogen	Androstenedione-11 α -Succinate-BSA
Host Species	Rabbit
Source	Serum
Purification	-
Form	Lyophilized form. (Whole serum) * Reconstitute with distilled water. No preservative is contained.
Volume	2000 tests ** 1 TEST: [250 μ l / Tube use] in RIA experiment
Label	Unlabeled
Specificity	Androstenedione-11 α
Cross reactivity	Human, Bovine, Equine, Other mammalian and Pisces Other species have not been tested.
Storage Conditions	Store at 4°C. (below -20°C for prolonged storage) After reconstitution, it should be divided into small quantity to avoid freezing and thawing.

Application notes • **RIA** : 2000 tests (250 μ l / tube use / 1 test)

*** Dilution buffer and Optimal dilutions/concentrations should be determined by the end user
Other applications have not been tested.

[ANTIBODY CHARACTERIZATION]

The Cross reactivities of Anti Androstenedione-11 α (Cat# 106) by RI Assay

Steroids	Cross Reaction
Androstenedione	100.00 %
Androstenetrione	200.00 %
5 α -Androstanedione	32.30 %
DHA	5.90 %
Progesterone	3.90 %
Testosterone	2.90 %
Cortisol	0.25 %
Corticosterone	0.10 %
Estradiol	0.01 %

Related products

Product Name	Maker	Cat. Num.	Quantity
Testosterone-3(E)-HRP	FKA	101	1000 TEST
Anti Testosterone-3(E) Polyclonal Antibody	FKA	102	2000 TEST
Anti Testosterone-3(E) Polyclonal Antibody (Purified)	FKA	102E	2000 TEST
Testosterone-11 α -HRP	FKA	103	1000 TEST
Anti Testosterone-11 α Polyclonal Antibody	FKA	104	2000 TEST
Anti Testosterone-11 α Polyclonal Antibody (Purified)	FKA	104E	2000 TEST
DHT-11 α -HRP	FKA	111	1000 TEST
Anti DHT-11 α Polyclonal Antibody	FKA	112	2000 TEST
Anti DHT-11 α Polyclonal Antibody (Purified)	FKA	112E	2000 TEST
11-Ketotestosterone (11-KT)	FKA	117ST	10 MG
11-Oxo-Testosterone-3-HRP	FKA	117	1000 TEST
Anti 11-Oxo-Testosterone-3 Polyclonal Antibody	FKA	118	2000 TEST
Anti 11-Oxo-Testosterone-3 Polyclonal Antibody (Purified)	FKA	118E	2000 TEST
11 β -OH-Testosterone-3-HRP	FKA	147	1000 TEST
Anti 11 β -OH-Testosterone-3 Polyclonal Antibody (Purified)	FKA	148E	2000 TEST
Anti Androstenedione-11 α Polyclonal Antibody	FKA	106	2000 TEST
Anti Androstenedione-11 α Polyclonal Antibody (Purified)	FKA	106E	2000 TEST
16 α -OH-Androstenedione-3-HRP	FKA	123	1000 TEST
Anti 16 α -OH-Androstenedione-3 Polyclonal Antibody (Purified)	FKA	124E	2000 TEST
Androstenedione-3-HRP	FKA	137	1000 TEST
Anti Androstenedione-3 Polyclonal Antibody	FKA	138	2000 TEST
Anti Androstenedione-3 Polyclonal Antibody (Purified)	FKA	138E	2000 TEST
11 β -OH-Androstenedione-HRP	FKA	139	1000 TEST
Anti 11 β -OH-Androstenedione Polyclonal Antibody (Purified)	FKA	140E	2000 TEST
4-Androstenetrione-3-HRP	FKA	149	1000 TEST
Anti 4-Androstenetrione-3 Polyclonal Antibody (Purified)	FKA	150E	2000 TEST
Anti 5 α -Androstane-3 α , 17 β -diol Polyclonal Antibody	FKA	130	2000 TEST
Anti 5 α -Androstane-3 α , 17 β -diol Polyclonal Antibody (Purified)	FKA	130E	2000 TEST
5 α -Androstane-3 α , 17 β -diol-3-Glucuronide-HRP	FKA	141	1000 TEST
Anti 5 α -Androstane-3 α , 17 β -diol-3-Glucuronide Polyclonal Antibody (Purified)	FKA	142E	2000 TEST
5 α -Androstane-3 α , 17 β -Diol-17-Glucuronide-HRP	FKA	143	1000 TEST
Anti 5 α -Androstane-3 α , 17 β -Diol-17-Glucuronide Polyclonal Antibody (Purified)	FKA	144E	2000 TEST
Anti 5 α -Androstane-3 α , 11 α , 17 β -Triol-17-Glucuronide-11-Succinate Polyclonal Antibody (Purified)	FKA	146E	2000 TEST
DHA-11 α -HRP	FKA	107	1000 TEST
Anti DHA-11 α Polyclonal Antibody	FKA	108	2000 TEST
Anti DHA-11 α Polyclonal Antibody (Purified)	FKA	108E	2000 TEST
DHA-3-HRP	FKA	109	1000 TEST
Anti DHA-3 Polyclonal Antibody	FKA	110	2000 TEST
Anti DHA-3 Polyclonal Antibody (Purified)	FKA	110E	2000 TEST
16 α -OH-DHA-3-HRP	FKA	125	1000 TEST
Anti 16 α -OH-DHA- Polyclonal Antibody	FKA	126	2000 TEST
Anti 16 α -OH-DHA- Polyclonal Antibody (Purified)	FKA	126E	2000 TEST

For research use only, Not for diagnostic use.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619