

Anti Dehydroepiandrosterone-3-Succinate polyclonal antibody

Product No. FKA-110

The Cross reactivities of anti Dehydroepiandrosterone-3-Succinate (FKA-110) by RI Assay

Antigen	Dehydroepiandrosterone-3-Succinate-BSA
Antibody[FKA-110]	lyophilized Serum 50μl / Vial
Dilution of Serum	2000test for RI Assay 250μl / tube use
Binding of Total	60%
DHA dose of B0/2	200 pg/ml
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

Steroid	Cross Reaction
Dehydroepiandrosterone (DHA)	100 %
DHA-3-Sulfate	41.50 %
4-Androstenedione	85.00 %
Androsterone	18.00 %
Etiocholanolone	1.90 %
Testosterone	0.80 %
5-Androstenediol	1.50 %
16Alpha-Hydroxy-DHA	6.50 %
Dihydrotestosterone(DHT)	3.80 %
Progesterone	0.30 %
Pregnenolone	0.27 %

For research use only; not for use as a diagnostic.