



Anti 6-Oxo-Estrone-6-CMO-BSA polyclonal antibody

Product No. FKA-202E

The Cross reactivities of anti 6-Oxo-Estrone-6-CMO-BSA (FKA-202E)
by ELISA with 6-Oxo-Estrone-6-CMO-HPR (FKA-201)

Antigen	6-Oxo-Estrone-6-CMO-BSA		
Antibody(FKA-202E)	Purified by Protein A	lyophilized	3μl / Vial
Dilution of IgG	2000test for ELISA	50μl/well use	
Label (FKA-201)	6-Oxo-Estrone-6-CMO-HPR		1ml / Vial
	1,000test for ELISA	50μl/well use	
O.D. (492/600nm) of B0	1.0		
Estrone dose of B0/2	5 ng/ml		
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.		

Steroid	Cross Reaction
Estrone(E1)	100 %
Estradiol or Estradiol-17Beta(E2)	5.00 %
Estriol(E3)	1.20 %
Estrone Sulfate	0.40 %
Estrone Glucuronide	0.15 %
2-Methoxyestrone	0.05 %
Pregnenolone	0.00 %
Progesterone	0.00 %
Cortisol	0.00 %
Cortisone	0.00 %
4-Androstenedione	0.00 %
Testosterone	0.00 %
DHA	0.00 %
DHT	0.00 %
Progesterone	<0.01 %
Coltisol	<0.01 %
Corticosterone	<0.01 %

For research use only; not for use as a diagnostic.

