



Anti Estradiol-3-carboxymethylether polyclonal antibody

Product No. FKA-236

The Cross reactivities of anti Estradiol-3-carboxymethylether (FKA-236) by RI Assay

Antigen	Estradiol-3-carboxymethylether-BSA	
Antibody[FKA-236]	lyophilized Serum	25μl / Vial
Dilution of Serum	2000test for RI Assay	250μl / tube use
Binding of Total	68%	
Estradiol dose of B0/2	200 pg	
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.	

Steroid	Cross Reaction
Estradiol(E2)	100 %
Estrone(E1)	5.20 %
Estriol(E3)	1.40 %
Estradiol-3-Sulfate	31.50 %
Estradiol-3-glucuronide	57.50 %
Estrone-3-Sulfate	2.70 %
Estrone-3-glucuronide	3.00 %
4-Androstenedione	0.60 %
Testosterone	0.48 %
Progesterone	0.11 %
DHA	0.07 %
17Alpha-OH-progesterone	0.02 %
Pregnenolone	0.01 %
Cortisol	0.01 %
Cortisone	0.01 %
DHT	0.01 %

For research use only; not for use as a diagnostic.

