

The RI Assay of anti Progesterone-3(E)-CMO-BSA according to FKA-302

The Cross reactivities of anti Progesterone-3(E)-CMO-BSA (FKA-302) by RI Assay with Progesterone(1.2.6.7-³H)

Antigen	Progesterone-3(E)-CMO-BSA		
Antibody	FKA-302	Serum	lyophilized serum of 25µl / Vial
Dilution of Serum	2000	test for RI Assay	250µl/tube use
Binding of Total	60%		
Progesterone of B0/2	50 pg		
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.		

Steroid	Cross Reaction
Progesterone	100 %
5 Alpa-pregnanedione	12.50 %
Pregnenolone	2.00 %
20 Alpa-hydroxyprogesterone	0.20 %
17 Alpa-hydroxyprogesterone	0.01 %
etc.	0.01 %

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20,TOYO 2CHOME,KOTO-KU,TOKYO 135-0016,JAPAN
 TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp