

The ELISA of anti Ecdysone-6-CMO-BSA serum according to FKA-612E and FKA-611

The Cross reactivities of anti Ecdysone-6-CMO-BSA serum (FKA-612E) by ELISA with Ecdysone-6-CMO-HRP (FKA-611)

Antigen	Ecdysone-6-CMO-BSA serum	
Antibody	FKA-612E	lyophilized serum of 0.2 µl / Vial
Dilution of IgG	2000test for ELISA	50µl/well use
Label (FKA-611)	Ecdysone-6-CMO-HRP	1ml / Vial
	1,000test for ELISA	50µl/well use
O.D. (492/600nm) of B0	about 1.0	
Ecdysone dose of B0/2	about 1 ng/ml	
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.	

Steroid

Cross Reaction

Ecdysone	100 %
20beta-OH-ecdysone	1.50 %
Cholesterol	0.00 %
Testosterone	0.00 %
Pregnenolone	0.00 %
Progesterone	0.00 %
Cortisol	0.00 %

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20,TOYO 2CHOME,KOTO-KU,TOKYO 135-0016,JAPAN
 TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp