



Code No. KAL-KR051-EX


For research use only

Anti Human Aurora-A Polyclonal Antibody, Rabbit

Aurora-A (Aurora-A/STK-15/Aik) encoding a centrosome-associated kinase, is amplified and overexpressed in many cancer cell types, and is involved in the induction of centrosome duplication-distribution abnormalities and aneuploidy in mammalian cells. So that, Aurora-A is a potential oncogene. Aurora family have been identified three types (Aurora-A, B, C) which has high homology domain in C terminal region. Aurora family is indispensable kinase that regulate the structure and function of centrosomes and spindle. Aurora-A could be one of targeting molecule of anti-cancer therapy.

This antibody was established from the purified serum of the immunized rabbit with N terminal region of Aurora-A and did not show any cross-reactivity to Aurora-B and Aurora-C. This antibody is useful for western blotting and immunohistochemistry.

Package Size	25 μ g (100 μ L/vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	Block Ace as a stabilizer, containing 0.1% Proclin as a bacteriostat
Storage	Below -20°C until needed.
Purification method	This antibody was purified from rabbit serum by affinity chromatography.
Working dilution	For western blotting ; 0.5~2 μ g/mL For Immunohistochemistry ; 1~2 μ g/mL



green : Aurora-A

blue : DNA

Immunocytochemistry

Sample : HeLa Cell

Preparation of antibodies and instruction

Dr.Saya H. and Dr.Marumoto T. at Department of Tumor Genetics and Biology, Kumamoto University School of medicine.



Anti Human Aurora-A Polyclonal Antibody, Rabbit

【Reference】

- 1.Giet R, Prigent C.:J Cell Sci 1999 Nov;112(Pt21):3591-601
- 2.Bischoff JR,et al.:EMBO J 1988 Jun 1;17(11):3052-65
- 3.Zhou H,et al.:Nat Genet 1988 Oct;20(2):189-93
- 4.Lim Shen Kiat,et al.:J Biol Chem 2002 Mar 22;277(12):10719-10726
- 5.Goepfert TM,et al.:Cancer Res 2002 Jul 15;62(14):4115-22

Distributor



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer

 Trans Genic Inc.

7-1-6 Minatojimmaminami-machi,
Chuo-ku, Kobe,650-0047 JAPAN
TEL : +81-78-306-0590
FAX : +81-78-306-0589
URL : <http://www.transgenic.co.jp/>
e-mail : techstaff@transgenic.co.jp