

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. 2ZCKHH3

Anti Human HGF Beta chain

Product type Primary antibodies

Host Mouse **Source** Ascites

Hybridoma cell line (Hyb- T7156) was grown in ascitic fluid of BALB/c Mouse

Form Liquid

Ammonium sulfate precipitation

Physiological salt solution containing 0.1% Sodium Azide

Volume 0.5 mL Concentration 0.5 mg/mL

Specificity This antibody specifically reacts with Human HGF Beta chain

Antigen Human HGF synthetic peptide antigen

Clone Hyb-T7156

The splenic lymphocytes from BALB/c mouse immunized with Human HGF synthetic peptide were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-T7156) were

screened.

Isotype IgG₁

Application notes ELISA, IP, WB

Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Maintain at 2-8°C (-20°C for long term storage) in undiluted aliquots.

Avoid repeated freeze/thaw cycles.

References



Cosmo Bio Co., Ltd.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: <u>export@cosmobio.co.jp</u> [Outside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: +81-3-5632-9610

FAX: +81-3-5632-9618 FAX: +81-3-5632-9619



モノクローナル抗体

For research use only. Not for clinical diagnosis.

Catalog No. 2ZCKHH3

抗ヒト HGF β 鎖

Host マウス

Source BALB/c マウス腹水

ハイブリドーマ(Hyb-T7156)を BALB/c に接種して得られた腹水

Form 硫酸アンモニウム塩析精製

0.1% アジ化ナトリウムを含む生理的食塩水溶液

Volume 0.5 mL Concentration 0.5 mg/mL

Specificity ヒト HGF β 鎖に特異的に反応 Antigen ヒト HGF の合成ペプチド抗原

Clone Hyb- T7156

ヒト HGF の合成ペプチドを免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞

(NS-1)を融合して得られたハイブリドーマ

Isotype IgG₁

Application notes ELISA, IP, WB

Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage 2~8°C 3ヶ月間

-20°C 3年間

(凍結融解の繰り返しは避けてください)

References



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: +81-3-5632-9610 FAX: +81-3-5632-9618 FAX: +81-3-5632-9619