

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. 2ZCP9

Anti HCV Core (CP9)

Product type	Primary antibodies
Host	Mouse
Source	Ascites
Form	Hybridoma cell line (Hyb-9380B) was grown in ascitic fluid of BALB/c Mouse Liquid Ammonium sulfate precipitation Physiological salt solution containing 0.1% Sodium Azide
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	This antibody specifically reacts with HCV core region (Amino acid 39-74)
Antigen	HCV core region synthetic peptide(CP9)
Clone	Hyb-9380B The splenic lymphocytes from BALB/c mouse immunized with HCV core region synthetic peptide (CP9) were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-9380B) were screened.
Isotype	IgG ₁

Application notes ELISA
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Maintain at 2-8 °C (-20 °C for long term storage) in undiluted aliquots.
Avoid repeated freeze/thaw cycles.

References



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619



モノクローナル抗体

For research use only. Not for clinical diagnosis.

Catalog No. 2ZCP9

抗 HCV コア (CP9)

Product type	一次抗体
Host	マウス
Source	BALB/c マウス腹水 ハイブリドーマ (Hyb-9380B) を BALB/c マウスに接種して得られた腹水
Form	硫酸アンモニウム塩析精製 0.1%アジ化ナトリウムを含む生理的食塩水溶液
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	HCV コア領域 (アミノ酸配列 39~74) に特異的に反応
Antigen	HCV コア領域合成ペプチド (CP9) 抗原
Clone	Hyb-9380B HCV のコア領域の合成ペプチド (CP9) を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞(NS-1)を融合して得られたハイブリドーマ
Isotype	IgG ₁

Application notes

ELISA

Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage

2~8°C 3 ヶ月間
-20°C 3 年間
(凍結融解の繰り返しは避けてください)

References



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619