

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHCMA3

Anti Human AFP (Antigenic determinant c)

Product type	Primary antibodies
Host	Mouse
Source	Ascites
Form	Hybridoma cell line (Hyb-2065) was grown in ascitic fluid BALB/c Mouse Liquid Ammonium sulfate precipitation Physiological salt solution containing 0.1% Sodium Azide
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	This antibody specifically reacts with Human AFP (determinant c)
Antigen	Purified human alpha fetoprotein (AFP) antigen
Clone	Hyb-2065 The splenic lymphocytes from BALB/c mouse immunized with purified human alpha fetoprotein (AFP) were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-2065) were screened.
Isotype	IgG ₁

Application notes ELISA
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Maintain at 2-8 °C (-20 °C for long term storage) in undiluted aliquots.
Avoid repeated freeze/thaw cycles.

References



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

モノクローナル抗体

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHCMA3

抗ヒト AFP (抗原決定基 c)

Product type	一次抗体
Host	マウス
Source	BALB/c マウス腹水 ハイブリドーマ (Hyb-2065) を BALB/c マウスに接種して得られた腹水
Form	液状 硫酸アンモニウム塩析精製 0.1%アジ化ナトリウムを含む生理的食塩水溶液
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	ヒト AFP (抗原決定基 c) に特異的に反応
Antigen	精製ヒト α フェトプロテイン (AFP) 抗原
Clone	Hyb-2065 精製ヒト α フェトプロテイン (AFP) を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得られたハイブリドーマ
Isotype	IgG ₁

Application notes ELISA
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage 2~8°C 3 ヶ月間
-20°C 3 年間
(凍結融解の繰り返しは避けてください)

References



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619