



Anti PHI (Rat) Serum **Cat. No. Y022** **Lot No. 025171011**

Description: This antiserum was raised in a rabbit by immunization with a carrier free synthetic PHI (rat) peptide. The product vial contains 50 µL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution . It can be used for immunoassay, immunohistochemistry or any other immunoreaction with PHI (rat).

Immunogen: Synthetic PHI (rat), carrier free **Host:** Rabbit

Amino Acid Sequence of PHI (rat)¹⁾: HADGVFTSDY SRLLGQISAK KYLES LI-NH₂

Product Form: Lyophilized unpurified serum **Size:** 50 µL

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon recon- stitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing- thawing should be avoided.

Suggested Working Dilution Range: 1:2,000-5,000 (final dilution ~1:35,000) for radioimmunoassay;
1: 1,000-4,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

Specificity (based on radioimmunoassay): PHI (rat) 100%, PHI (human) 80%, PHI (porcine) 75%, secretin (porcine) 0%, VIP (porcine) 0%, PACAP38 (human) 0%

Positive Control (immunohistochemistry): Rat colon.

Species Tested: Rat

REFERENCES:

1) M. Nishizawa, N. Yanaihara et al., Nucleotide sequence divergence and functional constraint in VIP precursor mRNA evolution between human and rat. FEBS Letter 183:55-59, 1985

FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



2480-1 Awakura, Fujinomiya-shi,
Shizuoka 418-0011, JAPAN
Phone: +81-544-22-2771
Fax: +81-544-22-2770
URL : www.yanaihara.co.jp
E-mail : ask@yanaihara.co.jp