



## Anti Gastrin 34 (1-15) (Human) Serum

Cat. No. YII-Y110-EX      Lot No. 023071109

**Description:** This antiserum was raised in a rabbit by immunization with a bovine serum albumin (BSA) conjugate of synthetic gastrin 34 (G34) (1-15) (human) peptide. The product vial contains 50 µL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with gastrin (human).

**Immunogen:** Synthetic gastrin 34 (1-15)(human)-BSA conjugate **Host:** Rabbit

**Amino Acid Sequence of Gastrin 34 (1-15) (human)<sup>1)</sup>:**

QLGPQGPPHL VADPSKKQGP WLEEEEEAYG WMDF

**Product Form:** Lyophilized unpurified serum **Size:** 50 µL

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN3 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

**Suggested Working Dilution Range:** 1:5,000 (final dilution~1:45,000) for radioimmunoassay; 1:500-2,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

**Specificity** (based on radioimmunoassay): G34 (1-15) (human) 100%, G34 (human) 43%, G34 (1-5) (human) 30%, G34 (6-15) (human) 0%

**Positive Control** (immunohistochemistry): Rat antrum

**Species Tested:** Rat

### REFERENCES:

1) K. Kato, Y. Hayashizaki et al., Molecular cloning of the human gastrin gene. Nucleic Acids Research 11: 8197-8205, 1983

### FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

#### Distributor



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

#### Manufacturer



2480-1 Awakura, Fujinomiya-shi,  
Shizuoka 418-0011, JAPAN  
Phone: +81-544-22-2771  
Fax: +81-544-22-2770  
URL : [www.yanaihara.co.jp](http://www.yanaihara.co.jp)  
E-mail : [ask@yanaihara.co.jp](mailto:ask@yanaihara.co.jp)